

Techno-Vision LED Illuminators

Make the switch to high-tech lighting! LED Radiology Viewbox with Dimmer and Film Activating Switch. This illuminator will hold two 14 x 17 x-ray films.



Cat# V2-LED
Double Bank with dimmer and Film Activated Switch



DIMENSIONS: FRAME	895 x 496 x 39
DIMENSIONS: VIEWING AREA	711 x 422 mm
AVERAGE BRIGHTNESS	Approx. 4,500 Cd/m ² + 10%
AVERAGE UNIFORMITY	Over 90%
LIGHT SOURCE	LED (Light Emitting Diode)
COLOR TEMPERATURE	9,600 K
LAMP LIFETIME	Approx. 100,000 Hours
INPUT POWER	110Vac (50Hz) - 240Vac (60Hz), Free Voltage
POWER CONSUMPTION	Max 25W
POWER SUPPLY	Built-In
OPERATING TEMPERATURE	10 degree C—40 degree C
STORAGE TEMPERATURE	-15 degree C—60 degree C
FILM CLIP	Rolling Bar
MOUNTING TYPE	Wall Mount
SAFETY AGENCY	C.E.
WEIGHT	12.3 lbs

- Other sizes available
- Desk stand available
- Vesa mount available

