



2-Step Platform

The unique design of the *2-Step Platform* permits efficient and accurate X-rays of knees, ankles, feet, and the lower leg (lateral-medial, medial-lateral, or anterior-posterior).

The platform places the patient in the optimal position for easy capture of the highest quality weight-bearing X-ray images. Safety grab bars provide stability for patients, both young and old.

- Designed for patients weighing up to 450 lbs / 204 kg
- Keeps patient's foot as close to the image plate as possible to ensure high-quality images
- Three imaging panel placement locations provide great versatility in posing the patient for optimal imaging
- Designed for use with bucky or fixed cassette
- **Tip-N-Tow** – Easily portable with built-in rollers and grab bars
- Easy to keep clean, the platform is made from germ-resistant, non-porous polyethylene, and stainless steel
- Models for all film, CR, and DR based systems



Free Standard Shipping
All Products - 48 States

Clear Image Devices LLC
5225 68th St., Suite 8 Caledonia, MI 49316
Phone: +1 734-474-6537
Fax: +1 734.661.0220



www.ClearImageDevices.com

copyright ©2016 Clear Image Devices LLC
All rights reserved.

2-STEP PLATFORM	
Item #	Description
24302-C50	For DR Systems: Canon CXDI-50, 60
24302-C55	For CR Systems and the following DR Systems: Canon CXDI-55, 70, 80, 401C, 501, 701C, 801C, Carestream DRX-1, Fuji FDR D-EVO, Konica Minolta Aero DR, AGFA
24302-SY	For DR Systems: Siemens Ysio, GE Flash Pad, Philips Wireless
24302-GE	For DR Systems: GE XR650, GE AMX 700
Dimensions: Top Platform	
Height: 16 in / 40 cm, Width: 19.5 in / 50 cm	
Length: 17 in / 43 cm	
Handle Height: 51 in / 130 cm	
50122	Side rail accessory

Clear Image Devices has become the #1 brand in imaging accessories through innovative product development in cooperation with the medical imaging community. We believe that value begins with quality; quality products and customer service. We are advancing medical imaging every day by bridging the gap between technology and people.