

3-Step Platform



The unique design of our series of step platforms permits efficient and accurate X-rays of knees, ankles, feet, and the lower leg (lateral-medial, medial-lateral, or anterior-posterior).

The platform places the patient in the optimal position for easy capture of the highest quality weight-bearing X-ray images. Safety grab bars provide stability for patients, both young and old.

- *Designed for patients weighing up to 450 lbs / 204 kg*
- *Keeps patient's foot as close to the image plate as possible to ensure high-quality images*
- *Three imaging panel placement locations provide great versatility in posing the patient for optimal imaging*
- *Designed for use with bucky or fixed cassette*
- *Tip-N-Tow – Easily portable with built-in rollers and grab bars*
- *Easy to keep clean, the platform is made from germ-resistant, non-porous polyethylene, and stainless steel*
- *Models for all film, CR, and DR based systems*

3-STEP PLATFORM

Item #	Description
24306-C50	For DR Systems: Canon CXDI-50, 60
24306-C55	For CR Systems and the following DR Systems: Canon CXDI-55, 70, 80, 401C, 501, 701C, 801C, Carestream DRX-1, Fuji FDR D-EVO, Konica Minolta Aero DR, AGFA
24306-SY	For DR Systems: Siemens Ysio, GE Flash Pad, Philips Wireless
24306-GE	For DR Systems: GE XR650, GE AMX 700
Dimensions: Top Platform	
Height: 24 in / 41 cm Width: 19.5 in / 50 cm	
Length: 17 in / 43 cm	
Handle Height: 59 in / 150 cm	

Clear Image Devices has become the #1 brand in imaging accessories through innovative product development in cooperation with the medical imaging community. We believe that value begins with quality; quality products and customer service. We are advancing medical imaging every day by bridging the gap between technology and people.



Free Standard Shipping
All Products - 48 States

