


2 Step Platform



Our weight bearing 2 Step Platform is designed to accommodate digital DR Panels, CR Cassettes and Analog.

This product is designed for easy mobility plus stability while allowing for easy weight bearing imaging on the foot, ankle and knee.

Our 2 Step Platform is constructed of high density ABS plastic and stainless steel.



RC Imaging
your radiology resource

2 Step Platform



- Anti-slip platform
- Comfortable textured platform
- Non-skid rubber feet
- Multiple detector slots for positioning
- Medical grade stainless steel rails for patient stability
- Rear mounted offset casters for easy relocation
- Simple cleaning and sterilization



RC Imaging
your radiology resource

800.528.6311 585.392.4336 585.392.5575
50 Old Hojack Lane, Hilton, New York 14468
www.rcimaging.com

