


Step Platform for U-Arm Systems



The *U-Arm Step* is designed for easy capture of the highest quality weight-bearing lateral and AP foot images. This positioning device is ideal to satisfy the imaging needs of orthopedic surgeons, sports medicine clinics, and general orthopedic practices.

Optimal patient positioning also helps acquire better diagnostic images for professional and collegiate sports physicians, with fewer retakes, thus reducing exposure to radiation.

- *Designed for patients weighing up to 500 lbs / 227 kg*
- *Keeps patient's foot as close to the image plate as possible to ensure high-quality images*
- *Easily portable with four full-locking swivel casters*
- *Easy to keep clean, the platform is made from germ-resistant, non-porous polyethylene, polycarbonate, and stainless steel*

STEP PLATFORM FOR U-ARM SYSTEMS

Item #	Description
24701	500 lb Capacity: Sedecal, Swissray, IMIX, IDC and other U-Arm Systems
Dimensions: Top Platform	
Height: 18 in / 45 cm, Width: 28 in / 71 cm	
Length: 32 in / 81 cm	
Handle Height: 53 in / 135 cm	
Docking Space:	
Height: 15.50 in / 39 cm	
Width: 28.25 in / 71.8 cm	
* Special order for use with CR cassettes or DR panels. Request a quote.	

Clear Image Devices has become the #1 brand in imaging accessories through innovative product development in cooperation with the medical imaging community. We believe that value begins with quality; quality products and customer service. We are advancing medical imaging every day by bridging the gap between technology and people.



Free Standard Shipping
All Products - 48 States

Clear Image Devices LLC
5225 68th St., Suite 8 Caledonia, MI 49316
Phone: +1 734-474-6537
Fax: +1 734.661.0220



www.ClearImageDevices.com

copyright ©2016 Clear Image Devices LLC
All rights reserved.