

# Lightweight. Flexible. Affordable.

### Here's an economical, space-saving solution for small X-ray rooms .

The CARESTREAM SPT2 Radiographic Table is specifically designed for smaller X-ray rooms – or, for any imaging facility with a smaller budget. This highly mobile table allows for easy, precise positioning of patients for efficient radiographic exams – and provides these benefits at a highly affordable price.

### Advantages At a Glance

- Lightweight, for easy portability\*
- Small size to easily position the table over wallstand Buckys
- Economical price to fit virtually any budget

\* Caution: The table is only intended to be used to support patients during radiographic examinations. Do not use the table to transport patients.



# SPT2 Radiographic Table

You can rely on the SPT2 Table for high performance at an exceptionally low cost. Its features and benefits include:

## Easy Handling

Directional wheel castors and a large handgrip allow for easy movement in any direction.

## Easy to Clean

The streamlined table design simplifies cleaning and disinfecting.

## Four Wheel Brakes

A brake on every wheel helps to ensure stability and patient safety, and all four wheels have a pedal-activated locking mechanism.

## Carbon-fiber Surface

The table's carbon-fiber top is composed of a superior radiolucent material.

## Durability

Its exceptionally sturdy design is intended to provide a ten-year service life.



## Product Specifications

### Table Dimension

Length: 83.5in (212.0 cm)

Width: 27.0 in (68.5 cm)

Height: 27.6 in (70.0 cm)

### Table Clearance from Floor

Height: 26.69 in (67.8 cm)

### Table Thickness

Thickness: 1.65 in (4.2 cm)

### Table Weight

Weight: 92.6 lb (42.0kg)

### Table Patient Weight Limit

650 lb (295 kg)

### X-ray Transparent Surface

Length: 74.8 in (190.0cm)

Width: 22.0 in (56.0 cm)

### X-ray Absorption Factor

Factor  $\leq$  0.7mm Aluminum Equivalents at 100 kVp

## Environment Requirements

### Transit and Storage

The receiving and storage area(s) must be dry and provide the proper humidity and temperature controls requested for this table.

### Storage Temperature

-9 to 150°F (-23 to 66°C)

### Relative Humidity

35% to 75% (non-condensing)

### Atmospheric Pressure

70 to 106kPa

### Operating Temperature

50 to 104°F (10 to 40°C)

[carestream.com](http://carestream.com)



"Rx only"

© Carestream Health, Inc., 2018. CARESTREAM is a trademark of Carestream Health. CAT 2000 254 9/18

**Carestream**