

December 2011

Low Price
Guarantee

MRIequip.com

Non-Magnetic Equipment

Now at your fingertips is a diverse line of non-ferrous and PVC products for your MRI environment



Many Items Now Available
with **GREEN** Upholstery

Sanitary Headset Covers

See page 63



NEW

Step Stools

See pages 17-18



Oxygen Equipment

See pages 34-39



Visit www.mriequip.com for many more products and options.



What can we do for you?

Upon initial contact we will customize a package based upon the procedures performed as well as the strength, weight capacity and brand of MR system(s) used.

We will take into account specific product requests and determine the most cost effective solution to our customers needs.

We present a quotation on a full custom package including freight, rather than a pre-packaged kit that is incomplete and contains items that may go unused.

We specialize in custom non-ferrous products that may require additional time for production. We will provide the necessary lead time and discuss a final purchase order date that will allow for delivery before a grand opening or re-opening.

Once a purchase order has been received, upon the customers request; we may hold all products to be delivered at one time or have it delivered as available. Products in stock at the time of order will be immediately tagged for that order and held until the release date. This eliminates the possibility of a stock item at order date becoming out of stock on the release date.

We supply MSDS sheets on pertinent products for our customer's records.

Two days before release date, we will contact our customer to verify they're ready for the delivery. If there is a delay in the facilities opening, we will hold the product until needed at no additional charge.

Above all else, we will take that extra step in providing our customers with knowledgeable advice, low prices, on-time delivery and excellent customer service. We recognize that each of these attributes are vital to a smooth transaction and a satisfied customer.

Before ordering, please read our company terms and policies on inside back cover.

Safety is our Business!

Our dedication to safety begins at 7:00 AM every day of the work week - it is our livelihood. All products in our warehouse are checked on entry and then re-checked on exit. Products drop shipped from our manufacturers directly to your door are brought in randomly and re-tested for safety verification. We continually work with all manufacturers to ensure safety and proper designation.

For more information on the ASTM Standard visit www.ASTM.org and search for F2503-05.

For more helpful information pertaining to MR Safety, you can visit

www.MRI-Planning.com

www.MRISafety.com

www.radiology.upmc.edu/MRsafety

The highly anticipated Standard Practice for Marking Medical Devices for Safety in the Magnetic Resonance Environment by American Society for Testing and Materials (ASTM) has arrived.

As recommended by ASTM Intl; the intent of this practice is to provide needed information about the safety of items in and near MR scanners using a compact and easily recognized set of symbols and terms. The terms MR safe and MR compatible as first defined in 1997 in the FDA draft guidance document, "A Primer on Medical Device Interactions with Magnetic Resonance Imaging Systems," have been used to describe the safety of devices in and near MR systems. However, there has been a great deal of confusion surrounding this terminology. An incorrect assumption is often made. Users often incorrectly assume that items labeled MR safe or MR compatible are safe or compatible for any MR environment. MR environments vary in terms of magnetic field strength and RF conditions. This new system eliminates much of that confusion.

Designation

MR Safe—an item that poses no known hazards in all MR environments.

MR Conditional—an item that has been demonstrated to pose no known hazards in a specified MR environment with specified conditions of use. Field conditions that define the specified MR environment include field strength, spatial gradient, dB/dt (time rate of change of the magnetic field), radio frequency (RF) fields, and specific absorption rate (SAR). Additional conditions, including specific configurations of the item, may be required.

MR Unsafe—an item that is known to pose hazards in all MR environments.

Recommended Icon

Catalog Designation

Labels available on Page 77



■ MR Safe



▲ MR Conditional



● Not MR Safe

Products **NOT** marked with a designation have been tested under the past methods and are currently being re-evaluated under the new ASTM Standards. Products will be updated continually.



Now, at your fingertips, is a diverse line of Non-Ferrous and PVC products for your MRI environment. MRlequip.com was established to provide the highest quality MRI compatible equipment with price tags significantly lower than the current marketplace.

In the past, it was a priority to clear the room of all Ferrous materials that pose an immediate risk. Now, a majority of those items can be replaced, in extent saving time and creating a friendlier environment for both the patients and the providers.

Beyond the standard items, such as Wheelchairs, I.V. Poles and Stools, a whole new world of products has opened up. Ranging from Sound Systems to Furniture, these products increase patient comfort while decreasing concerns for injury and liability.

At MRlequip.com, we are constantly researching new products and welcome any suggestions you might have. If you can't locate what you require, let us know. Our design teams are ready for their next challenge.

Index

Aluminum Stretcher & Acces.——1	Medication Cart Accessories——23	Stereotatic Table Drapes——83
Aluminum Wheelchair & Acc.——2	Medication Carts——21-22	Stethoscopes——31
Biopsy Tray——51	Mop Bucket——86-87	Stools——17-19
Blood Pressure——33	Multi Use Carts——25-26	Storage Bins——43
Breast Coil Draping——49-50	Multi-Modality Markers——41	Storage Cases——56
Ceiling Tile Art——57	Name Badge Accessories——80	Stretchers / Gurneys——9-13
Clipboard——56	Needle——51	Stretchers / Gurneys Acces.——10, 12
Coil Carts——48	Oxygen Carts——35	Suction——39-40
Coils——47	Oxygen Cylinder Carts——37-38	Supply Carts——45-46
Crash Cart——27	Oxygen Cylinder Stands——37-38	Titanium Forceps——52
CT Modality——81-82	Oxygen Cylinders——35	Titanium Tools——89-93
Ear Plugs / Muffs——65-66	Oxygen Flowmeters——34	Titanium Tweezers——53
Electrodes——31	Oxygen Mounting——36	TMJ Device——41
Equipment Stands——36	Oxygen Regulators——34	Tool Kits——90, 94
Ferromagnetic Detection——71	Padlocks——56	Tools——88-94
Fiducial Markers——42	Paper Clips——56	Transfer Belt——7
Fire Extinguishers——67-68	Patient Transfer Boards——13-15	Transport Wheelchairs——7
Floor Safety Cones——79	Pens——56	Trash Cans——93
Floor Signs——86	Positioners——59-60	Utility Carts——27-28
Furniture——20	Prescription Goggles——58	Utility Tables——16
GaussAlert Alarm System——70	Prism Glasses——58	Ventilator Stands——36
Glove Box Holders——44	Pulse Oximeter——32	Walkers——7
Hampers——24	Safety Carpets / Floor——74, 76	Water Basin——53
Hand Sanitizer——66	Safety Manual——70	Wheelchair Accessories——8
Instruments ,Disposable——52	Safety Signs——72-79	Wheelchair Parts——8
IV Clamp——55	Safety Stickers——72-78	Wheelchairs——3-7
IV Stands——36	Sand Bags——60	Wheelchairs, Bariatric——6
Janitorial Supplies——85-87	Sanitation——84-85	Zone Signs——73
Ladders——95	Scalpel——52	
Lanyards——80	Scissors——52	
Laryngoscope——54	Screening Magnet——69	
Laryngoscope Blades——54	Sleep Masks——58	
Lights——29	Sound System Headsets——63-64	
Locking Supply Cart——45	Sound Systems——61-62	
Magna Cleaning Wand——84	Sphygmomanometers——33	
Mayo Stand——53	Step Stools——17-18	

Specialized MRI Catalogs!

Call to request
your copy today or
e-mail your
request to
cs@mrlequip.com.

Breast Imaging Products



Pediatric Products



Disposable Products





Now Available with
GREEN Upholstery



Aluminum Non-Ferromagnetic Gurney

This Non-Ferromagnetic Gurney consists of materials such as Brass, Titanium, Aluminum and Plastic. It has been tested and certified to pass through the pre-screen and pre-entry door systems.

Features:

- Aluminum construction
- Drop down side rails
- 4" Ceramic Total Lock Casters
- Full surround rubber bumper
- Adjustable headrest
- 3" pad included

Specifications:

- Height 28½"
- Height with Rails up 40"
- Height of Rails 15"
- Length 77¼"
- Rail Length 39"
- Width: 27"
- Weight: 130 lbs.
- Back Adj. Range 0-75 Degrees
- Number of Positions: 11

Weight Capacity: 450 lbs

Warranty

2 years on frame
90 days on soft goods and moving parts

FERRALERT™ HALO
FERROMAGNETIC DETECTION SYSTEMS
Item Passes Ferromagnetic Detection
See Page 71 Non-Ferromagnetic
Detection System

Non-Ferromagnetic Aluminum IV Pole for Stretcher

Model	Description	Price
ST-5100	▲ Aluminum IV Pole	\$180.00 ea.



Non-Ferromagnetic Aluminum Oxygen Pole for Stretcher

Model	Description	Price
ST-5101	▲ Aluminum O2 Pole	\$199.00 ea.



Model	Description	Price
ST-5000	▲ Aluminum Gurney	\$3,950.00 ea.

■ MR Safe ▲ MR Conditional ● Not MR Safe



20" Wide Non-Ferromagnetic Aluminum Wheelchair

This Non-Ferromagnetic Wheelchair consists of materials such as Brass, Titanium, Aluminum and Plastic. It has been tested and certified to pass through the pre-screen and pre-entry door systems. This is the first of its kind to have 24" wheels, which allows the patient to move independently.

Flip-up Desk Length Padded Arms
Swing-away Footrests
Color: Burgundy

Specifications

Seat Width: 20"
Overall Width: 28"
Overall Depth: 36" w/footrest
Seat Depth: 16"
Overall Height: 36"
Back Height: 19"
Front Casters: 8" x 1" solid
Rear Wheel: 24" solid
Wheel Locks: push to lock
Weight: 31 lbs w/footrest

Weight Capacity: 300 lbs

Warranty

5 years on frame and crossbraces
6 months on soft goods and moving parts



FERRALERT HALO
FERROMAGNETIC DETECTION SYSTEMS
Item Passes Ferromagnetic Detection



Back of upholstery



FERRALERT HALO
FERROMAGNETIC DETECTION SYSTEMS
Item Passes Ferromagnetic Detection

Non-Ferromagnetic Aluminum IV Pole

Model	Description	Price
WC-2100 ▲	Aluminum IV Pole	\$144.50 ea.

Non-Ferromagnetic Aluminum IV/O2 Pole

Model	Description	Price
WC-2101 ▲	Aluminum IV/O2 Pole	\$203.00 ea.

Non-Ferromagnetic Aluminum Oxygen Pole

Model	Description	Price
WC-2102 ▲	Aluminum Oxygen Pole	\$110.45 ea.





Back View

18" Wide Wheelchair, With Fixed Footrest

Removable Full Length Padded Arms
 Fixed Footrests

Specifications

Seat Width: 18"
 Overall Width: 26"
 Overall Depth: 44" w/footrest
 Seat Depth: 16"
 Overall Height: 35"
 Back Height: 18"
 Front Casters: 8" x 1" solid
 Rear Wheel: 24" solid
 Wheel Locks: push to lock
 Weight: 44 lbs w/footrest

Weight Capacity: 250 lbs



Now Available with
GREEN Upholstery

Warranty

5 years on frame and crossbraces
 6 months on soft goods and moving parts

Model	Color / Description	Price
WC-1000	▲ Burgundy Chair With Footrests	\$1,398.00 ea.
WC-5000	▲ Green Chair With Footrests	\$1,498.00 ea.

See page 8 for wheelchair parts and accessories



Chair With
 Footrests



Back View

Chair With Legrests

18" Wide Wheelchair, With Flip-Up Arms

Flip-Up Desk Length Padded Arms
 Swingaway, Detachable Footrests or
 Swingaway, Elevating, Detachable Legrests

Specifications

Seat Width: 18"
 Overall Width: 26"
 Overall Depth: 31" w/o footrest
 Seat Depth: 16"
 Overall Height: 35"
 Back Height: 18"
 Front Casters: 8" x 1" solid
 Rear Wheel: 24" solid
 Wheel Locks: push to lock
 Weight: 40 lbs w/o footrest

Weight Capacity: 250 lbs



Now Available with
GREEN Upholstery

Warranty

5 years on frame and crossbraces
 6 months on soft goods and moving parts

Model	Color / Description	Price
WC-1001	▲ Burgundy Chair With Footrests	\$1,570.00 ea.
WC-1021	▲ Burgundy Chair With Legrests	\$1,670.00 ea.
WC-5001	▲ Green Chair With Footrests	\$1,670.00 ea.
WC-5021	▲ Green Chair With Legrests	\$1,770.00 ea.

See page 8 for wheelchair parts and accessories



20" Wide Wheelchair, With Full Length Arms

Swingaway, Detachable Footrests or
Swingaway, Elevating, Detachable Legrests

Specifications

Seat Width: 20"
Overall Width: 28"
Overall Depth: 33" w/o footrest
Seat Depth: 16"
Overall Height: 38"
Back Height: 18"
Front Casters: 8" x 1" solid
Rear Wheel: 24" solid
Wheel Locks: push to lock
Weight: 46 lbs w/o footrest

Weight Capacity: 250 lbs

Warranty

5 years on frame and crossbraces
6 months on soft goods and moving parts



Back View

Model	Color / Description	Price
WC-1004	▲ Burgundy Chair With Footrests	\$1,575.00 ea.
WC-1022	▲ Burgundy Chair With Legrests	\$1,670.00 ea.
WC-5004	▲ Green Chair With Footrests	\$1,675.00 ea.
WC-5022	▲ Green Chair With Legrests	\$1,770.00 ea.

See page 8 for wheelchair parts and accessories

22" Wide Wheelchair, With Desk Length Arms

Swingaway, Detachable Footrests or
Swingaway, Elevating, Detachable Legrests

Specifications

Seat Width: 22"
Overall Width: 30"
Overall Depth: 33" w/o footrest
Seat Depth: 18"
Overall Height: 38"
Back Height: 18"
Front Casters: 8" x 1" solid
Rear Wheel: 24" solid
Wheel Locks: push to lock
Weight: 46 lbs w/o footrest

Weight Capacity: 350 lbs

Warranty

5 years on frame and crossbraces
6 months on soft goods and moving parts



Back View

Model	Color / Description	Price
WC-1002	▲ Burgundy Chair With Footrests	\$1,998.00 ea.
WC-1023	▲ Burgundy Chair With Legrests	\$2,198.00 ea.
WC-5002	▲ Green Chair With Footrests	\$2,098.00 ea.
WC-5023	▲ Green Chair With Legrests	\$2,298.00 ea.

See page 8 for wheelchair parts and accessories

■ MR Safe ▲ MR Conditional ● Not MR Safe

Chair With Footrests



Chair With Legrests



Chair With Footrests



Chair With Legrests





Chair With Footrests



Chair With Legrests



Back View

24" Wide Wheelchair, With Desk Length Arms

Swingaway, Detachable Footrests or
 Swingaway, Elevating, Detachable Legrests

Specifications

Seat Width: 24"
 Overall Width: 32"
 Overall Depth: 33" w/o footrest
 Seat Depth: 18"
 Overall Height: 38"
 Back Height: 18"
 Front Casters: 8" x 1" solid
 Rear Wheel: 24" solid
 Wheel Locks: push to lock
 Weight: 46 lbs w/o footrest

Weight Capacity: 350 lbs

Warranty

5 years on frame and crossbraces
 6 months on soft goods and moving parts



Now Available with
GREEN Upholstery

Model	Color / Description	Price
WC-1024	▲ Burgundy Chair With Footrests	\$2,148.00 ea.
WC-1025	▲ Burgundy Chair With Legrests	\$2,348.00 ea.
WC-5024	▲ Green Chair With Footrests	\$2,248.00 ea.
WC-5025	▲ Green Chair With Legrests	\$2,448.00 ea.

See page 8 for wheelchair parts and accessories

Specialized MRI Catalogs!

Pediatric Products



Call to request your copy today or
 e-mail your request to cs@mrlequip.com.



24" Wide Bariatric Wheelchair, With Desk Length Arms

Swingaway, Detachable Footrests or
Swingaway, Elevating, Detachable Legrests

Specifications

Seat Width: 24"
Overall Width: 33"
Seat Depth: 20"
Overall Height: 38"
Back Height: 18"
Front Casters: 8" x 2" solid
Rear Wheel: 24" solid
Wheel Locks: push to lock
Weight: 58 lbs w/o footrest

Weight Capacity: 550 lbs

Warranty

5 years on frame and crossbraces
6 months on soft goods and moving parts



Now Available with
GREEN Upholstery



Back View



Chair With
Footrests

Chair With
Legrests



Model		Color / Description	Price
WC-1029	▲	Burgundy Chair With Footrests	\$2,595.00 ea.
WC-1030	▲	Burgundy Chair With Legrests	\$2,795.00 ea.
WC-5029	▲	Green Chair With Footrests	\$2,695.00 ea.
WC-5030	▲	Green Chair With Legrests	\$2,895.00 ea.

See page 8 for wheelchair parts and accessories

26" Wide Bariatric Wheelchair, With Desk Length Arms

Swingaway, Detachable Footrests or
Swingaway, Elevating, Detachable Legrests

Specifications

Seat Width: 26"
Overall Width: 35"
Overall Depth: 36" w/o footrest
Seat Depth: 20"
Overall Height: 38"
Back Height: 18"
Front Casters: 8" x 2" solid
Rear Wheel: 24" solid
Wheel Locks: push to lock
Weight: 61 lbs w/o footrest

Weight Capacity: 550 lbs

Warranty

5 years on frame and crossbraces
6 months on soft goods and moving parts



Now Available with
GREEN Upholstery



Back View



Chair With
Footrests

Chair With
Legrests



Model		Color / Description	Price
WC-1003	▲	Burgundy Chair With Footrests	\$2,595.00 ea.
WC-1026	▲	Burgundy Chair With Legrests	\$2,795.00 ea.
WC-5003	▲	Green Chair With Footrests	\$2,695.00 ea.
WC-5026	▲	Green Chair With Legrests	\$2,895.00 ea.

See page 8 for wheelchair parts and accessories

■ MR Safe ▲ MR Conditional ● Not MR Safe



**Chair With
 Footrests**



Back View

**Chair With
 Legrests**



20" Wide Transport Wheelchair

Removable Full Length Padded Arms, Swingaway, Detachable Footrests or Swingaway, Elevating, Detachable Legrests

Specifications

Seat Width: 20"
 Overall Width: 25"
 Overall Depth: 27" w/o footrest
 Seat Depth: 16"
 Overall Height: 35.75"
 Back Height: 18"
 Front Casters: 8" x 1" solid
 Rear Wheel: 12.50" x 2.50"
 Wheel Locks: pull to lock
 Weight: 45lbs

Weight Capacity: 250 lbs

Warranty

5 years on frame and crossbraces
 6 months on soft goods and moving parts



Model		Color / Description	Price
WC-1027	▲	Burgundy Chair With Footrests	\$1,250.00 ea.
WC-1028	▲	Burgundy Chair With Legrests	\$1,375.00 ea.
WC-5027	▲	Green Chair With Footrests	\$1,350.00 ea.
WC-5028	▲	Green Chair With Legrests	\$1,475.00 ea.

See page 8 for wheelchair parts and accessories

MRI Folding Walkers

18" Wide, Depth 17"

250 lb weight capacity

Warranty

3 years - excludes tips

Model		Description	Price
WA-1000	▲	Adult - 32" to 36"	\$195.00 ea.
WA-1001	▲	Junior - 28" to 32"	\$195.00 ea.

MRI Heavy Duty / Extra Wide Folding Walker

24" Wide, Depth 21"

500 lb weight capacity

Warranty

3 years - excludes tips

Model		Description	Price
WA-1002	▲	Extra Wide - 32" to 38½"	\$227.50 ea.

5" Adjustable Walker Wheels

Converts folding walker into wheeled walker. 8 height adjustments

350 lb weight capacity

Model		Description	Price
WA-1010	▲	5" Walker Wheels	\$49.75 pr.

MRI Transfer / Gait Belt

Model		Description	Price
GB-1000	■	2" x 48" Transfer/Gait Belt	\$18.00 ea.

■ MR Safe ▲ MR Conditional ● Not MR Safe




MRI Wheelchair Accessories



H Custom Stenciling, add your company name and/or other info on back of wheelchair

Fig	Model	Description	Price
A	RP-1053	▲ Anti-Tippers	\$44.00 pr.
B	RP-1085	▲ Elevating Legrest for 18" wide flip back arms, 20" wide, 22" wide & 24" wide wheelchair	\$225.00 pr.
B	RP-1086	▲ Elevating Legrest for 26" wide wheelchair	\$300.00 pr.
C	WC-1102	▲ Oxygen Holder	\$90.45 ea.
D	WC-1100	▲ I.V. Pole	\$101.85 ea.
E	WC-1101	▲ I.V. & Oxygen Holder	\$163.00 ea.
F	RP-1054	▲ Anti-Theft Bar for 18" Wheelchair	\$101.85 ea.
F	RP-1077	▲ Anti-Theft Bar for 20" Wheelchair	\$153.50 ea.
F	RP-1055	▲ Anti-Theft Bar for 22" Wheelchair	\$153.50 ea.
F	RP-1078	▲ Anti-Theft Bar for 24" Wheelchair	\$153.50 ea.
G	RP-1114	▲ Anti-Theft Clamps for 26" Wheelchair, set of 4	\$74.50 set
H	WC-1110	Custom Stenciling (Company Name, etc.)	\$175.00 ea.



MRI Wheelchair Parts - MR Conditional (See web site for PDF part breakdowns on all wheelchairs)

Model	Description	Cost	Model	Description	Cost
RP-1000	Footrest for 18" wide flip back arm and 20" wide chair, left	\$69.25 ea.	RP-1018	Bearing for rear wheel - 7/16" Axle	\$8.95 ea.
RP-1001	Footrest for 18" wide flip back arm and 20" wide chair, right	\$69.25 ea.	RP-1019	Bearing for rear wheel - 5/8" Axle	\$8.95 ea.
RP-1002	Footrest for 22" & 24" wide chair, Left	\$119.10 ea.	RP-1025	Brake for all MRI chairs, left	\$26.75 ea.
RP-1003	Footrest for 22" & 24" wide chair, right	\$119.10 ea.	RP-1026	Brake for all MRI chairs, right	\$26.75 ea.
RP-1004	Footrest for 18" wide fixed chair, left	\$63.50 ea.	RP-1028	Seat uph, for 18" wide chair, Burgundy	\$64.50 ea.
RP-1005	Footrest for 18" wide fixed chair, right	\$63.50 ea.	RP-1029	Seat uph, for 22" wide chair, Burgundy	\$85.00 ea.
RP-1006	Front roller bumper for footrest	\$9.75 ea.	RP-1030	Back uph, for 18" std chair, Burgundy	\$75.00 ea.
RP-1007	Heel loop and hardware for footrest	\$13.85 ea.	RP-1031	Back uph, for 18" flip back arm, Burgundy	\$77.40 ea.
RP-1009	Front fork assy. fork, wheel, bearings, Hardware	\$65.45 ea.	RP-1032	Back uph, for 22" wide chair, Burgundy	\$77.40 ea.
RP-1010	8" Wheel assy. with bearings	\$18.25 ea.	RP-1033	Seat uph, for 18" wide chair, Dark Blue	\$64.50 ea.
RP-1011	Bearing for front caster wheel 5/16"	\$7.95 ea.	RP-1034	Back uph, for 18" std chair, Dark Blue	\$75.00 ea.
RP-1012	8" front fork with nut	\$41.75 ea.	RP-1035	Padded armrest, full length, Burgundy	\$12.40 ea.
RP-1013	Fork stem bearing	\$7.95 ea.	RP-1036	Padded armrest, full length, Dk Blue	\$12.40 ea.
RP-1016	24" complete rear wheel, with 7/16" bearings	\$155.00 ea.	RP-1037	Plastic armrest, Black for flip back arm chair	\$7.95 ea.
RP-1017	24" complete rear wheel, with 5/8" bearings	\$155.00 ea.	RP-1038	Padded armrest, desk length, Burgundy	\$12.40 ea.
			RP-1046	Button lock for removable arms 22" wide	\$16.25 ea.
			RP-1115	Locking kits for foot/legrest	\$14.95 pr.
			RP-1117	Vice Tire Setter for Snap on Tires	\$29.95 ea.

■ MR Safe ▲ MR Conditional ● Not MR Safe



Fixed Gurney

NEW

Now Available with
GREEN Upholstery



ST-2000



ST-2002



ST-2001

Fixed Gurney

Stainless Steel construction
Drop down side rails
4-5" Total Lock Casters or 2-8" swivel & 2-8" locking casters
Grounding chain
Storage tray
Full surround rubber bumper
Adjustable headrest
3" pad included

Specifications:

Height	29 1/2"
Height with Rails up	38 1/2"
Height of Rails	9"
Length	77"
Rail Length	57 1/2"
Width:	30"
Storage Tray	22" x 31"
Weight:	100 lbs.
Back Adj. Range	0-75 Degrees
Number of Positions:	11

Weight Capacity: 600 lbs

Warranty

2 years on frame
90 days on soft goods and moving parts



**Optional Oxygen
Holder Available
See page 10**



Model	Description	Price
ST-2000	▲ Fixed Gurney, 5" Total Lock Casters, Black Pad	\$3,219.00 ea.
ST-2001	▲ Fixed Gurney, 2-8" Swivel & 2-8" Locking Casters, Black Pad	\$3,508.00 ea.
ST-2002	▲ Fixed Gurney, 5" Total Lock Casters, Green Pad	\$3,219.00 ea.
ST-2003	▲ Fixed Gurney, 2-8" Swivel & 2-8" Locking Casters, Green Pad	\$3,508.00 ea.



MRI Stretcher/Gurney Accessories for ST-2000, ST-2001, ST-2002, ST-2003

Model	Description	Price
PV-3000	■ Contoured Shifter	\$229.00 ea.
ST-2006	▲ I.V. Pole	\$175.00 ea.
ST-2007	■ Side Rail Pads, Burgundy	\$224.00 pr.
ST-2008	■ Safety Straps	\$82.50 set
ST-2005	■ Replacement Pad, Black	\$240.00 ea.
ST-2010	■ Replacement Pad, Green	\$240.00 ea.
ST-2009	▲ Oxygen Holder	\$179.50 ea.



Contoured Shifter

Specifications:

Width 22"
 Length 72"
 Thick 3/16"
 See page 14 for more details

Side Rail Pads

Fits the ST-2000, ST-2001, ST-2002 or ST-2003 Only

Specifications:

Width 8½"
 Length 58"

I.V. Pole

Fits the ST-2000, ST-2001, ST-2002 or ST-2003 Only

Can be attached at any of the four corners

Multiple I.V. Poles can be attached to a single gurney

Warranty: 2 years

Safety Straps

Fits the ST-2000, ST-2001, ST-2002 or ST-2003 Only

Set of 3 Straps

Replacement Pad

Specifications:

Height 3"
 Width 24½"
 Length 73½"

Oxygen Holder



Fits the ST-2000, ST-2001, ST-2002 or ST-2003 Only

Warranty: 2 years

Other Replacement Stretcher Pads

Replacement stretcher pads covered with Sure-Check Fusion III material. Available with square corners or angle cut corners. Stretcher pads are **Non-Returnable**. Custom sizes available



Square Corner Pad



Cut Corners Pad

Model	Description	Price
ST-4000	■ 24"W x 75"L x 2" T	\$209.00 ea.
ST-4001	■ 24"W x 75"L x 3" T	\$234.00 ea.
ST-4002	■ 24"W x 75"L x 4" T	\$259.00 ea.
ST-4003	■ 26"W x 75"L x 2" T	\$199.70 ea.
ST-4004	■ 26"W x 75"L x 3" T	\$234.00 ea.
ST-4005	■ 26"W x 75"L x 4" T	\$255.00 ea.
ST-4006	■ 28"W x 75"L x 2" T	\$199.70 ea.
ST-4007	■ 28"W x 75"L x 3" T	\$234.00 ea.
ST-4008	■ 28"W x 75"L x 4" T	\$255.00 ea.

Model	Description	Price
ST-4009	■ 24"W x 75"L x 2" T	\$199.70 ea.
ST-4010	■ 24"W x 75"L x 3" T	\$234.00 ea.
ST-4011	■ 24"W x 75"L x 4" T	\$255.00 ea.
ST-4012	■ 26"W x 75"L x 2" T	\$199.70 ea.
ST-4013	■ 26"W x 75"L x 3" T	\$234.00 ea.
ST-4014	■ 26"W x 75"L x 4" T	\$255.00 ea.
ST-4015	■ 28"W x 75"L x 2" T	\$199.70 ea.
ST-4016	■ 28"W x 75"L x 3" T	\$234.00 ea.
ST-4017	■ 28"W x 75"L x 4" T	\$255.00 ea.



Ferno MRI Gurney 400 lb Capacity

Model 30NM The "Standard" and the lightest lift-in cot in the industry. Ferno's Model 30 Multi-Level, Lift-in Cot served as the industry standard for many years and is still used by many ambulance companies and medical facilities. The Ferno Model 30 cot adjusts to six height levels, enabling patient transfer from most nursing home beds or hospital treatment tables. To operate the height adjustment mechanism, the EMTs must assume the patients weight prior to raising or lowering the cot-bed frame. Because the Model 30 must be lifted to be loaded into the vehicle, it can be used with ambulances that have various floor heights. The Model 30 cot will work well in station wagons, vans or modular ambulances. The telescoping head-end main frame section allows the Model 30 cot to be shortened for maneuverability around tight corners and spaces. The Model 30NM is the Non-Magnetic version for medical treatment areas where the cot must be non-magnetized.



Contoured
Style Pad

Specifications:

Max Height: 32"
Min Height: 12"
Length: 78"
Width: 23"
Weight: 56 lbs
Number of Positions: 6
Casters: All Locking
Weight Capacity: 400 lbs

Warranty

2 Years on frame
90 Days on pad and straps

Model	Description	Price
ST-1000	▲ With Contoured Pad	\$3,999.00 ea.
ST-1001	▲ With Flat Style Pad	\$3,999.00 ea.

Ferno MRI Gurney 500 lb Capacity

Model 35ANM The lightest cot available with a 500 lb weight capacity. The Ferno Model 35ANM Mobile Transporter Multi-Level Roll-In Cot combines the benefits of a multi-level cot with the easy loading feature of a roll-in style cot. Designed for a variety of ambulance floor heights, the Mobile Transporter will work well with station wagons, vans or modular ambulances. EMTs no longer have to lift from the side to load the cot. One EMT is positioned at the foot of the cot, while the other EMT simply raises the undercarriage during loading. Eight height positions enable easy patient transfer from any couch, bed or treatment table. The Mobile Transporter should be rolled in a level position. The two inclined positions are designed as loading positions for either vans or modular ambulances. The Model 35ANM is the Non-Magnetic version for medical treatment areas where the cot must be non-magnetized.



Contoured
Style Pad

Specifications:

Max Height: 40"
Min Height: 9"
Max Length: 63" -79"
Width: 24"
Weight: 74 lbs
Number of Positions: 8
Casters: All Locking
Weight Capacity: 500 lbs

Warranty

2 Years on frame
90 Days on pad and straps

Model	Description	Price
ST-1015	▲ With Contoured Pad	\$3,950.00 ea.
ST-1017	▲ With Flat Style Pad	\$3,950.00 ea.



Ferno Gurney Parts & Accessories

Model	Description	Price
ST-1004	▲ IV Pole	\$354.60 ea.
ST-1009	■ Safety Belts - Set of 3	\$305.00 ea.
ST-1007	▲ Locking Caster	\$220.00 ea.
ST-1008	▲ Non-Locking Caster 400lb Capacity Gurney	\$163.00 ea.
ST-1020	▲ Non-Locking Caster 500lb Capacity Gurney	\$284.00 ea.
ST-1005	▲ Contoured Pad	\$246.30 ea.
ST-1006	▲ Flat Style Pad	\$222.50 ea.



I.V. Pole
For Ferno Gurneys
Only

**Replacement
Safety Belts**



**Replacement
Casters**



Replacement Pads
Flat Style Pad



Shower Gurney with 3 Position Elevating Headrest

Specifications:

79" Long x 28" Wide x 33-1/2" High (rails down) 43"
High (rails up)

350 lb weight capacity

Heavy Duty Specifications:

79" Long x 32" Wide x 33-1/2" High (rails down) 43"
High (rails up)

600 lb weight capacity



Custom sizes available upon request

Warranty

1 Year on frame

6 Months on casters and soft goods

Model	Description	Price
PV-1096	▲ With 2" Vinyl Foam Pad Cushion (pictured)	\$1,104.98 ea.
PV-1098	▲ With Mesh Cushion	\$1,115.79 ea.
PV-1129	▲ Heavy Duty With 2" Vinyl Foam Pad Cushion (pictured)	\$1,209.01 ea.
PV-1180	▲ Heavy Duty With 2" Closed Cell Water Proof Pad	\$1,270.00 ea.
PV-1181	■ Replacement 2" Vinyl Foam Pad Cushion (pictured on stretcher)	\$256.50 ea.
PV-1169	■ Replacement 2" Closed Cell Water Proof Pad (pictured at right)	\$266.35 ea.

Closed Cell Water Proof Pad



More options and sizes of PVC Stretchers / Gurneys available online at www.MRlequip.com



Patient Wheelchair Transfer Board

30" Long x 12.5" Wide

The transfer board bridges the gap between two transfer points. The board's anti-slip tape enables it to remain immobile during transfer procedures. The upper side of the board is slick to allow the patient to slide easily.

Flexible plastic.

The Board is Radiotranslucent.

2 Year Warranty



Model	Description	Price
TR-1010	■ Transfer Board	\$97.00 ea.



Shockless Shifter

Won't build up a static charge, radiolucent, can't damage electronics, 100% shock protection. Surprise zaps aren't just painful and dangerous for patients, they can jeopardize expensive electronic medical equipment and cause power outages.

Shifter is made with carbon-filled plastic that actually conducts electricity, not simply treated with an antistatic product to reduce static. Static charges can't build up on a conductive surface, so there are no more painful shocks ever, GUARANTEED!

350lb Weight Limit

Model	Description	Price
TR-1000	■ Shifter 22" wide x 72" Long, Black	\$375.00 ea.
Disc.	■ Shifter 22" wide x 72" Long, White	Disc.
TR-1003	■ Shifter 18" wide x 72" Long, Aqua	\$262.00 ea.

PVC Body Transfer Board

Specifications:

71" Long x 23-1/2" Wide x 1/4" Thick

Made of high density polyethylene

Weight capacity: 400 lb

Warranty: 6 Months

PVC Body Transfer Board Storage Brackets

Set of 3 brackets that allow for vertical or horizontal storage on a wall.



Model	Description	Price
PV-1121	■ PVC Body Transfer Board	\$171.39 ea.
TR-1009	■ PVC Transfer Board Storage Brackets	\$101.80 set





Contoured Shifter - Patient Mover

Slides easily on a bed sheet or stretcher pad thanks to the nearly friction-free surface. Special slip coating prevents the patient from sticking to the mover during transfers and results in less static buildup. No more wrenching pulls; no more back strains. Uneven transfers are a breeze, the semi rigid board bridges uneven gaps. Radiolucent and artifact free, the shifter can be left in place during X-ray procedures and minor surgery.

Features

New contoured shape for easier patient loading and unloading, especially in tight quarters

Built in ionic silver antimicrobial agents to control levels of bacteria - active for the life of the board

Antistatic inhibitors (static **reducing** properties)

Anti-pinch hand holds

Ergonomic design to reduce healthcare workers stress and injuries

X-Ray translucent

Artifact free

Specifications:

Length 72"
Thick 3/16"
Width 22"

500 lb Weight capacity, can move 300 lb

Please call for certificate of performance for MRSA / s.aureus



MT-1651



Wall Hanger Hooks

NOT MR SAFE

This product is magnetic

Model	Description	Price
PV-3000	■ 22" Wide Patient Mover	\$229.00 ea.
MT-1651	■ Replacement Jumbo Hooks, 2 ea	\$18.95 set
TR-1004	● Vertical Wall Hanger Hooks	\$101.80 ea.
TR-1005	● Horizontal Wall Hanger Hooks	\$159.35 ea.

Bariatric Patient Shifter

Bariatric Shifter makes lateral transfers for patients up to 600 lbs. easier and safer. Reduce caretaker back strain by minimizing friction on transfer services. Elongated shape brings fourteen handholds to ergonomically correct position. Extra-wide design better supports excess tissue. It even supports patient's head. Flexible plastic is great for uneven surfaces.

Recommend using 1 caretaker per 100 lbs. 28"W x 72"L x 3/16"H.

Weight capacity: 600 lbs. Purple.

Must ship by truck

Features:

Helps reduce back strain

Ergonomic design

Supports up to 600 lbs.

Extra wide for better support of excess tissue



Model	Description	Price
TR-1600	■ Bariatric Patient Shifter	\$358.90 ea.

■ MR Safe ▲ MR Conditional ● Not MR Safe



Load Grips, w/Rubberized Pre-Molded Handle

Features:

Use to lift backboards and other heavy and awkward gear.
 Compatible with all backboards.
 Quick and easy application.
 Durable and long lasting with normal cleaning protocol.
 2" poly strap with rubberized pre-molded handle for a no-slip grip.

Model	Description	Price
LG-1000-CL	Load Grips - set of 4, Red or Black	\$18.00 set

Safe Sit Manual Handling Aid for CT, MRI and PET



The Safe Sit has been designed to be used after all the procedures and imaging has occurred. It is a simple and easy to apply sheet which is placed onto the bottom of the table for assisting the patient to sit up.

The product is a disposable sheet - made from a laminated paper and as such would pose no issues for magnetic procedures.

Safe Sit™ is a single-patient use manual handling aid designed to be used by the patient several times before disposal thereby helping to minimise crossinfection within the healthcare environment.

Safe Sit™ has been specifically designed to be quick and easy to use providing patients with additional help when sitting up after a clinical procedure. Safe Sit™ minimizes potential injury to the healthcare professional assisting them. The use of Safe Sit™ may help to safely improve efficiency within radiotherapy and imaging

departments. In addition to use in radiotherapy suites Safe Sit™ can be used on CT, MRI and PET scanners, trolleys, stretchers, beds, examination couches and even operating tables after day case surgery or investigations.

Model	Description	Price
SS-4000	Safe Sit	\$54.00 ea.



MRI Non-Magnetic PVC Crib

Specifications

Height: 51" Width: 31"
 Length: 54" Bed to Top: 29"
 Rails slide down to mattress level
 2.5" space between rails
 6 qty 3" twin plastic casters with non magnetic stainless steel stems
 Mattress not included, but comes with an Antimicrobial 3" padded cushion, Weight Capacity: 300 lbs



Model	Description	Price
PV-1190	White PVC Crib	Call
PV-1191	Yellow PVC Crib	Call





Tables/Carts

Constructed of Stainless Steel, 2" Dual Wheel Casters
 Without rails 34" high, With rails 36¾" high
 Assembly Required, Can be assembled for additional \$50.00

18" x 24"

Model	Description	Price
TA-5001	▲ With top shelf and rails	\$483.30 ea.
TA-5007	▲ With top shelf, no rails	\$439.50 ea.
TA-5013	▲ With two shelves and rails	\$651.95 ea.
TA-5019	▲ With two shelves, no rails	\$609.75 ea.

18" x 30"

Model	Description	Price
TA-5002	▲ With top shelf and rails	\$519.50 ea.
TA-5008	▲ With top shelf, no rails	\$471.88 ea.
TA-5014	▲ With two shelves and rails	\$715.40 ea.
TA-5020	▲ With two shelves, no rails	\$670.65 ea.

18" x 36"

Model	Description	Price
TA-5003	▲ With top shelf and rails	\$559.50 ea.
TA-5009	▲ With top shelf, no rails	\$507.65 ea.
TA-5015	▲ With two shelves and rails	\$780.95 ea.
TA-5021	▲ With two shelves, no rails	\$731.35 ea.

24" x 36"

Model	Description	Price
TA-5004	▲ With top shelf and rails	\$621.40 ea.
TA-5010	▲ With top shelf, no rails	\$562.55 ea.
TA-5016	▲ With two shelves and rails	\$890.50 ea.
TA-5022	▲ With two shelves, no rails	\$831.15 ea.

24" x 48"

Model	Description	Price
TA-5005	▲ With top shelf and rails	\$697.60 ea.
TA-5011	▲ With top shelf, no rails	\$628.60 ea.
TA-5017	▲ With two shelves and rails	\$1,025.45 ea.
TA-5023	▲ With two shelves, no rails	\$953.55 ea.

24" x 72" Must Ship via Truck

Model	Description	Price
TA-5006	▲ With top shelf and rails	\$849.35 ea.
TA-5012	▲ With top shelf, no rails	\$765.60 ea.
TA-5018	▲ With two shelves and rails	\$1,294.45 ea.
TA-5024	▲ With two shelves, no rails	\$1,207.75 ea.





MRI Step Stools

Specifications:

16" x 12" Base
 8 1/2" Step Height
 Made of Heavy Duty Stainless Steel
500 lb Weight Capacity

Model	Description	Price
SS-1000	Without 41" Handrail	\$177.45 ea.
SS-1001	With 41" Handrail	\$219.90 ea.



MRI Bariatric Step Stool

Specifications:

30" x 12" Base
 8 1/2" Step Height
 Made of Heavy Duty Stainless Steel
500 lb Weight Capacity



Model	Description	Price
SS-1733	Without 41" Handrail	\$250.00 ea.
SS-1734	With Single 41" Handrail	\$342.50 ea.
SS-1005	With Dual 41" Handrail	\$374.00 ea.



MRI Double Step Stool With 41" Handrail

Specifications:

36" L X 12" W X 8" H Top Step
 36" L X 24" W X 8" H Base
 Weight: 47 lbs
500 lb Weight Capacity

Model	Description	Price
SS-1002	With Dual 41" Handrail	\$715.00 ea.

More options and sizes of Step Stools available online at www.MRIequip.com, including PVC options with casters and higher weight capacities.



MRI Narrow Double Step Stool

Specifications:

16" L X 12" W X 8" H Top Step

16" L X 12" W X 8" H Base

Weight: 20 lbs

500 lb Weight Capacity



Model	Description	Price
SS-1007	▲ Narrow Step Stool with out Handle	\$449.00 ea.



MRI Narrow Double Step Stool with Handle

Specifications:

16" L X 12" W X 8" H Top Step

16" L X 12" W X 8" H Base

Handrail: 47" Tall x 26" Deep

Weight: 20 lbs

500 lb Weight Capacity



Model	Description	Price
SS-1008	▲ Narrow Step Stool with Handle	\$499.00 ea.



MRI Stools

Stainless steel stools with solid non-slip material mounted into recess of stool seat. Seat diameter is 14". Comes with rubber tips or casters. Choose from a fixed height stackable stool or heights that adjust with threaded stem, from 15" to 21" or 21" to 27".



CH-8058



CH-8052



CH-8056



CH-8057

Model	Description	Price
CH-8058	▲ Fixed Height Stool w/ Rubber Tips, 18"	\$276.75 ea.
CH-8050	▲ Adjustable Stool w/ Rubber Tips, 15"-21"	\$405.00 ea.
CH-8054	▲ Adjustable Stool w/ 2" Casters, 15"-21"	\$433.00 ea.
CH-8051	▲ Adjustable Stool w/ Rubber Tips, 15"-21" & Back	\$545.00 ea.
CH-8055	▲ Adjustable Stool w/ 2" Casters, 15"-21" & Back	\$573.00 ea.
CH-8052	▲ Adjustable Stool w/ Rubber Tips, 21"-27"	\$410.70 ea.
CH-8056	▲ Adjustable Stool w/ 2" Casters, 21"-27"	\$459.26 ea.
CH-8053	▲ Adjustable Stool w/ Rubber Tips, 21"-27" & Back	\$555.00 ea.
CH-8057	▲ Adjustable Stool w/ 2" Casters, 21"-27" & Back	\$580.00 ea.



CH-8030

Stools

4" Padded Seat
Burgundy Naugahyde Upholstery
Stainless Steel Construction
5 Leg Design for Maximum Stability
Padded Armrests (on selected models)
2" Dual Wheel Casters (on selected models)
250 lb Weight Capacity

Warranty

2 years on frame
6 months on soft goods and moving parts



Adjusts 15" to 21" with Rubber Tips

Burgundy	Green	Description	Price
CH-8030	CH-9030	▲ Stool w/ Rubber Tips	\$388.00 ea.
CH-8031	CH-9031	▲ Stool w/ Rubber Tips & Back	\$458.00 ea.
CH-8032	CH-9032	▲ Stool w/ Rubber Tips, Back & Arms	\$498.00 ea.



CH-8032

Adjusts 21" to 27" with Rubber Tips

Burgundy	Green	Description	Price
CH-8033	CH-9033	▲ Stool w/ Rubber Tips	\$408.00 ea.
CH-8034	CH-9034	▲ Stool w/ Rubber Tips & Back	\$478.00 ea.
CH-8035	CH-9035	▲ Stool w/ Rubber Tips, Back & Arms	\$518.00 ea.

Adjusts 16" to 22" with Casters

Burgundy	Green	Description	Price
CH-8036	CH-9036	▲ Stool w/ Rubber Tips	\$438.00 ea.
CH-8037	CH-9037	▲ Stool w/ Rubber Tips & Back	\$508.00 ea.
CH-8038	CH-9038	▲ Stool w/ Rubber Tips, Back & Arms	\$549.00 ea.



CH-8038

Adjusts 22" to 28" with Casters

Burgundy	Green	Description	Price
CH-8039	CH-9039	▲ Stool w/ Rubber Tips	\$444.00 ea.
CH-8040	CH-9040	▲ Stool w/ Rubber Tips & Back	\$513.00 ea.
CH-8041	CH-9041	▲ Stool w/ Rubber Tips, Back & Arms	\$575.00 ea.



Custom Stenciling Option

Customize your doctor stool. Add your facility name, department name or any information you would like.

Model	Description	Price
CH-8914	Custom Stenciling Option	\$75.00 ea.



MRI Furniture

The MRI furniture line is built completely with non-magnetic supplies. The furniture is made to go directly into the MR environment. Available now, Chair, Bariatric Chair, Sofa, End Table, and Coffee Table. Custom items may be available upon request. Choose from 4 fabric choices or you may supply your own fabric at NO extra charge (would need 3 yards of fabric for love seat and 1½ yards of fabric for chair).

Chair




Specifications:

24" W x 22 1/2" D
Seat Height: 20"
Arm Height: 28"

Bariatric Chair Specifications:

26" W x 22 1/2" D
Seat Height: 20"
Arm Height: 28"
500 lb Weight Capacity

Sofa



Specifications:

53" W x 24 1/2" D
Seat Height: 20"
Arm Height: 28"



Coffee Table

Specifications:

24" W x 48" L
18" High



Pediatric Chair (no arms)

NEW

Specifications:

16" W x 14" D
Seat Height: 10"
NO Arms
Burgundy Nagahyde
Weight Capacity: 250 lbs



End Table

Specifications:

22 1/2" W x 22 1/2" D
24" High



Pediatric Table Solid Oak

NEW

Specifications:

21" L x 25 1/2" W x 18" H




Purchase 1 Pediatric Table and 2 Pediatric Chairs SAVE \$50.00

Model	Description	Price
FR-1005	▲ Pediatric Chair, Burgundy Nagahyde	\$225.00 ea.
FR-1006	▲ Pediatric Table	\$280.00 ea.
FR-1000	▲ Sofa, 350 lb Weight Capacity	\$1,251.00 ea.
FR-1001	▲ Chair, 250 lb Weight Capacity	\$660.45 ea.
FR-1004	▲ Bariatric Chair, 500 lb Weight Capacity	\$835.45 ea.
FR-1002	▲ End Table	\$441.20 ea.
FR-1003	▲ Coffee Table	\$626.20 ea.

You may supply your own fabric at NO extra charge, except Pediatric Chair

■ MR Safe ▲ MR Conditional ● Not MR Safe

Fabric Choices





All Carts Feature:

- Aluminum
- Replaceable Plastic Top
- Stainless Steel Pull-Out Shelf
- Ball Bearing Full Extension Drawer Slides

Dimensions:

45.75" High
30.50" Wide
23.75" Deep
126 Lbs.



Dimensions:

45.75" High
32" Wide
25" Deep
136 Lbs.



Dimensions:

45.75" High
32" Wide
25" Deep
136 Lbs.



Breakaway Lock Carts

6 Drawers

One 4", Two 3.25", Two 6.5" & One 9.75" Drawers

7 Drawers

One 4", Four 3.25", One 6.5" & One 9.75" Drawers

Warranty

5 years on carts
1 year on locks & moving parts

Model	Description	Price
CA-1000	▲ 6 Drawer Breakaway Lock Cart	\$2,460.70 ea.
CA-1002	▲ 7 Drawer Breakaway Lock Cart	\$2,704.25 ea.

Keyed Lock Carts

6 Drawers

One 4", Two 3.25", Two 6.5" & One 9.75" Drawers

7 Drawers

One 4", Four 3.25", One 6.5" & One 9.75" Drawers

Two Keys come with Keyed Lock carts Extra Keys Available (see page 23)

Warranty

5 years on carts
1 year on locks & moving parts

Model	Description	Price
CA-1004	▲ 6 Drawer Keyed Lock Cart	\$2,909.75 ea.
CA-1006	▲ 7 Drawer Keyed Lock Cart	\$3,181.41 ea.

Keyed Lock Carts with Rails and Organizers

Includes:

- 3 Tier raised back rail system
- 4 Side rails
- 2 Full drawer trays with organizers

Warranty

5 years on carts
1 year on locks & moving parts

6 Drawers

One 4", Two 3.25", Two 6.5" & One 9.75" Drawers

7 Drawers

One 4", Four 3.25", One 6.5" & One 9.75" Drawers

Two Keys come with Keyed Lock carts Extra Keys Available (see page 23)

Model	Description	Price
CA-1005	▲ 6 Drawer Keyed Lock Cart	\$3,459.41 ea.
CA-1007	▲ 7 Drawer Keyed Lock Cart	\$3,737.00 ea.

Breakaway Lock Carts with Emergency Package

Includes:

I.V. Pole, Cardiac board, Oxygen tank holder, 4 utility hooks, 2 Drawer trays

6 Drawers

One 4", Two 3.25", Two 6.5" & One 9.75" Drawers

7 Drawers

One 4", Four 3.25", One 6.5" & One 9.75" Drawers

Warranty

5 years on carts
1 year on locks & moving parts

Model	Description	Price
CA-1001	▲ 6 Drawer Breakaway Lock Cart	\$2,913.87 ea.
CA-1003	▲ 7 Drawer Breakaway Lock Cart	\$3,000.00 ea.

Visit www.MRlequip.com or call customer service for more carts and options

All Carts Ship Red Unless Color Specified from page 23

■ MR Safe ▲ MR Conditional ● Not MR Safe



Narrow Design Breakaway Lock Carts

6 Drawers

One 4", Two 3.25", Two 6.5" & One 9.75" Drawers

7 Drawers

One 4", Four 3.25", One 6.5" & One 9.75" Drawers

Warranty

5 years on carts

1 year on locks & moving parts

Model	Description	Price
CA-1050	▲ 6 Drawer Breakaway Lock Cart	\$2,312.35 ea.
CA-1052	▲ 7 Drawer Breakaway Lock Cart	\$2,463.42 ea.

All Carts Feature:

Aluminum, Replaceable
Plastic Top, Stainless
Steel Pull-Out Shelf,
Ball Bearing Full
Extension Drawer Slides

Narrow Design Breakaway Lock Emergency Carts

Includes:

I.V. Pole, Cardiac board, Oxygen tank holder, 4 utility hooks, 2 Drawer trays

6 Drawers

One 4", Two 3.25", Two 6.5" & One 9.75" Drawers

7 Drawers

One 4", Four 3.25", One 6.5" & One 9.75" Drawers

Warranty

5 years on carts

1 year on locks & moving parts

Model	Description	Price
CA-1051	▲ 6 Drawer Breakaway Lock Cart	\$2,703.75 ea.
CA-1053	▲ 7 Drawer Breakaway Lock Cart	\$2,866.82 ea.

Narrow Design Keyed Lock Carts

6 Drawers

One 4", Two 3.25", Two 6.5" & One 9.75" Drawers

7 Drawers

One 4", Four 3.25", One 6.5" & One 9.75" Drawers

Two Keys come with Keyed Lock carts Extra Keys Available (see page 23)

Warranty

5 years on carts

1 year on locks & moving parts

Model	Description	Price
CA-1054	▲ 6 Drawer Keyed Lock Cart	\$2,708.04 ea.
CA-1056	▲ 7 Drawer Keyed Lock Cart	\$2,981.85 ea.

Narrow Design Keyed Lock Anesthesia Carts

6 Drawers

One 4", Two 3.25", Two 6.5" & One 9.75" Drawers

7 Drawers

One 4", Four 3.25", One 6.5" & One 9.75" Drawers

Includes:

3 Tier raised back rail system

4 Side rails

2 Full drawer trays with organizers

Two Keys come with Keyed Lock carts Extra Keys Available (see page 23)

Warranty

5 years on carts

1 year on locks & moving parts

Model	Description	Price
CA-1055	▲ 6 Drawer Keyed Lock Cart	\$3,239.35 ea.
CA-1057	▲ 7 Drawer Keyed Lock Cart	\$3,496.85 ea.



All Carts Ship Red Unless Color Specified from page 23



Optional Cart Accessories



CA-1026



CA-1020



CA-1033



CA-1013

CA-1022



CA-1010

CA-1012



CA-1017



CA-1018

CA-1019

Available Colors



Hammertone Colors:



Model	Description	Price
CA-1008	▲ 3 Tier Raised Back Rail	\$281.95 ea.
CA-1009	▲ 2 Tier Raised Back Rail	\$247.97 ea.
CA-1010	▲ Side Rail	\$81.06 ea.
CA-1011	▲ Cardiac Board & Mounting Brackets	\$163.15 ea.
CA-1012	▲ Rail Hook	\$21.52 ea.
CA-1028	▲ Key for Keyed Carts	\$19.10 ea.
CA-1029	▲ Breakaway Key for Carts (package of 50)	\$46.29 pk.
CA-1013	▲ I.V. Pole	\$136.37 ea.
CA-1033	▲ Oxygen Tank Holder	\$84.40 ea.
CA-1015	▲ Replacement Pull-Out Shelf	\$189.10 ea.
CA-1016	▲ Set of 4 Utility Hooks	\$48.92 ea.
CA-1017	▲ 4 Tilt Bin Organizer for Raised Back Rail System	\$152.44 ea.
CA-1018	▲ 5 Tilt Bin Organizer for Raised Back Rail System	\$118.34 ea.
CA-1019	▲ 6 Tilt Bin Organizer for Raised Back Rail System	\$89.61 ea.
CA-1020	▲ 3 Gallon Waste Basket	\$79.41 ea.
CA-1022	▲ 3" Drawer Organizer (does half of shallow drawer)	\$64.07 ea.
CA-1023	▲ 6" Drawer Organizer (does half of 6.5" drawer)	\$79.16 ea.
CA-1026	▲ 3 Gallon Sharps Container	\$79.41 ea.



Stainless Steel Hampers

All Stainless Hampers Do Not Include a bag

Specifications:

Height: 36"
Inside Diameter: 16 1/2"
Outside Diameter: 18"
Casters: 2" Dual Wheel

Bag is NOT included, it must be ordered separately.

Model		Description	Price
HA-2000	▲	Hamper - No Lid	\$278.00 ea.
HA-2001	▲	Hamper with Lid	\$378.00 ea.
HA-2002	▲	Hamper with Foot Operated Lid	\$459.00 ea.
HA-2003	■	Laundry Bag	\$38.50 ea.



PVC Hampers

Warranty

1 year on frame
6 months on casters and soft goods

All mesh hampers sent black unless color specified.

All leak proof hampers sent white unless color specified. No zipper on Leak Proof Bags

Model		# of Bins	Bag Type	Length	Width	Height	Price
PV-1000	▲	1	28 Gallon Mesh	18"	18"	37"	\$334.65 ea.
PV-1001	▲	1	22 Gallon Leak Proof	18"	18"	37"	\$375.69 ea.
PV-1007	▲	1	42 Gallon Mesh	26 1/2"	22 1/2"	35 1/2"	\$349.73 ea.
PV-1002	▲	1	28 Gallon Mesh	23 3/4"	18 3/4"	38"	\$445.52 ea.
PV-1003	▲	2	22 Gallon Leak Proof	38 1/4"	18"	37"	\$520.09 ea.
PV-1004	▲	3	28 Gallon Mesh	58 1/4"	18"	37"	\$547.60 ea.
PV-1005	▲	3	22 Gallon Leak Proof	58 1/4"	22 1/2"	35 1/2"	\$772.86 ea.



PV-1001
Casters are different than pictured.

Available Colors



Leakproof Colors



PV-1006



PV-1003
Casters are different than pictured.



PV-1005

Foot Pedal for Hampers

Foot pedal available for hampers, order one per bin, must be ordered with hamper.

Model		Description	Price
PV-1006	▲	Foot Pedal for Hampers	\$24.72 ea.

Visit www.MRIequip.com or call customer service for more laundry baskets and linen carts

Casters are different than pictured.

■ MR Safe ▲ MR Conditional ● Not MR Safe



PV-1042

Casters are different than pictured.



PV-1038

Casters are different than pictured.



PV-1034

Casters are different than pictured.

PVC Linen / Multi-Use Carts

Features:

MR Compatible

Totally enclosed carts help prevent dust build-up on expensive coils.

All carts come with solid vinyl fabric, mesh can be substituted.

Carts are shipped with black vinyl, unless color specified.

Velcro closures

Caster sockets are double wall reinforced with schedule 40 one inch PVC

Can upgrade to 5" casters for additional charge

Warranty

1 year on frame

6 months on casters and soft goods

All multi-use carts are shown with the front cover open.

Available Colors



Yellow



Tan



Mauve



Hunter Green



Red



Navy



Black



Gray



Light Blue



White



Royal Blue



Mint

Model		# of Shelves	Shelf Size	Height Between Shelves	Height with Casters	Weight Capacity per Shelf	Price
PV-1031	▲	2	25" x 20"	12"	37"	75 lb	\$402.36 ea.
PV-1033	▲	3	25" x 20"	12"	51 ¼"	75 lb	\$451.39 ea.
PV-1034	▲	4	25" x 20"	12"	65 ½"	75 lb	\$535.80 ea.
PV-1036	▲	2	32" x 20"	12"	37"	100 lb	\$423.94 ea.
PV-1038	▲	3	32" x 20"	12"	51 ¼"	100 lb	\$516.18 ea.
PV-1039	▲	4	32" x 20"	12"	65 ½"	100 lb	\$592.76 ea.
PV-1042	▲	2	45" x 20"	16"	46"	125 lb	\$592.76 ea.
PV-1043	▲	3	45" x 20"	16"	65"	125 lb	\$705.55 ea.
PV-1044	▲	4	45" x 20"	16"	77"	125 lb	\$775.28 ea.
PV-1045	▲	2	50" x 20"	16"	45 ¼"	150 lb	\$633.96 ea.
PV-1046	▲	3	50" x 20"	16"	65"	150 lb	\$749.73 ea.
PV-1047	▲	4	50" x 20"	16"	77"	150 lb	\$885.18 ea.



MR Safe



MR Conditional



Not MR Safe



3 Drawer Universal Cart

Length: 25"
 Depth: 16"
 Height with Casters: 48"
 Shelf Size: 25" x 16"
 Drawers: 3 Slide Out
 Drawer Size: 19.125" W x 14" D x 6.5" H
 Can upgrade to 5" casters for additional charge

Warranty

1 year on frame
 6 months on casters and soft goods

Model	Description	Price
PV-1149 ▲	3 Drawer Cart w/ Shelf	\$508.30 ea.

6 Drawer Universal Cart

Length: 27"
 Depth: 16"
 Height with Casters: 57"
 Drawers: 6 Slide Out
 Drawer Size: 19.125" W x 14" D x 6.5" H
 Can upgrade to 5" casters for additional charge

Model	Description	Price
PV-1150 ▲	6 Drawer Cart w/o Shelf	\$632.00 ea.

10 Drawer Universal Cart

Length: 28"
 Depth: 14"
 Height with Casters: 49"
 Drawers: 10 Slide Out
 Drawer Size: 10" W x 13.5" D x 3" H
 (Top Eight Drawers)
 10" W x 13.5" D x 7" H
 (Bottom Two Drawers)
 Can upgrade to 5" casters for additional charge

Warranty

1 year on frame
 6 months on casters and soft goods

Model	Description	Price
PV-1151 ▲	10 Drawer Cart w/ Shelf	\$669.29 ea.

10 Drawer Locking Universal Cart

Locking Enclosure
 Length: 28" Depth: 14"
 Height with Casters: 49"
 Drawers: 10 Slide Out
 Drawer Size: 10" W x 13.5" D x 3" H (Top Eight Drawers)
 10" W x 13.5" D x 7" H (Bottom Two Drawers)
 Cart is lockable (Pad-Lock is not included -PL-1000 or PL-1001)
 Not approved for medications or treatments
 Can upgrade to 5" casters for additional charge

Model	Description	Price
PV-1152 ▲	10 Drawer Locking Cart	\$896.00 ea.





Crash Cart with Bins Only
Oxygen Tank & Suction
Pump NOT included.

PVC Crash Cart and Optional Accessories

Specifications:

Length: 31-1/2"
Width: 20"
Height: 34.5"
Comes with 3" stainless steel casters

Warranty

1 year on frame
6 months on casters and soft goods

Accessories shipped in black vinyl unless color specified.

Model	Description	Price
PV-1103	▲ Crash Cart	\$473.80 ea.
PV-1106	▲ CPR Board, 20" x 20"	\$100.57 ea.
PV-1107	▲ Bag for CPR Board	\$42.12 ea.
PV-1104	▲ Hood for Top Section	\$98.46 ea.
PV-1105	▲ Side Panel Covers Choose Mesh or Solid vinyl	\$70.40 ea.

Custom sizes available upon request

Available Colors



PV-1106 & PV-1107



PV-1105



PV-1104

PVC Utility Carts

Specifications:

3 shelves
Length with handle: 31-1/2"
Width between shelves: 12"
Comes with 3" stainless steel casters
Height: 36"
Shelf Weight Capacity: 75 lb

Warranty

1 year on frame
6 months on casters and soft goods

Model	Description	Price
PV-1030	▲ 20" x 25" Shelf Size	\$359.05 ea.
PV-1172	▲ 24" x 25" Shelf Size	\$383.50 ea.

Custom sizes available upon request



Casters are different than
pictured.



MRI Non-Magnetic 3 Shelf Tub Cart

24"W X 18" D X 38 1/2" H

Tubs are 2 3/4" Deep. Shelves and legs are made from high density plastic that will not stain, scratch, rust or dent. Push handle molded into the top shelf/tub. All shelves are reinforced with 2 aluminum bars.

Weights 29 lbs

NEW



Model	Description	Price
PV-4004	▲ Blue, 3 Shelf	\$223.85 ea.
PV-4005	▲ Black, 3 Shelf	\$223.85 ea.
PV-4006	▲ Gray, 3 Shelf	\$223.85 ea.
PV-4007	▲ Putty, 3 Shelf	\$223.85 ea.

Non-Magnetic Utility Carts

Ideal mobility solution for medical equipment that is operated near MRI machines. The universal logos displayed on the wheels clearly show the medical facility that these items are approved for usage around MRI machines. Also for non-corrosive applications where water and moisture can present potential issues.

- Holds 100 lbs. per shelf.
- MRI room compatible.
- Three shelf unit will hold 300lbs. evenly distributed



2 Shelf Cart, Choose Color

Specifications: 18"W x 24"D x 24 1/2"H, 22lbs.
Space between Shelves is 15.5"

Model	Description	Price
PV-4000	▲ 2 Shelf Cart	\$248.00 ea.

3 Shelf Cart, Choose Color

Specifications: 18"W x 24"D x 42"H, 29 lbs.
Space between shelves is 15.5"

Model	Description	Price
PV-4003	▲ 3 Shelf Cart	\$272.00 ea.

3 Shelf Cart, Choose Color

Specifications: 18"W x 24"D x 34"H, 28 lbs.
Space between Shelves is 12" from bottom to middle shelf and 12" from middle to top.

Model	Description	Price
PV-4002	▲ 3 Shelf Cart	\$264.00 ea.

3 Shelf Cart, Choose Color

Specifications: 18"W x 24"D x 34"H, 28lbs.
Space between Shelves is 16" from bottom to middle shelf and 8" from middle to top.

Model	Description	Price
PV-4001	▲ 3 Shelf Cart	\$259.00 ea.





Mobile Light

39" heavy duty stainless steel dual joint arm, 50 watt head. Mounted on a 25" diameter base with five 2" casters, 14 lb. stainless steel weight added for stability and a 39" tall pole. Light has a 30 foot cord to transformer box (see warning below) 9.5 foot cord from transformer to wall outlet.

Warranty: 3 years

Model	Description	Price
EX-1000	▲ Mobile Light with Power Box	\$2,800.00 ea.
EX-1013	▲ Mobile Light w/o Power Box	\$1,655.25 ea.
EX-1004	▲ Power Box	\$1,144.75 ea.
EX-1011	▲ Replacement Bulb for EX-1000	\$61.00 ea.

Warning: The installation guidelines must be strictly followed. The electrical box contains ferrous materials and must **NEVER** be placed close to the magnet. If the electrical box is moved closer to the magnet than the 100 gauss line, electromagnetic interference may be seen, to avoid this **NEVER** move the electrical box close to the magnet.

- **NEVER** operate the light while the electrical box is unsecured.
- The lamp may slightly vibrate



More Lights

Mobile and Dual Versions at
www.mriequip.com



Ceiling Mount and Mobile Lights

29" extension arm, 720° rotation & 69" reach to beam center, 359° lamp head/ balance arm rotation in horizontal and vertical planes, Unlimited beam placement and Balanced arm rotation. Light has a 30 foot cord to transformer box (see warning above) 9.5 foot cord from transformer to wall outlet.

Model	Description	Price
EX-1002	▲ Single Pod	\$4,700.00 ea.
EX-1005	▲ Dual Pod (not pictured)	\$9,500.00 ea.
EX-1006	▲ Mobile Single Pod	\$4,500.00 ea.
EX-1008	▲ 3 Pack of Replacement Bulbs	\$183.15 pk.
EX-1012	▲ 35 Watt Replacement Bulb	\$72.03 ea.

■ MR Safe ▲ MR Conditional ● Not MR Safe



SCHILLER

A M E R I C A



MAGLIFE Lite
Basic Vital Signs



MAGSCREEN
Remote Control



MAGLIFE Serenity
Patient Monitor

**Special Limited
Time Offer!**

Contract your purchase for your MAGLIFE Serenity
MRI Patient Monitor before July 31st, 2011, and receive a
Free 2-Channel Invasive Blood Pressure Upgrade!
A \$5,000 Value! (Upgrade # 1-094-2380)

MAGLIFE Serenity - Perfection in MRI Monitoring

The new MAGLIFE Serenity guarantees highest ECG quality during magnetic resonance imaging (MRI) scanning - even under strongest gradient influence. It monitors all vital parameters during anesthesia in an MRI environment and is specifically developed for adults, children and neonates.

The Advantages of MAGELIFE at a glance:

- Highest ECG quality even under strongest gradient influence (MRI)
- Several gating outputs for optimal MRI triggering
- Core and skin temperature measurement by use of an optical sensor
- State-of-the-art gas analysis during anesthesia
- Network compatibility (wireless LAN, fiber cable LAN)
- Save penetration of the Faraday cage
- Accessories suitable for MRI
- Flexibility in the configuration of options



**Call Schiller America Today for a
Quotation or to schedule your
demonstration!**

Toll Free: 1 (888) 845-8455


or e-mail:

sales@schilleramerica.com



Stethoscope

Model	Description	Price
BP-1000	Stethoscope	\$299.00 ea.



Disposable Stethoscopes

Stethoscopes feature a plastic binaural and ultra sensitive plastic chestpiece with remarkable acoustical performance. Designed for single-patient use to prevent cross contamination in infectious areas. Case includes 10 each.

Model	Description	Price
BP-1012	■ Yellow Stethoscope	\$33.00 cs.
BP-1013	■ Blue Stethoscope	\$33.00 cs.



Radiotranslucent Electrodes

Perforated clear tape for comfort
Radiotranslucent for X-Ray/MRI
Clear tape does not absorb fluids
Gel pocket stabilizes the electrode
Clear tape for skin monitoring
Case includes 30 ea.

Model	Description	Price
BP-1014	■ Electrodes	\$60.00 cs.



Accusensor 38mm Radiotranslucent Foam Electrodes

Cardiac Cath Lab and Special Procedures Stress Testing Holter Monitoring.
Disposable 38mm foam single use silver/silver chloride radiotranslucent ECG electrode.
60 per bag, 600 per case

Model	Description	Price
BP-1400	■ Foam Electrodes	\$175.00 cs.



Foam Wet Gel Radiotranslucent Carbon ECG Electrodes

Adult radio-transparent medical ECG electrodes are designed for use during cath lab, MRI, X-Ray and thallium procedures.

Features:

Carbon snap will not interfere with Nuclear Imaging Systems
Silver/silver chloride sensing system
Exceeds AAMI specifications
All Electrodes are latex free.
300 per case (10 pouches of 30)

Model	Description	Price
BP-1039	■ Foam Wet Gel Electrodes	\$54.00 cs.



Desktop/ Portable Pulse Oximeter

The Nonin 7500FO Pulse Oximeter incorporates a Fiber Optic sensor, meeting the unique needs of MRI users for patient safety and image quality. Product highlights include:

Safe-to-use - Fiber optic sensors eliminate the possibility of RF burns in the MRI environment

Convenient - Non-magnetic pole mount system for easy access

Economical - Long-life internal rechargeable battery system and light-weight design

Customizable - pulse tone with variable pitch; audible and visual alarms for pulse rate & SpO2

Expandable - Optional 18 hour memory for software review of SpO2 & pulse rate

Flexible - Multiple interface & printing options

Pole Mount - Fits poles with 5/8" to 1 1/2" Diameter Poles sold separately below.

Warranty

3 years on pulse oximeter


1 year on internal battery

90 days on fiber optic cords and other components



Fiber Optic Sensors must be ordered separately, they are not included with the pulse oximeter.

Model	Description	Price
PO-1020	▲ Pulse Oximeter	\$1,663.89 ea.
PO-1001	▲ Infant Fiber Optic Sensor w/ 20' Cable	\$1,060.00 ea.
PO-1002	▲ Infant Fiber Optic Sensor w/ 30' Cable	\$1,165.00 ea.
PO-1003	▲ Adult Fiber Optic Sensor w/ 20' Cable	\$1,060.00 ea.
PO-1004	▲ Adult Fiber Optic Sensor w/ 30' Cable	\$1,165.00 ea.
PO-1021	● AC Adaptor	\$145.00 ea.
PO-1022	● Power cord for Adaptor	\$52.95 ea.
PO-1010	● nVision Software-Data Mgmt Cable Required	\$435.00 ea.
PO-1023	● Cable for nVision Software	\$145.00 ea.
PO-1024	● Black carrying case for pulse Oximeter	\$114.95 ea.
OX-1061	▲ 48" High Equipment Pole	\$258.00 ea.
OX-1062	▲ 60" High Equipment Pole	\$258.00 ea.



Pole mount hardware included on OX-1061 & OX-1062

NeoPulse Disposable Pulse Oximeter Wrap

100 Per Case

Features:

Made with unique, patented NeoFoam™

Soft and ultra skin friendly

Easily adjustable to fit any size

Lightweight and comfortable

Provides stable, secure attachment that won't slip

Breathable material has no unpleasant scent

Wicks away perspiration

Latex and phthalate (DEHP) free

Model	Description	Price
PO-1025	Small, 3/4" x 3"	\$101.25 cs.
PO-1026	Large, 1" x 5"	\$101.25 cs.





Latex-Free Sphygmomanometer

1.5 T Safe

Model	Description	Price
BP-1017	Adult	\$249.00 ea.
BP-1027	Large Adult	\$269.00 ea.
BP-1028	Child	\$265.00 ea.
BP-1029	Pediatric	\$262.00 ea.

Replacement Latex-Free Bulb and Valve Assembly

Model	Description	Price
BP-2020	Bulb and Valve Assembly	\$17.95 ea.

Replacement Latex-Free Cuffs and Bladders

One and Two tube bladders
Black nylon cover

Model	Description	Price
BP-2004	■ 2 Tubes, Thigh	\$45.51 ea.
BP-2009	■ 1 Tube, Thigh	\$45.51 ea.
BP-2003	■ 2 Tubes, Large Adult	\$32.45 ea.
BP-2008	■ 1 Tube, Large Adult	\$32.45 ea.
BP-2000	■ 2 Tubes, Adult	\$23.40 ea.
BP-2005	■ 1 Tube, Adult	\$23.40 ea.
BP-2002	■ 2 Tubes, Child	\$23.40 ea.
BP-2007	■ 1 Tube, Child	\$23.40 ea.
BP-2001	■ 2 Tubes, Pediatric	\$23.10 ea.
BP-2006	■ 1 Tube, Pediatric	\$23.10 ea.

SmartTourniquet

A set of high quality thigh blood pressure cuffs for use in the MRI environment.

Each SmartTourniquet kit contains two each of thigh cuffs, extension tubing, and manometer.

Latex Free

Re-Useable


3 T Conditional

Model	Description	Price
BP-1016	▲ MRI Smart Tourniquet	\$499.00 set






Oxygen Regulators / Flowmeter




Model	Description	Price
OX-1021	Nut & Nipple Style, 50 PSI Oxygen Regulator (M cylinder)	\$325.00 ea.
OX-1020	Yolk Style, 50 PSI Oxygen Regulator	\$325.00 ea.
OX-1022	Yolk Style (Most E-Cylinders), 0-15 LPM Regulator w/ Flow Control Knob	\$325.00 ea.
OX-1023	Nut & Nipple Style (Most M-Cylinders), 0-15 LPM Regulator w/ Flow Control Knob	\$332.00 ea.
OX-2002	Yolk Style, 50 PSI Preset Oxygen Regulator	\$164.95 ea.
OX-2003	Yolk Style, 0-15 LPM Regulator w/ Flow Control Knob	\$273.95 ea.
OX-6000	Click Style - Dome Handle, 0-15 LPM Regulator	\$265.35 ea.
OX-6001	Click Style - Dome Handle, 0-15 LPM Mini Regulator	\$245.00 ea.



Oxygen & Air Flow Meters

Model	Description	Price
OX-1045	Quick Connect for Ohmeda	\$33.93 ea.
OX-1046	Quick Connect for Chemtron	\$31.93 ea.
OX-1047	Quick Connect for Puritan	\$16.70 ea.
OX-1025	Oxygen Flowmeter with DISS Connection	\$79.20 ea.
OX-2007	Oxygen Flowmeter with Female Nut & Nipple DISS Connection for M Cylinders	\$72.50 ea.
OX-1053	Air Flowmeter, 0-15 LPM, with DISS Inlet Connection	\$94.50 ea.
OX-1026	Oxygen Flowmeter with Ohmeda or Hillrom Connection	\$88.95 ea.
OX-1055	Air Flowmeter, 0-15 LPM with Ohmeda Quick Connect Inlet	\$106.50 ea.
OX-1027	Oxygen Flowmeter with Chemetron Connection	\$92.50 ea.
OX-1054	Aluminum Air Flowmeter, 0-15 LPM with Chemetron Connection	\$106.50 ea.
OX-1040	Oxygen Flowmeter with Puritan Bennett Connection	\$88.95 ea.
OX-1056	Aluminum Air Flowmeter, 0-15 LPM with Puritan Bennett Connection	\$106.50 ea.
OX-1017	Manifold with Dual Flowmeters & DISS Connection	\$215.55 ea.
OX-1018	Manifold with Dual Flowmeters & Ohmeda Connection	\$215.55 ea.
OX-1019	Manifold with Dual Flowmeters & Chemetron Connection	\$215.55 ea.
OX-1041	Manifold with Dual Flowmeters & Puritan Bennett Connection	\$215.55 ea.



Complete Oxygen Cart

Includes:

E Cylinder
Cart
Cylinder Wrench
Regulator with 0-15 Flowmeter

Model	Fig	Description	Price
OX-5000	A	Round Gauge	\$559.65 ea.
OX-5001	B	Tube Gauge	\$489.00 ea.
OX-5002	C	Click Style	\$459.00 ea.
OX-5003	C	Mini Click Style	\$399.00 ea.



Oxygen Labels

Oxygen warning stickers, package of 10.

Model	Description	Price
MT-4014	Oxygen Labels	\$13.90 set

Oxygen Cylinders



Model # OX-1028 E Cylinder



Model # OX-1029 M Cylinder

CGA 540 valve is ok for
MRI use and accepts nut &
nipple regulators

Cylinder Wrenches



Model #
OX-1031



Model #
OX-1057


Oxygen Cylinders & Wrenches

Model	Description	Price
OX-1028	Oxygen E Cylinder, 25" x 4½", Accepts Yoke Regulators	\$114.26 ea.
OX-1058	Air E Cylinder, 25" x 4½", Accepts Yoke Regulators	\$208.99 ea.
OX-1029	Oxygen M Cylinder, 36¼" x 8", Accepts Nut & Nipple Regulators	\$355.00 ea.
OX-1031	Small E Cylinder Wrench	\$8.95 ea.
OX-1057	E & M Cylinder Wrench	\$16.75 ea.

■ MR Safe ▲ MR Conditional ● Not MR Safe



I.V., Ventilator & Equipment Stands



Model	Description	Price
IV-1000	2 Hook I.V. Pole, Plastic Base w/ Five Casters, Adjustable Height 59"-88"	\$248.00 ea.
IV-1003	4 Hook I.V. Pole, Plastic Base w/ Five Casters, Adjustable Height 59"-88"	\$258.00 ea.
IV-1006	2 Hook I.V. Pole, Aluminum Base w/ Five Casters, Adjustable Height 59"-88"	\$269.00 ea.
IV-1007	4 Hook I.V. Pole, Aluminum Base w/ Five Casters, Adjustable Height 59"-88"	\$279.00 ea.
OX-1007	Magellan & Tau Omni Ventilator Stand w/ Dual D/E Cylinder Holder & Top Tray, Aluminum Base, 1 1/4" Diameter Pole	\$325.00 ea.
OX-1008	Ventilator Stand w/ Dual D/E Cylinder Holders, Aluminum Base, 60" High, 1 1/4" Diameter Pole	\$337.62 ea.
OX-1059	Equipment Stand Only, 48" High, 1 1/4" Diameter Pole, Casters	\$250.00 ea.
OX-1060	Equipment Stand Only, 60" High, 1 1/4" Diameter Pole, Casters	\$250.00 ea.
OX-1061	Equipment Stand Only, 48" High, 1 1/4" Diameter Pole, Rubber Tips	\$258.00 ea.
OX-1062	Equipment Stand Only, 60" High, 1 1/4" Diameter Pole, Rubber Tips	\$258.00 ea.

Mounting Brackets

Model	Description	Price
OX-1011	Bio Medical Device Blender or Bird Blender, Wall Mount	\$37.60 ea.
OX-1012	Bio Medical Device Blender or Bird Blender, Pole Mount	\$120.35 ea.
OX-1014	Check Valve Manifold on Wall Mount Bracket 1 Input - 4 Output	\$103.20 ea.

Wall Mount
Bracket



OX-1011

Pole Mount
Bracket



OX-1012

Check Valve
Manifold



OX-1014




Aluminum Oxygen Cylinder Carts & Stands

Model	Description	Price
OX-8000	Aluminum Oxygen Cart on Wheels for D or E Cylinders, 38" H x 12" W x 10" L	\$193.25 ea.
OX-8001	Aluminum Dual Oxygen Cart on Wheels for D or E Cylinders, 38.5" High	\$279.50 ea.
OX-8002	4 Tank Aluminum Oxygen Cart for D or E Cylinders, 38.5" H x 13.75" W x 13.25" L	\$278.50 ea.
OX-8003	6 Tank Aluminum Oxygen Cart for D or E Cylinders, 40.25" H x 13.75" W x 16.75" L	\$319.00 ea.
OX-8004	Deluxe Aluminum Oxygen Cart for M/M60 Cylinders	\$369.00 ea.
OX-8005	Aluminum Wall Mount Cylinder Holder for E Cylinders	\$160.00 ea.
OX-8006	Aluminum Oxygen Stand for all Size Cylinders, Free Standing with Securing Strap	\$99.75 ea.
OX-8007	Aluminum Oxygen Cart on Wheels for M Cylinders, (8" diameter Cylinders) 38" H x 14" W x 13" L	\$252.50 ea.




Model # OX-8003
Oxygen Stand
6 Tank



Model # OX-8007
Oxygen Cart 8" Dia.
Cylinders

Model # OX-8006
Free Standing
Cylinder Holder



Model # OX-8005
Wall Mount
Cylinder Holder



Model # OX-8004
Deluxe Oxygen Cart, M/M60 Cylinders



Model # OX-8001
Oxygen Cart - Dual Tank



Model # OX-8002
Aluminum Cart for 4 Cylinders



Model # OX-8000
Aluminum Cylinder Cart



Oxygen Cylinder Carts & Stands

Model	Description	Price
OX-1000	▲ Stainless Steel Oxygen Cart for D or E Cylinders, Adjustable Handle 38½" to 44", 6" White Plastic Wheels	\$205.00 ea.
OX-1001	▲ Oxygen Cart with Heavy Duty Bed Hook for D or E Cylinders, White Powder Coat Finish	\$295.00 ea.
OX-1002	▲ Dual Tank Oxygen Cart with Heavy Duty Bed Hook for D or E Cylinders, White Powder Coat Finish	\$310.00 ea.
OX-1050	▲ 6 Tank Oxygen Stand for D or E Cylinders	\$340.70 ea.
OX-1003	▲ 4 Tank Aluminum Oxygen Stand for D or E Cylinders, White Powdered Coat Finish	\$284.40 ea.
OX-1004	▲ Oxygen Cart with Heavy Duty Bed Hook for M/M60 Cylinders, White Powder Coat Finish	\$244.00 ea.
OX-1005	▲ Plastic Floor Cylinder Stand for D or E Cylinders, Impact Resistant, Green	\$47.95 ea.
OX-1006	▲ Wall Mount Cylinder Holder for D or E Cylinders, White Epoxy Powder Coated Aluminum	\$208.00 ea.



Model # OX-1050
Oxygen Stand
6 Tank



Model # OX-1003
Oxygen Stand
4 Tank


Model # OX-1005
Cylinder Floor Stand



Model # OX-1006
Wall Mount
Cylinder Holder



Model # OX-1004
Oxygen Cart - Large Cylinder



Model # OX-1002
Oxygen Cart - Dual Tank



Model # OX-1001
Heavy Duty Cylinder Cart
with Bed Hook



Model # OX-1000
Stainless Steel Cylinder Cart



Du-O-Vac Plus, Suction System

0-350 mm hg of suction, 2 each 50 PSI power take off, 500 ml collection canister with built in filter. Powered by 50 PSI gas source

Model	Description	Price
OX-1049	Hill-Rom Wall Connection	\$589.00 ea.
OX-1032	DISS Wall Connection	\$589.00 ea.
OX-1033	Ohmeda Wall Connection	\$589.00 ea.
OX-1039	Chemtron Wall Connection	\$589.00 ea.
OX-1034	Puritan Bennett Wall Connection	\$589.00 ea.
OX-1038	10 pack Replacement Canisters, w/ frosted writing strip with 100 ml increments to 400 ml	\$89.00 pk.

Puritan



Connector

Chemtron



Connector

Ohmeda



Connector

Replacement Canisters



Du-O-Vac Plus with Regulator and Flowmeter, Suction System

Includes: main body, collection canister, O2 flowmeter, oxygen regulator.
(Does not include cylinder and stand)

Model	Description	Price
OX-1037	With Regulator & Flowmeter	\$840.00 ea.

Vacuum Slide Wall Bracket for Disposable Suction Canisters

Accepts the plastic rings that hold the disposable suction canisters.

Model	Description	Price
OX-1072	Wall Bracket	\$42.58 ea.



Du-O-Vac Plus Code Cart, Suction System

Includes: Du-O Vac Plus, Regulator, Cylinder Cart, E Cylinder
Option, 0-15 lpm Flowmeter

Model	Description	Price
OX-1035	With E Cylinder	\$1,428.35 ea.
OX-1036	With Out E Cylinder	\$1,378.35 ea.



3 Mode Continuous Vacuum Suction Regulator

It's Field Serviceable ...The unique modular design makes repair a breeze for any biomed or central supply department. The unit is made up of six (6) easy-to-replace modular components.

It's Compact ...As a result of its modular design, the unit is 25-50% smaller than other suction regulators making it the product of choice for crowded wall and rail areas.

It's Sturdy ...The case is made of strong ABS plastic. It features a smooth surface with rounded corners making cleaning easy. The case also houses the gauge for extra protection.

It's User-Friendly ...The large control knob and easy-to-read settings make vacuum adjustments simple, even in emergency situations.

Specifications:

Gauge Reading

0-200 mmHg and "Full Vacuum"

Vacuum Operating Range

0-200 mHg and "Full Vacuum"

Modes

Line-Off-Reg



Model	Description	Price
OX-7500	Tubing Nipple and NO Wall Connection (1/8" N.P.T.F.)	\$346.60 ea.
OX-7501	Tubing Nipple and DISS Hex Nut	\$355.66 ea.
OX-7502	Tubing Nipple and DISS Hand Tight	\$357.14 ea.
OX-7503	Tubing Nipple and Ohmeda Wall Connection	\$373.30 ea.
OX-7504	Tubing Nipple and Chemetron Wall Connection	\$373.30 ea.
OX-7505	Tubing Nipple and Oxequip (old style) Wall Connection	\$362.14 ea.
OX-7506	Tubing Nipple and Oxequip Medstar Wall Connection	\$415.04 ea.
OX-7507	Tubing Nipple and Puritan-Bennett Wall Connection	\$362.42 ea.
OX-7508	DISS Male and NO Wall Connection (1/8" N.P.T.F.)	\$351.96 ea.
OX-7509	DISS Male and DISS Hex Nut Wall Connection	\$361.03 ea.
OX-7510	DISS Male and DISS Hand Tight Wall Connection	\$362.51 ea.
OX-7511	DISS Male and Ohmeda Wall Connection	\$378.69 ea.
OX-7512	DISS Male and Chemetron Wall Connection	\$378.69 ea.
OX-7513	DISS Male and Oxequip (old style) Wall Connection	\$367.50 ea.
OX-7514	DISS Male and Oxequip Medstar Wall Connection	\$425.13 ea.
OX-7515	DISS Male and Puritan-Bennett Wall Connection	\$367.78 ea.
OX-7516	Ohmeda Style Trap Connection and NO Wall Connection (1/8" N.P.T.F.)	\$350.30 ea.
OX-7517	Ohmeda Style Trap Connection and DISS Hex Nut Wall Connection	\$359.36 ea.
OX-7518	Ohmeda Style Trap Connection and DISS Hand Tight Wall Connection	\$360.84 ea.
OX-7519	Ohmeda Style Trap Connection and Ohmeda Wall Connection	\$377.03 ea.
OX-7520	Ohmeda Style Trap Connection and Chemetron Wall Connection	\$377.03 ea.
OX-7521	Ohmeda Style Trap Connection and Oxequip (old style) Wall Connection	\$370.55 ea.
OX-7522	Ohmeda Style Trap Connection and Oxequip Medstar Wall Connection	\$418.75 ea.
OX-7523	Ohmeda Style Trap Connection and Puritan-Bennett Wall Connection	\$366.11 ea.



Multi-Modality Markers

Product Benefits:

Provides clear reference points
Shows on all sequences
Does not indent into soft tissue
No artifacts
50 Per Container

For further details visit www.mriequip.com or call customer service

Model	Description	Price
BP-1015	■ Multi-Modality Markers	\$100.50 cs.
BP-1037	■ Multi-Modality Markers, Solid	\$105.60 cs.
BP-1030	■ Multi-Modality Markers for Nuclear Medicine/Pet	\$122.10 cs.

"Great Markers! Perfect size and show on both T1 and T2 weighted images. Also, the markers do not indent the skin." - Open MRI Specialists of Stockbridge



Burnett BiDirectional TMJ Device

ACQUIRE DYNAMIC MR IMAGES

Key Features

Precisely opens and closes the jaws.
Audible click indicates 1mm movement of mouthpiece.
Greater bite area
Improved strength and tracking method
Scale indicating distance between mouthpieces up to 55mm
Disposable mouthpieces for single patient use

Benefits

Simple one-handed operation
Reduced motion artifacts
Virtually eliminates device flex for better images

Warranty
1 Year

Model	Description	Price
MO-1000	Device with 3 Sets of Disposable Mouthpieces	\$1,018.40 ea.
MO-1001	10 Sets of Disposable Mouthpieces	\$56.20 cs.



SmartTourniquet

A set of high quality thigh blood pressure cuffs for use in the MRI environment.
Each SmartTourniquet kit contains two each of thigh cuffs, extension tubing, and manometer.
Latex Free
Re-Useable
3T Conditional

Model	Description	Price
BP-1016	▲ MRI Smart Tourniquet	\$499.00 set



Fiducial Markers

A powerful new tool, engineered to provide features never before available on an elongated fiducial marker. Designed for higher visibility and greater stability in tissue, the Fiducial Marker lets you target tumors with confidence, whether using CT, MR, Portal, or Ultrasound imaging. Elongated gold fiducial markers that provide accuracy, stability in soft tissue, and greater visibility in MRI, CT Portal and Ultrasound. Features a 99.95% gold construction and unique snake-like design. Its structure enables it to flex with soft tissue and it will not unravel like other coil-style markers.


Since it is longer than a standard fiducial, fewer markers need to be used, thus minimizing the number of needle sticks necessary. Additionally, the innovative design greatly reduces bending problems associated with the injection process as compared to other linear soft tissue markers. Sterile, 1 per pack

Preloaded 20 cm Needle, 17 Gauge

Model	Description	Price
BP-1050	1.15 mm x 1 cm Fiducial	\$162.50 ea.
BP-1051	1.15 mm x 2 cm Fiducial	\$175.00 ea.
BP-1052	1.15 mm x 3 cm Fiducial	\$187.50 ea.

Preloaded 20 cm Needle, 18 Gauge

Model	Description	Price
BP-1053	0.85 mm x 1 cm Fiducial	\$162.50 ea.
BP-1054	0.85 mm x 2 cm Fiducial	\$175.00 ea.
BP-1055	0.85 mm x 3 cm Fiducial	\$187.50 ea.



Apparel / Hairnet Dispensers

Keep protective and sanitary apparel in an easy-to-reach location.

Features:

Can be used for anything from headset covers (**RA-1006**) and shoe covers to wet wipes and goggles
Mounts easily wherever you most need it with tape
Top-loading for easy refilling
Clear acrylic to easily see contents
Center hole for pulling items out

Model	Description	Price
PC-1027	Small - 8.5" x 11" x 4"	\$54.84 ea.
PC-1103	Medium - 15" x 17" x 6"	\$92.02 ea.



Large 3-in-1 Acrylic Apparel Dispenser

Keep disposable protective wear where it is needed most.

Features:

Ear plugs (**RA-1009**), headset covers (**RA-1006**) and eye shades (**RA-6000**) at lab entry way in a wall-mountable dispenser
Acrylic front for easy viewing

Model	Description	Price
PC-1086	20 inches wide, 18" tall and 6" deep	\$212.86 ea.





Eight Compartment Storage Bin

Storage bins are an easy and affordable way to organize your lab. Bins are ready to use; no assembly required. Attractive design instantly improves lab appearance. Contains eight bins. This supply bin is 24" long, 12" tall and 10" deep with 8 compartments each 4.5" wide. Optional, hinged lids come in multiple colors, choose from below. Ships ready to use, mounting tape included and holds 8-10lbs

Model	Description	Price
PC-1203	Storage Bin	\$110.61 ea.
PC-2203-CL	Storage Bin with Colored Lid	\$138.50 ea.



Clear, smoke, neon red, orange, green and blue



Sixteen Adjustable Compartment Storage Bin

Storage bins are an easy and affordable way to organize your lab. Bins are ready to use; no assembly required. Attractive design instantly improves lab appearance. 16 adjustable dividers to create the right sized bins for your supplies. This supply bin is 12" long, 11.75" tall and 7.25" deep. Bright white PVC plastic makes the area look clean and supplies stand out for easy use. Durable PVC construction with smooth, non-scratching surface. Mounting tape included and holds 8-10lbs. Optional, hinged lids come in multiple colors, choose from below.

Model	Description	Price
PC-1012	Storage Bin	\$101.32 ea.
PC-2012-CL	Storage Bin with Colored Lid	\$115.26 ea.



Clear, smoke, neon red, orange, green and blue



High Capacity Storage Bin with Shelves

Storage bins are an easy and affordable way to organize your lab. Bins are ready to use; no assembly required. Attractive design instantly improves lab appearance. 14 adjustable dividers to create the right sized bins for your supplies. Two shelves to store bottle, gloves or other larger items. This supply bin is 24" long, 20" tall and 10" deep. Bright white PVC plastic makes the area look clean and supplies stand out for easy use. Durable PVC construction with smooth, non-scratching surface. Mounting tape included and holds 8-10lbs. Optional, hinged lids come in multiple colors, choose from below.

Model	Description	Price
PC-1017	Storage Bin	\$166.38 ea.
PC-2017-CL	Storage Bin with Colored Lid	\$204.49 ea.



Clear, smoke, neon red, orange, green and blue



Supply Boxes

Walking in and out of the lab, or dropping off samples, these all purpose lab boxes make a great addition to the entry. Clear acrylic sides allow for easy viewing. PVC and acrylic; mounts to wall with tape. Use at lab entrance for safety glasses, samples, or supplies. Choose from clear or colored acrylic for color coding. Available in two sizes and multiple colors, choose below. Convenient quick safety glass or sample storage while walking in and out of the lab, or dropping off samples. Use at lab entry and exit. Mounting tape included and holds 8-10lbs

Model	Description	Price
PC-2056-CL	Small - 8.85" W x 5.875" T x 6.1" D	\$36.25 ea.
PC-2162-CL	Large - 8.85" W x 9" T x 6.1" D	\$45.55 ea.



Clear, smoke, neon red, orange, green and blue



Glove Box Holders

Made of white styrene plastic construction that is tough and durable. Affordable and patient friendly way to encourage safety within the MRI environment. Wall mountable with double sided tape not included. These are only single glove box holders, but double and triple holders are also available in all designs please call for pricing on them. Size: 10.75" x 5.5" x 3.88"

Model	Description	Price
PC-1520	Sea Turtle Glove Box Holder	\$12.86 ea.
PC-1521	Pelican Glove Box Holder	\$12.86 ea.
PC-1522	Seagull Glove Box Holder	\$12.86 ea.
PC-1523	Cat Lovers Glove Box Holder	\$12.86 ea.
PC-1524	Dog Lovers Glove Box Holder	\$12.86 ea.
PC-1525	Dolphin Glove Box Holder	\$14.75 ea.



Side Loading Acrylic Glove Box Holders

Right or left filling. These glove box holders are perfect for your organized lab. Available in three different configurations. The compartments are designed to fit most types of glove boxes. Prefabricated holes make mounting this holder in your lab easy. Mounting tape included and holds 8-10lbs

Model	Description	Price
PC-1504	Single, 1 Compartment 5.6" x 10.4" x 4.4"	\$43.69 ea.
PC-1503	Double, 2 Compartments 10" x 10.4" x 4.4"	\$66.92 ea.
PC-1499	Triple, 3 Compartments 16.4" x 10.4" x 4.4"	\$70.64 ea.





MRI CORE DX Cart

High capacity Deluxe Cart offers versatility along with more than 17ft of moveable storage. This non-ferrous, non-magnetic cart is designed especially for the MRI Lab to easily move supplies in and out of patient areas. Four side shelves provide roomy storage for supplies along with built-in glove and Kimwipe dispensers. A multi compartment tool holder on the back holds lab books, patient records, pens and other useful items. Deluxe MRI Lab Cart with three drawers, four shelves, two glove and wiper dispensers and one large tool and log book holder. One large compartment in the bottom. Two locking casters for braking. High Density Polyethylene (HDPE) top and bottom; hinged acrylic door. 29" Long x 17" Wide x 34.5" High

Model	Description	Price
PC-2470-CL	Core DX Cart	\$870.94 ea.



Clear, smoke, neon red, orange, green and blue



Caymen Lab Island

The most flexible laboratory work space available. Caymen Lab Island™ combines a large work surface on wheels with maximum useful storage compartments for an ideal blend of covered storage and movable work surface. Made to handle the tough laboratory environment, the high density polyethylene construction gives years of valuable service. Lab Islands™ keeps your lab neat and clean on even the busiest days. Easy-to-use, no-swing-out doors improves safety and saves space. Secure lock protects contents. Lockable- locking mechanism is made of brass. 4" casters. Large capacity, private transportation. 4' Long x 2' Wide x 34.75" High

Model	Description	Price
PC-1475	Caymen Lab Island	\$1,924.07 ea.



Locking Roll-Top Cart

A durable high density polyethylene frame made just for the harsh laboratory environment and easy-to-use roll top door make this cart a great choice. Smooth rolling tambour door features plastic slates to travel smooth up and down the tracks. Top is also HDPE for chemical resistance so it gives years of good looking service. Casters provide easy walking speeds to help move around the laboratory. Reliable push-brakes secure cart when necessary. Built-in handle for easy moving. Lockable cart with disappearing roll top door. Six removable drawers for use and storage. Lockable- locking mechanism is made of brass. 11.5" Long x 16.5" Wide x 35" High

Model	Description	Price
PC-1471	Locking Roll-Top Cart	\$774.27 ea.

More Carts Available at
www.MRlequip.com



Space-Saving Cart

Our Space Saver Lab Cart is only 17" wide, allowing it to easily navigate through narrow openings and makes excellent use of its 35 inches of height. Use it to store your favorite lab supplies and roll it along in this space saving cart with the convenient built in handle. Our Space Saver Cart has easy-rolling casters (two are locking) and has three deep drawers and one large shielded compartment in the bottom, great for bulk storage. The sides are made from ABS, and the top and bottom HDPE. 17" Long x 17" Wide x 34.5" High

Model	Description	Price
PC-2468-CL	Space-Saving Cart	\$1,092.34 ea.



Clear, smoke, neon red, orange, green and blue



MRI Non-Magnetic Light Duty Cart

If you need a light-weight cart with heavy-weight function the MRI Light Duty Cart is for you. This super easy to move, ultra-light cart, makes easy work of moving your supplies to your patient's side. This light duty cart has four removable drawers, four side shelves, built-in clipboard and a bulk storage area.

Features:

- Nonferrous light duty cart
- 4 removable drawers
- 4 side shelves
- Built-in clipboard
- 1 bulk storage area
- 2 3/8 inch non-ferrous braking casters
- 17 x 33 1/8 x 15 1/2 inches (width, height, depth)
- 19 pounds

Model	Description	Price
PC-1526	▲ Light Duty Cart	\$478.69 ea.



Access / Egress Organizer

Storage for personal gear or for entry or exit. Space for gloves, safety glasses, protective clothing and changing area supplies. Light, white, room-for-everything and built for maximum ease in and out of the lab. Hinged acrylic doors allow easy stocking and viewing of supplies. Fabricated design to hold lab gear and keep your clothes and personal items off the floor. Entire organizer is wall mountable and easy to install with built-in key holes, so you'll go from the office to the lab in no time.

Features one shelf, four safety glass cubbies, 3 covered compartments, a glove box dispenser and three hooks. Mounting hardware NOT included. 18.125" Long x 7.2" Wide x 19.875" High

Model	Description	Price
PC-1137	Access / Egress Organizer	\$250.04 ea.





Small Extremity General Purpose Coil

Primary Applications:

Small extremities
Wrist
Hands, Fingers
Forearms

The coil is ideal for imaging a variety of maladies to include, but are not limited to, carpal tunnel, nerve compression, ligament injuries, tendon abnormalities, tumor detection, bone injuries and scarring within the wrist.

Coil Coverage (FOV) \approx 10 cm (3.9 inches)

Coil Dimensions:

Length: 23.5 cm (9.25 inches)
Width: 14 cm (5.51 inches)
Height: 15 cm (5.90 inches)
Weight: 1361 grams (3 lbs)
Cable Length: 175 cm (69 inches)

1.5T

1.5T Small Extremity General Purpose Coil is a Transmit/Receive quadrature coil that has been optimized for specialty high-resolution musculoskeletal imaging. This design is notable for its demonstrated long-term reliability; low transmit radio frequency (RF) power, high image uniformity, and high image Signal-to-Noise Ratio (SNR).

System Requirements:

Compatible on the following MR systems (GE Signa 1.5T system's that include HD, HDx, HDi, VIBRANT, MR450 and MR450w (70cm wide bore), also includes systems that have been upgraded to this level.

Regulatory Compliance: (Approved)

The 1.5T single channel wrist coil is an FDA cleared and a CE-compliant device that satisfies regulations regarding Electro-Magnetic Compatibility (EMC) and Electro-Magnetic Interference (EMI), pursuant to IEC-601.

3T

3.0T Small Extremity General Purpose Coil is a Transmit/Receive quadrature coil that has been optimized for specialty high-resolution musculoskeletal imaging. This design is notable for its demonstrated long-term reliability; low transmit radio frequency (RF) power, high image uniformity, and high image Signal-to-Noise Ratio (SNR).

System Requirements:

Compatible on the following MR systems (GE Signa 3.0T system's that include HD, HDx, HDxt MR750. Also includes GE MR systems that have been upgraded to this level.

Regulatory Compliance: (Approved)

The 3.0T single channel wrist coil is an FDA cleared and a CE-compliant device that satisfies regulations regarding Electro-Magnetic Compatibility (EMC) and Electro-Magnetic Interference (EMI), pursuant to IEC-601

Warranty:

1 yr from date of shipment

Model	Description	Price
CO-1000	1.5 T Coil	\$13,500.00 ea.
CO-1001	3.0 T Coil	\$22,500.00 ea.



MRI Universal Oak Wood Coil Cart

Solid oak coil cart with 4" casters (except 8' carts) must add 4" to height for all carts, casters not available on carts 8' and over. Can be custom made to any size and number of shelves. Burgundy nylon cover (must be ordered separately) attaches to cart with hook and loop velcro. Top is considered as one of the shelves.



2 Shelf



5 Shelf



3 Shelf

Model	Size D x H x L	# of Shelves	Price No Cover	Add Cover	Add Handles
FR-2011	▲ 2' x 3' x 2'	2	\$736.00 ea.	\$145.00 ea.	\$33.00 ea.
FR-2012	▲ 2' x 3' x 3'	2	\$763.00 ea.	\$145.00 ea.	\$33.00 ea.
FR-2013	▲ 2' x 3' x 4'	2	\$853.00 ea.	\$145.00 ea.	\$33.00 ea.
FR-2000	▲ 2' x 4' x 2'	3	\$778.00 ea.	\$145.00 ea.	\$33.00 ea.
FR-2001	▲ 2' x 4' x 3'	3	\$845.00 ea.	\$145.00 ea.	\$33.00 ea.
FR-2002	▲ 2' x 4' x 4'	3	\$895.00 ea.	\$145.00 ea.	\$33.00 ea.
FR-2003	▲ 2' x 6' x 2'	4	\$956.00 ea.	\$155.00 ea.	\$33.00 ea.
FR-2004	▲ 2' x 6' x 3'	4	\$1,028.00 ea.	\$155.00 ea.	\$33.00 ea.
FR-2005	▲ 2' x 6' x 4'	4	\$1,082.00 ea.	\$155.00 ea.	\$33.00 ea.
FR-2006	▲ 2' x 8' x 2'	5	\$1,000.00 ea.	\$155.00 ea.	\$33.00 ea.
FR-2007	▲ 2' x 8' x 3'	5	\$1,083.00 ea.	\$155.00 ea.	\$33.00 ea.
FR-2008	▲ 2' x 8' x 4'	5	\$1,166.00 ea.	\$155.00 ea.	\$33.00 ea.



MRlequip.com proudly introduces a revolutionary patent pending draping system for MRI breast coils.

These drapes will minimize the daily cleaning routine of coils as well as ease the mind of clientele. The drapes provide comfort to patients, while providing protection from contagious diseases such as MRSA. Finally, your staff will be delighted with your decision to simplify clean up between patients, thus achieving a greater throughput.

AccuFit MRI Sentinelle Breast Coil Drape

Provides complete coverage and protection of all Sentinelle Breast Coil Products. Total Coverage - added an extra 60" in length now covering the entire table which eliminates costly usage of sheets and provides quick procedure turn around.



Sentinelle MRI Breast Coil Disposable Drape

25 per case

Compatible with:

Siemens MAGNETOM® Avanto
 MAGNETOM Espree
 MAGNETOM Symphony
 GE Signa® HDx or HDi MRI
 Toshiba 1.5T Vantage Atlas™
 Vantage Titan™

May fit coils other than those listed



Model	Description	Price
BS-1011	■ Sentinelle Breast Coil Drape (25 per case)	\$149.75 cs.



AccuFit Sentinelle Large Headrest Cover

This product fits the large 5 inch high headrest utilized with the Sentinelle Breast Coil Table and stand alone coils.

100 per case

Model	Description	Price
BS-1012	■ Sentinelle Headrest cover (100 per case)	\$149.00 cs.

Please Call for a Free Sample.



AccuFit Headrest Disposable Covers

100 per case



Model	Description	Price
BS-1003	■ AccuFit Headrest Covers (100 per case)	\$149.00 cs.


AccuFit MRI Breast Coil Disposable Drape Style One

50 per case

Compatible with:

Philips Achieva 1.5T
Philips Intera 1.5T
Philips Achieva 3.0T
Philips Intera 3.0T
GE Excite 1.5T
GE Excite HD 1.5T
GE Excite HDx 1.5T
GE Signa LX 1.5T

May fit coils other than those listed



Model	Description	Price
BS-1004	■ AccuFit Breast Coil Drape (50 per case)	\$89.50 cs.

AccuFit MRI Breast Coil Disposable Cover Style Two

Compatible with:

Siemens 1.5T Avanto/Espree 7 Channel
Siemens 1.5T Symphony 7 Channel
Philips 1.5T Achieva 7 Channel
Philips 3.0T SENSE
GE 8 Channel coil
GE 1.5T Excite 7 Channel
GE 3.0T 7 Channel Breast
CP Siemens Breast Array
Insight

May fit coils other than those listed



Model	Description	Price
BS-1007	■ Breast Coil Drape Style Two (100 per case)	\$369.00 cs.
BS-1008	■ Breast Coil Drape Style Two (20 per case)	\$89.50 cs.

Please Call for a Free Sample.



Are your current MRI supplies certified?

MRlequip.com proudly introduces the first of its kind non-ferromagnetic MRI certified biopsy tray to include both a scalpel and needle. The first biopsy tray significantly reduces your fine institution of liability, when used in the prescribed fashion. Please call for MRI certification for future purchases. Thanks to your radiology family, larger production runs have lowered the price to \$399.00 per case, that's just **\$39.90** per tray.



Non-Ferromagnetic Certified Scalpel sold separately on page 54.

MRI Non-Ferromagnetic "Certified" Biopsy Tray

Featuring a Non-Ferromagnetic Scalpel and Non-Ferromagnetic Needle

SAFETY: MRlequip.com is committed to bring to you the only MRI Non-Ferromagnetic Certified Biopsy Tray. Designed and certified within the guidelines of the new ASTM standards in effect so your facility no longer has to assume the liability for any safety accidents due to items made with ferrous materials.

USAGE FLEXIBILITY WITH LESS WASTE: Based on customer feedback for increased ordering and procedure flexibility you'll note that we packaged the main biopsy components separately as well as in a MRI biopsy tray. Trays can also be used as a MRI localization tray.

UNIQUE DESIGN: MRlequip.com offers the only MRI Non-Ferromagnetic certified disposable scalpel and MRI Non-Ferromagnetic certified needle.

COST EFFECTIVE: This product is a first of its kind. Developed by technologists who perform the procedures. It is cost effective, practical and safe to use. **The trays are completely sterile and disposable.**

MR Conditional to: 3-Tesla

MRI Biopsy Tray is SAFE in the MR environment when used in the prescribed fashion.

The trays have a 2 year expiration date.



Tray Contents:

- 1 - 5cc syringe luer lock
- 2 - 10cc syringe luer lock
- 1 - 20cc syringe luer lock
- 1 - 22-gauge non-ferromagnetic injection needle 9cm
- 1 - # 11 blade non-ferromagnetic scalpel with plastic handle
- 1 - plastic forcep fine tip
- 6 - 4x4's 12 ply
- 2 - 2x2's 12ply
- 1 - ChloraPrep® 3ml
- 1 - Betadine swabs (pkg 3)
- 1 - sterile pen with ruler (to write on labels)
- 1 - sheet pre-labeled labels
- 1 - Alcohol prep
- 1 - Steri-strip ½ x 4 inch (3 per card)
- 1 - 4oz specimen cup
- 1 - towel 17 x 20 white
- 1 - prep tray 10x5x2

Model	Description	Price
BS-1000	▲ MRI Certified Biopsy Tray	\$399.00 cs.

MRI Non-Magnetic 20 G Needle

SAFETY: MRI Equip.com is committed to bring to you the only MRI Non-Ferromagnetic Certified Needle. Designed and certified within the guidelines of the new ASTM standards in effect so your facility No Longer has to assume the liability for any safety accidents due to items made with ferrous materials. Even the "BEST" employees can make a human error. It only takes one incident to hurt someone. **USAGE:** This MRI Non-Ferromagnetic Certified Needle is for utilization as a puncture or injection needle. Ideal for interventional procedures. **MR Conditional to: 3-Tesla.** MRI Non-Ferromagnetic Certified Needle is SAFE in the MR environment when used in the prescribed fashion for all Tesla Strengths



Model	Description	Price
BS-2000	▲ 20 Gauge 7 cm Needle	\$94.48 ea.
BS-2001	▲ 20 Gauge 9 cm Needle	\$94.48 ea.
BS-2002	▲ 20 Gauge 5 cm Needle (Qty 5 Min.)	Call

■ MR Safe ▲ MR Conditional ● Not MR Safe



#11 MRI Non-Ferromagnetic "Certified" Disposable Scalpel

This scalpel is safe to use for any interventional procedure for the initial puncture of the skin. Sterile and disposable.

SAFETY: MRlequip.com is committed to bring to you the only MRI Disposable Scalpel. Designed and certified within the guidelines of the new ASTM standards in effect so your facility no longer has to assume the liability for any safety accidents due to items made with ferrous materials.

Model	Description	Price
BS-1001	▲ MRI Certified Scalpel / 4 per case	\$95.00 cs.

MRI Non-Magnetic Titanium Round Bayonet Forceps

Model	Description	Price
TW-7021	▲ 8.25" Smooth Straight 1mm tip	Call


Scissors

Model	Description	Price
MT-6001	▲ 1¾" Blade	\$95.00 ea.

Disposable Instruments

Fig	Model	Description	Price
A	TW-5006	■ 5" Dressing Forceps, 10 Per Case	\$19.75 cs.
A	TW-5014	■ Sterilized 5" Dressing Forceps, 10 Per Case	\$19.90 cs.
A	TW-5021	■ Sterilized 5" Dressing Forceps, 100 Per Case	\$92.80 cs.
B	TW-5000	▲ 9½" Sponge Forceps, 10 Per Case	\$43.45 cs.
B	TW-5008	▲ Sterilized 9½" Sponge Forceps, 10 Per Case	\$70.65 cs.
B	TW-5015	▲ Sterilized 9½" Sponge Forceps, 100 Per Case	\$537.55 cs.
C	TW-5003	▲ 5¼" Straight Hemostat, 10 Per Case	\$37.95 cs.
C	TW-5011	▲ Sterilized 5¼" Straight Hemostat, 10 Per Case	\$32.75 cs.
C	TW-5018	▲ Sterilized 5¼" Straight Hemostat, 100 Per Case	\$230.00 cs.
D	TW-5004	■ Foam Nose Clip, 10 Per Case	\$21.95 cs.
D	TW-5012	■ Sterilized Foam Nose Clip, 10 Per Case	\$19.95 cs.
D	TW-5019	■ Sterilized Foam Nose Clip, 100 Per Case	\$67.20 cs.
E	TW-5005	■ 4" Narrow Jaw Utility Forceps, 10 Per Case	\$15.95 cs.
E	TW-5013	■ Sterilized 4" Narrow Jaw Utility Forceps, 10 Per Case	\$19.50 cs.
E	TW-5020	■ Sterilized 4" Narrow Jaw Utility Forceps, 100 Per Case	\$76.80 cs.
F	TW-5002	▲ 5¼" Curved Hemostat, 10 Per Case	\$34.20 cs.
F	TW-5010	▲ Sterilized 5¼" Curved Hemostat, 10 Per Case	\$32.75 cs.
F	TW-5017	▲ Sterilized 5¼" Curved Hemostat, 100 Per Case	\$220.50 cs.
G	TW-5001	▲ 5" Non-Perforating Towel Clamp, 10 Per Case	\$31.95 cs.
G	TW-5009	▲ Sterilized 5" Non-Perforating Towel Clamp, 10 Per Case	\$48.80 cs.
G	TW-5016	▲ Sterilized 5" Non-Perforating Towel Clamp, 100 Per Case	\$363.20 cs.

■ MR Safe ▲ MR Conditional ● Not MR Safe






TW-2000



TW-2001



TW-2002



TW-2003



TW-2020



TW-2021



TW-2022



Titanium Tweezers

These tweezers are 100% anti-magnetic & manufactured with Titanium material which is 100% non-corrosive and therefore are immune to rust. Titanium tweezers are 45% lighter than those made of stainless steel. Titanium as compared to stainless steel is much harder therefore the utility of the tweezers is greatly enhanced.

Model	Description	Price
TW-2000	▲ High Precision Tweezer, 4 1/4" Long	\$136.00 ea.
TW-2001	▲ High Precision Tweezer, 4 1/2" Long	\$136.00 ea.
TW-2002	▲ High Precision Tweezer, 4 3/4" Long	\$136.00 ea.
TW-2003	▲ High Precision Tweezer, Curved, 4 3/4" Long	\$136.00 ea.
TW-2020	▲ Micro Suture Tying Forcep, 4" Long	\$159.50 ea.
TW-2021	▲ Micro Point Corneal Forcep, 1x2 Teeth, 3 3/4" Long	\$159.50 ea.
TW-2022	▲ Micro McPherson Forcep, 3 1/2" Long	\$159.50 ea.
TW-2030	▲ Jacobson Needle Holder with Lock, 6 1/4" Long	\$209.00 ea.
TW-2031	▲ Micro Needle Holder, 4 1/4" Long	\$209.00 ea.
TW-2032	▲ Castroviejo Needle Holder without Lock, 6 1/2" Long	\$214.75 ea.
TW-2033	▲ DeBakey Atrauma Tissue Forcep, 8" Long	\$207.75 ea.

TW-2030




TW-2031



TW-2032



TW-2033



MRI Single Basin

Stainless steel water basin with 2" plastic twin wheels. 37" tall, 14" water bowl



Model	Description	Price
BA-1000	▲ Single Basin	\$303.05 ea.

Mayo Stand

Stainless Steel
Adjusts from 31" to 51"
2" twin plastic wheels



Model	Description	Price
TA-1000	▲ 13" x 19" Tray Size	\$585.00 ea.
TA-1001	▲ 16" x 21" Tray Size	\$640.00 ea.



Laryngoscope

Introducing the first fiber-optic incubating laryngoscope set designed specifically for use within the MR room.

Long Life - Laryngoscope handle developed from selected low-magnetic materials. Proprietary low-magnetic signature lithium battery delivers 4 times the useful life of a standard alkaline battery.

High Intensity - A 3.5V halogen lamp yields a high-intensity, fiber-optic delivered light beam.

Compatible - Works with a broad range of Fiber-optic Laryngoscope Blades.

Model		Description	Price
MO-1100	▲	Laryngoscope Handle, Adult	\$394.30 ea.
MO-1122	▲	Laryngoscope Handle, Pediatric	\$394.30 ea.
MO-1101	▲	Lithium Battery, 2 pack	\$337.00 pk.
MO-1102	▲	American Macintosh 0 Blade	\$194.25 ea.
MO-1103	▲	American Macintosh 1 Blade	\$194.25 ea.
MO-1104	▲	American Macintosh 2 Blade	\$194.25 ea.
MO-1105	▲	American Macintosh 3 Blade	\$194.25 ea.
MO-1106	▲	American Macintosh 3.5 Blade	\$194.25 ea.
MO-1107	▲	American Macintosh 4 Blade	\$194.25 ea.
MO-1108	▲	English Macintosh 1 Blade	\$194.25 ea.
MO-1109	▲	English Macintosh 2 Blade	\$194.25 ea.
MO-1110	▲	English Macintosh 3 Blade	\$194.25 ea.
MO-1111	▲	English Macintosh 4 Blade	\$194.25 ea.
MO-1112	▲	Miller 0 Blade	\$194.25 ea.
MO-1113	▲	Miller 1 Blade	\$194.25 ea.
MO-1114	▲	Miller 2 Blade	\$194.25 ea.
MO-1115	▲	Miller 3 Blade	\$194.25 ea.
MO-1116	▲	Miller 4 Blade	\$194.25 ea.
MO-1117	▲	Wisconsin 2 Blade	\$194.25 ea.
MO-1118	▲	Wisconsin 3 Blade	\$194.25 ea.
MO-1119	▲	Wisconsin 4 Blade	\$194.25 ea.
MO-1120	▲	C Halogen Light Bulbs, (2 pack)	\$73.25 pk.



**Handle Does NOT
Include Battery**

Warranty
1 year



Macintosh & Miller Disposable Blades

The all-new disposable laryngoscope blade which feature the Classic blade shape. The disposable blade avoids cross-infection from one patient to the next. Latex free. Compatible with handles above - MO-1100 and MO-1122.

Model		Description	Price
MO-1200	■	Macintosh #1 (25 per case)	\$192.00 cs.
MO-1201	■	Macintosh #2 (25 per case)	\$192.00 cs.
MO-1202	■	Macintosh #3 (25 per case)	\$192.00 cs.
MO-1203	■	Macintosh #4 (25 per case)	\$192.00 cs.
MO-1204	■	Miller #0 (25 per case)	\$192.00 cs.
MO-1205	■	Miller #1 (25 per case)	\$192.00 cs.
MO-1206	■	Miller #2 (25 per case)	\$192.00 cs.
MO-1207	■	Miller #3 (25 per case)	\$192.00 cs.
MO-1208	■	Miller #4 (25 per case)	\$192.00 cs.





Accu-Draw IV Clamp

MRI ACCU-DRAW™ provides one-handed access to non-sterile vials while maintaining sterility. ACCU-DRAW keeps vial labels visible at all times, and angles vials so they can be completely emptied. Made of brass, plastic, and aluminum. How it Works: The durable plastic holder accommodates two single-dose and three multi-dose vials at a time. Clamps onto any table surface. Vials that run out can be replaced with fresh ones at any time during the procedure. You'll know vials are secure because they click and lock into place. Holds vials at the perfect angle, allowing complete emptying and reducing medication waste.

Model	Description	Price
FO-1117	Accu-Draw IV Clamp	\$257.85 ea.



MRI Non-Ferromagnetic 3 Gal. Sharps Container

Has a dual-purpose lid, which allows the container to remain open for code situations or in a free-fall position for normal, one-handed operation. Lid features a needle key for safe removal of needles when necessary. Safely stands on table top or cart.



MRI Non-Ferromagnetic 2 Gal. Sharps Container

Provides open access for disposal of a variety of shapes and sizes of sharps, (i.e. needles and syringes, vials, test tubes and glass slides). Lid features a needle key for safe removal of needles when necessary. Can be wall mounted or safely stands on table top or cart.



MRI Non-Ferromagnetic 3 Qt. Sharps Container

The rotating cylinder top permits convenient, one-handed disposal of sharps. The rotating lid automatically drops needles and syringes in a horizontal position to maximize capacity. Lid remains closed when not in use. A built-in needle key permits the safe removal of needles when necessary. Units lock onto a wall bracket and are removed with a master bracket key.

Wall Safe Bracket for SH-1014 with Key



Wall Safe Bracket for SH-1016 with Key



Model	Description	Price
SH-1015	■ 3 Gallon, 6 Per Case	\$83.50 cs.
SH-1014	■ 2 Gallon, 5 Per Case	\$49.95 cs.
SH-1016	■ 3 Quart, 6 Per Case	\$48.90 cs.
SH-1019	▲ Wall Safe with Key for SH-1014	\$17.31 ea.
SH-1018	▲ Wall Safe with Key for SH-1016	\$13.60 ea.



Clipboard

100% Aluminum clipboards are tough as nails and last for years in any environment. Holds 5 Sheets of paper.

Set of 3

Model	Description	Price
CL-1000	▲ Clipboard (set of 3)	\$45.00 set

Non-Magnetic Paper Clips

Comes in assorted colors

Medium - 1", 500 per box

Large - 1 3/8", 200 per box

X-Large - 2 3/8", 50 per box

Model	Description	Price
MT-1702	■ Medium	\$17.75 box
MT-1701	■ Large	\$13.50 box
MT-1700	■ X-Large	\$11.55 box

Bic Round Stic Pen

Medium point

10 per box

Have your name custom imprinted on pen with 500 minimum purchase and \$50.00 set-up fee.

Model	Description	Price
MT-1801	▲ Pens	\$8.90 box

Non-Magnetic Safety Padlock

Plastic shackle helps prevent electrocution and arc flashing, also prevents electrical current from shackle to key

All components are non-sparking and non-magnetic - plastic, copper, brass 1 1/2" wide 1 3/4" tall Xenoy green or red body, plastic shackle w/1 1/2" clearance

Designed exclusively for Lockout/Tagout applications

Durable, lightweight, non-conductive Xenoy lock body


Compliance with OSHA "one employee, one lock, one key" directive

Model	Description	Price
PL-1000	Green	\$18.95 ea.
PL-1001	Red	\$18.95 ea.

Storage and Transport Cases

Holes on each side of latch opening allow for padlocks or other securing devices. Top storage area for reaching smaller items without having to open the case. Overall Outside Box Dimensions: 15 1/4"W x 10 1/8"H x 7 5/8"D. Removable top tray. Tray Outside Dimensions: 13 5/8"W x 3"H x 7 1/8"D. Large carry handle. Special rubber sealing gasket keeps contents dry and clean.

Model	Description	Price
MT-1609	Orange	\$36.81 ea.





Ceiling Tile Art and Decorative Wall Murals

How can you have a blue sky in a room with no windows? How can you make your ceilings and walls come to life? With dropped ceiling tile art and decorative wall murals!

Ceiling Art is scientifically proven to:

- Lower Blood Pressure
- Reduce Stress/Anxiety
- Decrease Muscle Tension
- Reduce Pain
- Re-Direct Negative Thought



There have been hundreds of research reports showing the positive effects of incorporating healing art and biophilic friendly products into healthcare design. Biophilic design accompanies green building perfectly because its products are sustainable, blurring the boundaries between building and landscape to relax and soothe the soul. Even though people may be happy in a suburban bungalow or urban high-rise, we are most at peace when walking in a park, looking at the ocean or hiking in the woods. Biophilic spaces can reduce unnecessary stress and depression, enhance positive social response, speed up recovery times, provide medical patient distractions, increase staff productivity, and even make tight, claustrophobic spaces seem wide open. Best of all, our products are maintenance free and will provide years of enjoyment for customers and staff alike. With a ceiling mural overhead, you will immediately SEE the advantages.



Hospital Lobby Ceilings & Wall Murals

Custom, dropped ceiling tile, light diffuser panels and vinyl wall murals that easily integrate into drop ceilings and walls of hospital lobbies and waiting rooms. These no maintenance, patient friendly dropped ceiling tile installations, are perceived instantly upon viewing as a sign of caring and deeply appreciated by the patient as well as the caregivers themselves. Sky ceilings are becoming more popular in facility and office design because they are visually therapeutic and create an open feeling.

Imaging Suite Ceilings & Wall Murals

The future of radiology departments is not only new imaging equipment, it is focused on becoming more patient-care oriented. A typical MRI or CT Scan can terrify the bravest souls and cause not only unnecessary stress, but delayed recovery times. A relaxed patient is more receptive, therefore they are more responsive to treatment. State-of-the-art, user friendly environments are a positive marketing tool and a strong facility competitive edge. Viewing nature imagery reduces systolic blood pressure and pulse; helps redirect negative thought and decreases boredom.

Acoustic Ceiling Tiles

Our professional-grade acoustical ceiling tiles are engineered for use in highly demanding environments such as hospitals, department stores, corridors, theaters, studios, gyms, offices, schools, MRI rooms, restaurants, basements and places of worship. Our acoustic ceiling tiles can be easily installed into existing drop and suspended ceiling grids. Our acoustic tiles combine exceptional impact and scratch resistance along with the latest sound absorption technology.

Wall Wraps

Our wall wraps are printed on wallpaper-like material, specially formulated for our unique printing process. Our wall wraps can be applied to interior or exterior walls like traditional wallpaper. We also carry vinyl wraps that are applied to walls with no glue or water and can be taken down and reapplied over and over again.

Patient Treatment Room Ceilings & Wall Murals

Installing healing art around intimidating medical equipment creates a positive, patient friendly environment that is an instant sign of caring. Although the patient experience is enhanced through these attractive surroundings, it also helps the caregiver (nurse, technologist, doctor, etc.) because patients are more cooperative. Healing art on dropped ceiling tiles and wall murals in healthcare settings has been proven to lower blood pressure, speed up recovery times and reduce overall stress and anxiety in patients.

Light Lens/Light Diffuser Panels

Our translucent polycarbonate sheet stands out above all other products by providing unmatched light transmittance, glass-like clarity and impact resistance all in one versatile product. Its virtually unbreakable strength allows it to be used in a variety of applications, but can easily be cut to form a custom fit with any common hand tools. It is strongly recommended that full spectrum bulbs are used to illuminate our backlit panels.

Mobile Wraps

Everyday your mobile MRI or medical coach is seen by thousands of people and potential clients. With today's technological breakthroughs, you no longer have to settle for cut letters or simple decals. Using digital illustration in combination with photography, MRlequip.com can create brilliant, eye-catching vehicle graphics for your mobile MRI or medical coach. Through our large format digital printing, we can turn any company vehicle into 24 hour advertising. Customers view companies that take pride in their appearance as being more professional, offer a higher quality product and to be more trustworthy. Whether you have an entire fleet of vehicles you need wrapped or just one, contact us today to get your project started. Life expectancy of wrapped vehicles is 7-10 years.

Please contact customer service, 866-529-3526, for more information or to place an order!



Non-Magnetic Prism Glasses

Non-magnetic Prism Glasses allow the patient to see the room outside the bore of the magnet during their scan, allowing a companion or scenic picture to be viewed. A simple and cost effective comfort measure to be used alone, or in conjunction with a sound system for ultimate relaxation.

The new Deluxe Prism Glasses have absolutely no metal so they are a significant improvement because they will work with all scan protocols and will not produce any image artifacts because they are totally inert. Deluxe Prism Glasses may be used in two different positions to reduce claustrophobia, prisms can be slid out of their mount, rotated and slid back in 180° in reverse. This makes it possible for patients to be scanned **Head-First or Feet-First** with an identical view outside the MRI magnet. The Deluxe Prism Glasses include an elastic adjustable band for use with children and a carrying case (carrying case is not MR compatible). Patent Pending.



Mirrors are adjustable



Mirrors are not adjustable

Model	Description	Price
RA-2012	■ Deluxe Prism Glasses, Adjustable Mirrors, ▲ Carrying Case Included, MR Conditional	\$180.00 ea.
RA-1011	▲ Prism Glasses, Mirrors Not Adjustable	\$99.00 ea.

Sleep Mask

MRI Safe eyeshade sleep mask. Made from satin soft polyester, this is the same mask provided by major airlines and hotels. Available in Dark Blue Only (inside is black). 2 elastic straps.

Model	Description	Price
RA-6000	■ 100 Nylon Masks	\$127.00 cs.
RA-6004	■ 100 Cotton Masks	\$160.00 cs.
RA-6005	■ 1000 Cotton Masks	\$1,550.00 cs.



Goggles

Plastic goggles with interchangeable prescription eyeglass inserts in a wide range of Rx powers. Provides your patients with simple corrective eye wear while in the MRI environment. You simply interchange the Rx inserts based on the Rx needs of your patient. Rx same for both eyes left and right)

Kit Includes:

Two goggles
24 lenses ranging from +6 to -6 in ½ diopter increments



Model	Description	Price
GL-1001	Goggle Kit	\$1,358.00 ea.
GL-1000	Goggle Only, 1 each	\$50.90 ea.



Comfort Pillow Positioners

This innovative pillow was designed to gently elevate and support the head, thus relieving stress on the shoulders during diagnostic MRI and CT procedures where the patient's arms are held over head. Because it is constructed of high-density foam, the pillow provides excellent patient support. Two hand straps allow the patient to hold-on

during the procedure, thus increasing their security. As a result of unique design, the patient's arms can be kept in a straight and more comfortable position, for safe contrast injections.


25" Long x 16" Wide x 8¼" High



Model	Description	Price
FO-1009	■ Flat Bottom	\$130.00 ea.
FO-1010	■ Contoured Bottom	\$140.00 ea.

Leg Rest with Knee Cutouts

Provides leg and knee support while maintaining hip flexion. Helps relieve lower back stress and pain. Stain resistant scancoating is easily cleaned.



Model	Description	Price
FO-1033	■ Small 5" x 18" x 15" w/o scan coating	\$54.60 ea.
FO-1034	■ Medium 7" x 13" x 17.5" w/o scan coating	\$70.00 ea.
FO-1035	■ Large 10" x 24" x 18" w/o scan coating	\$94.00 ea.
FO-1090	■ Small 5" x 18" x 15", w/ scan coating	\$152.30 ea.
FO-1091	■ Medium 7" x 13" x 17.5", w/ scan coating	\$200.30 ea.
FO-1092	■ Large 10" x 24" x 18", w/ scan coating	\$233.54 ea.

Radiolucent Abdominal Immobilizer for MRI & CT

14"W hook and loop -secured band wraps around patient and narrow, contoured tabletop. Patient is secured onto the table, eliminating motion and ensuring clearer image results. Attractive, washable, sterilizable gray canvas with yellow and black hook and loop straps. One-Piece system wraps around table without side rails. Use Two-Piece system for tables with side rails.



Model	Description	Price
FO-1118	1 Piece 75" Long	\$162.00 ea.
FO-1119	1 Piece 87" Long	\$206.55 ea.
FO-1120	2 Piece 36" each side	\$162.00 ea.
FO-1121	2 Piece 60" each side	\$206.55 ea.

Converta-Curve Table Pad for MRI & CT

"Converts" a standard, radiused CT table top into a rigid, flat table top
Ideal for radiation therapy planning and for imaging Bariatric patients. Firm foam offers little compression under heavy loads. Heavy duty vinyl cover is durable and cleanable. Surface of the pad and is slick to allow easy movement of patients. Can be used in CT and MRI tables



Model	Description	Price
FO-3000	Siemens Somatom Tables	\$918.50 ea.
FO-3001	Toshiba Aquilion Tables	\$918.50 ea.
FO-3002	GE Lightspeed Tables	\$918.50 ea.
FO-3003	Phillips Brilliance Tables	\$918.50 ea.



Sand Bag Positioners

Colored outer bag is nylon reinforced vinyl, and dielectrically sealed to prevent leakage or dust. Inner bag is electrically sealed for added protection. Covers are nylon reinforced, color coded vinyl, labeled with weight. Complete with convenient handle.

Fig	Model	Description	Price
A	FO-1105	3 lb, Yellow, 9½" x 4½"	\$24.50 ea.
B	FO-1106	5 lb, Orange, 12" x 5½"	\$26.95 ea.
C	FO-1107	7 lb, Red, 15½" x 6½"	\$29.45 ea.
D	FO-1108	10 lb, Blue, 17" x 6½"	\$31.45 ea.



Artifact Free Sandbags

Will not affect quality of scans. Vinyl double wall construction. Strong, durable and leak-resistant. Nonmagnetic, radiolucent and artifact free.

Fig	Model	Description	Price
E	FO-1113	3½ lb, Black, 9" x 7"	\$60.45 ea.
F	FO-1114	7½ lb, Black, 17" x 5"	\$64.15 ea.
G	FO-1115	11½ lb, Black, 21" x 6"	\$67.80 ea.



Patient Comfort System General 12 Piece Positioner Kit

The Patient Comfort System was developed over a 2-year period in close cooperation with leading imaging centers to maximize patient comfort as well as conform to the specific requirements for compatibility and safety. This viscoelastic material used for medical imaging has been specifically formulated to maximize patient comfort and maintain a thin form factor while optimizing the available free space for the patient when inside the MRI bore.

These mattresses and positioners are compatible with existing tables and are designed to simply replace the current pads. The basic system consists of four rectangular mattress pads covered with advanced breathable space-age waterproof material which resists staining and makes cleaning easier. A knee wedge and neck roll are also included to maximize patient comfort. To increase safety, all pads have an anti-slip undersurface to keep the pads in place while patients are getting on and off the table.



Model	Description	Price
FO-5000	General Positioning Kit	\$1,291.70 kit

Wide Universal Armboard

Accommodates both arm and supplies. Slippery bottom allows for easy insertion under patient's shoulder with minimal patient movement. Patient's weight holds it in place. Cleanable with any type of disinfectant. Hanging hole for easy storage off the floor. ¾" thick. 31½"L. 11½"W (tapers to 6").

Model	Description	Price
FO-1122	Wide Universal Armboard	\$236.93 ea.



For a **LARGE** selection of foam positioning, visit www.MRIequip.com.

■ MR Safe ▲ MR Conditional ● Not MR Safe



Micro Audio System

Equipped for both patient and operator comfort, the Full Stereo system provides everything you need to relax your patients with music during their MRI scan.

Includes:

A powerful compact AM/FM CD system. This system includes 2-way bass reflex speakers with 300 watts from amp to interface box, 14 watts per channel, 5.1 inch woofers and 1.15 inch soft-dome tweeter which add up to superior acoustic performance. MP3 compatible
Stereo Transducer
Separate Patient Volume Control
MRI non-magnetic 25 db noise attenuation headset (fits into most head coils)
Starter package of sanitary headset covers
Patient sound tubing



Warranty
1 Year

Model	Description	Price
RA-1700 ▲	Full Stereo System	\$2,895.00 ea.
RA-1701 ▲	Microphone Option	\$199.00 ea.

Stereo Build Kit

This Kit can be used with an existing stereo system by utilizing the "Line Out" feature (most systems have this).

Push-to-Talk Desktop Bogen® Microphone (Technologist Can Speak to the Patient Headphone)
Integrated Amplifier - Controls the Music Levels for the Patient
iPod® and MP3-Player Adapter Enables Patients to Hear their iPod® Music During the Scan
2 Plastic Pneumatic Headphones for Use in all MRI Coils and all MRI Scanning
1 Bag of Eartips - This Initial Amount will Accommodate Approximately 100 Patients
Custom Cut Cables* - All Cables are Labeled for an Easy Out-of-the-Box Installation
Custom Installation Manual with Graphic Illustrations and Photographs Making it quite easy to install
Plastic Headphone Hanger Clips Keep the Plastic Headphones Handy for the Next Patient Scan



Warranty
1 Year

Model	Description	Price
RA-2005 ▲	Stereo Build Kit	\$3,284.00 ea.

I-Sonic Stereo System

I-Sonic Stereo System with 4 built in speakers-AM/FM HD Radio
Push-to-Talk Desktop Bogen® Microphone (Technologist Can Speak to the Patient Headphone)
iPod and MP3-Player Adapter Enables Patients to hear their iPod music during the scan
The FM Radio is HD Radio Compatible, it can pick up HD Radio Stations with exceptional clarity.
2 Plastic Pneumatic Headphones for Use in all MRI Coils and all MRI Scanning
1 Bag of Eartips - This Initial Amount will Accommodate Approximately 100 Patients
Custom Cut Cables* - All Cables are Labeled for an Easy Out-of-the-Box Installation
Custom Installation Manual with Graphic Illustrations and Photographs Making it Quite Easy to Install
Plastic Headphone Hanger Clips Keep the Plastic Headphones Handy for the Next Patient Scan



Warranty
1 Year

Model	Description	Price
RA-2001 ▲	I-Sonic Stereo System	\$3,695.00 ea.
RA-2010 ▲	CD Player Option	\$199.00 ea.



RA-2010



SI MRI Stereo Sound System

Features: Easy-to-use integrated controls, Amazing sound quality, Sleek appearance, and Headphones equipped with comfort-soft foam cushions

QUALITY MULTIMEDIA SYSTEM:

The SI MRI Stereo Sound System includes speakers, media player, and headphones. The advanced multimedia PC speakers let you enjoy a deeper bass with minimal distortion. These speakers have several plug-in options; they can plug into a computer system or connect to your MP3 player or CD player. The SI MRI Stereo Sound System media player, which plays both CDs and DVDs, also has an iPod dock for your multimedia content. It connects via Component or Composite video, and has a coaxial digital audio output as well.

The media player supports universal docking brackets, is equipped with an on-screen display, and includes a remote control that contains all of the basic iPod commands. The dock plays and charges your iPod at the same time.

The SI MRI Stereo Sound System Ear Muffs are fully dielectric with an extra padded headband for a comfortable fit. The oversized vinyl cushions distribute pressure evenly for user comfort, and can be worn with eye protection or spectacles. Create a more relaxing environment for your patients today!

SPECIFICATIONS:

- DVD player CD player,
- Progressive scan NTSC/PAL On-screen display for iPod and DVD Supports universal docking brackets,
- Basic iPod controls from remote,
- Coaxial digital audio output,
- Remote control Soft eject iPod docking drawer
- Earmuffs offering top-of-the-line protection and comfort. The dielectric design is safe for electrical environments.

Feature:

1. Non-deforming headband
2. Ventilated inner headband
3. Super soft ear cushion
4. Quick Click height adjustment
5. Dielectric
6. 30 NRR
7. 33 dB

- 50 Watt amplifier, Impedance 4/8 ohms, for better sound quality
- Technician stereo speakers feature dynamic, real-time bass equalization to maximize low frequency response while minimizing distortion. Customize your audio balance with easy-to-reach tone and volume controls.
- Push-button desktop microphone, which allows technician talk to patient
- Transducer which transfers audio to patient using twelve feet of Patient Sound Tubing
- Can be upgraded to video system (basic, call customer service for information)

A single set of cables extends from the transducer (attached to the penetration panel's exterior) to the patient head phones, and another set of cables extends between the technologist and the penetration panel. Please choose from the drop down list which lengths of cables are needed, there is no additional charge for different cable lengths.

For an extra charge we can install your system, please call for more information on installation.



Model	Description	Price
RA-3050 ▲	SI MRI Stereo Sound System (Model 1)-with Basic CD/MP3/IPOD Capability, AM and FM Radio with 60 Watt Audio Amplifier, Microphone	\$2,295.00 ea.
RA-3051 ▲	SI MRI Stereo Sound System (Model 2)- Advanced Non Upgradable DVD/ IPOD Player/Pandora/Youtube/Other	\$2,795.00 ea.
RA-3052 ▲	SI MRI Stereo Sound System (Model 3)- Advanced Upgradable for use with the SAMM System-DVD/IPOD Player/Pandora/Youtube/MP3/Other	\$3,795.00 ea.



These headsets have
this style plug...



Other types and
Brands available,
please call



Sanitary Head Set Covers

Large fits models RA-1013, RA-1001, RA-1002, RA-2003 and RA-1057, Small fits model RA-1003
1000 per box

Model	Description	Price
RA-1007	■ Small, White - 2-3" Diameter	\$125.00 box
RA-1006	■ Large, White - 5-6" Diameter	\$125.00 box
RA-1058	■ Large, Black - 5-6" Diameter	\$125.00 box

The Noise Guard Headset

Provides a full 29 db of noise attenuation. The large profile of this headset does not allow it to fit in some head coils.

Model	Description	Price
RA-1001	■ Noise Guard	\$386.00 ea.
RA-1027	■ Replacement Head Band	\$38.95 ea.

The Slimline Noise Guard Headset

Provides 25 db of noise attenuation. The slim profile of this headset fits into most head coils.

Model	Description	Price
RA-1002	■ Slimline Noise Guard	\$299.00 ea.
RA-1027	■ Replacement Head Band	\$38.95 ea.

3 Prong Headset Adapter

For sound systems requiring an alternate connection to our standard two prong male connection, this adapter will work with three prong systems.
Must be ordered with headset and goes with RA-1002 or RA-1004.

Model	Description	Price
RA-1028	Screw Top	\$141.65 ea.
RA-1029	Push to Lock	\$141.65 ea.

Headphones

Cup the ear and are invisible to Scanning. Inner ear tips fit into the ear, providing direct sound. Provides approximately 55 db of noise reduction.

Model	Description	Price
RA-2003	■ Headphones	\$395.00 ea.
RA-2011	■ Replacement Headband	\$37.50 ea.

Siemens Headphone Adaptor

Model	Description	Price
RA-1017	■ Adaptor	\$44.80 ea.



The Non-Headset

Designed for the tightest of head coils. It has two 1" foam ear muffs which can be positioned right up to the patients ear in the head coil, compressing down to 1/2". These headsets are not rated at this time.

Model	Description	Price
RA-1004	■ Non-Headset	\$139.00 ea.

The Foam Ear Piece Headset

Fits easily into any head coil on adults and children alike, and is completely pneumatic.

Model	Description	Price
RA-1005	■ Foam Ear Piece Headset	\$15.00 ea.
RA-1026	■ Case of 10 Foam Ear Piece Headsets	\$120.00 cs.
RA-2020	■ Replacement Foam Ear Pieces for RA-1005, RA-1026, 200 Per Box	\$55.00 bx.
RA-1008	■ Replacement Foam Ear Pieces for RA-1005, RA-1026, 1000 Per Box	\$125.00 bx.



Relaxation and Success Imagery Audio Wellness CD

Using The Power of The Mind-Body Connection, 31 Minutes Long Relaxation and Success Imagery is designed to help teens and adults obtain the relaxation response. Learning to relax is a skill that results in a variety of health benefits. You will be guided through a progression of suggestions to relax all the muscles in your body. Once relaxed you will be asked to repeat positive self-statements and visualize yourself being successful in any goal behavior of your choice. This type of imagery is a form of mental rehearsal and is a proven technique to enhance the attainment of behavior you imagine. These research-based techniques improve learning, enhance athletic performance, promote healing and reduce the negative effects of stress. Relaxation and imagery provide the foundation for positive change and relief from stress. The script, accompanied by tranquil acoustic guitar, offers a choice of male voice on track 1 (15:16) or female voice on track 2 (15:53).



Model	Description	Price
SF-1014	Relaxation Audio Wellness CD	\$99.95 ea.



Non-MRI Equipment

For those situations where non-magnetic items are not required

Call and request your free catalog today

Local: 1-218-828-7788
Toll Free: 1-800-328-5343
Fax Free: 1-866-829-9744
www.ocelco.com

Equipment Catalog

There are over 1,000 new items in this catalog. Products that range from ambulatory aids, beds, bedroom accessories, wheelchairs, scales, diagnostic and respiratory equipment to bariatric products, and much more. From such manufacturers as Drive, Invacare, Sunrise Medical, Lumex plus many more.

Parts & Accessories Catalog

There are a wide variety of parts and accessories in this catalog. Parts that range from bearings, armrest, upholstery, brakes, ramps, pads and cushions, casters, geri-chairs, recliners, patient lifts, bed and safety rails to grab bars. Part breakdowns for patient lifts and geri-chairs. Parts for Drive, Invacare, Everest & Jennings, Gendron, Quickie, Hoyer, Theradyne and many more.

■ MR Safe ▲ MR Conditional ● Not MR Safe





Neonatal Noise Guards

A must have for premie and newborn patients, Neonatal Noise Guards protect these little one's sensitive ears in the MRI environment, and provide peace and quiet in the NICU. Neonatal Noise Guards reduce the surrounding environment noise level by at least 7db, with a reduction in sound pressure of more than 50%. Disposable, latex free.

Box includes 6 pairs

Model	Description	Price
RA-1016	■ Neonatal Guards (6 pairs per box)	\$75.00 box



Reusable Earplugs

Case includes 50 each

Easy grip handle allows easy insertion.

Air bubble in the tip provides cushioned comfort.

Can be washed and reused.

Each pair comes with a carrying case to store plugs between uses.

Slips easily into the pocket.

Bright blue detectable model is ideal for food processing use.

Noise Reduction Rating 26

Model	Description	Price
RA-4000	■ Reusable Earplugs (case of 50)	\$64.30 cs.



Swirl Earplugs

Fun colored to spark up the hearing protection program.


Different colors in every package.

Tapered for easy insertion in any size ear canal.

Non-irritating, smooth surface for added comfort.

Noise Reduction Rating 32

Model	Description	Price
RA-1009	■ Without cord, 1000 Count / 500 Pair	\$160.00 cs.
RA-4006	■ With cord, 200 Count / 100 Pair	\$66.32 cs.



Pure Fit - Long Earplugs

Longer length for easy fitting and removal.

Bright green color for easy compliance checks.

Tapered for easy insertion in any size ear canal.

Non-irritating, smooth surface for added comfort.

Noise Reduction Rating 32

Model	Description	Price
RA-4004	■ Without cord, 1000 Count / 500 Pair	\$170.00 cs.
RA-4005	■ With cord, 200 Count / 100 Pair	\$71.30 cs.



Plug Dispenser

Clicking mechanism for easier use and dispensing. Chute delivers plugs right into the palm of your hand. Affordable enough to be placed near every work area. Eliminates the mess of individual paper/polybag packaging. Each dispenser is sealed and completely recyclable. Available in 250 or 500 count per Plug Dispenser with either Swirl Plugs or Pure Fit foam earplugs. Mounting bracket with hardware included with every Plug Dispenser.



Containers are Non-Refillable

Model	Description	Price
RA-4009	▲ Swirl, 250 Count / 125 Pair	\$53.30 ea.
RA-4010	▲ Swirl, 500 Count / 250 Pair	\$93.25 ea.
RA-4011	▲ Pure Fit, 250 Count / 125 Pair	\$66.32 ea.
RA-4012	▲ Pure Fit, 500 Count / 250 Pair	\$102.24 ea.

Ear Muffs

For sound protection.

Model	Description	Price
RA-1057	■ 33 db Noise Attenuation	\$99.50 ea.
RA-1013	■ 25 db Noise Attenuation	\$126.20 ea.
RA-1027	■ Replacement Head Band	\$38.95 ea.



RA-1057



RA-1013



RA-1027

Large Sanitary Head Set Covers

Fits models RA-1013, RA-1001, RA-1002, RA-2003 and RA-1057
1000 per box

Model	Description	Price
RA-1006	■ White (1000 per box)	\$125.00 box
RA-1058	■ Black (1000 per box)	\$125.00 box



Hand Sanitizer in Green Case with Plastic Clip

Features: Hand sanitizer contains 62% Ethyl Alcohol as an active ingredient/antiseptic. Black heavy duty clip attaches to belt loop, backpack or beach bag. Handy neoprene sleeve holds sanitizer bottle securely. Purell is as safe on hands as regular soap and water and you can use it as often as you need. Purell Instant Hand Sanitizer is a refreshing alcohol-based instant hand sanitizer that kills 99.99% of most common germs that may cause illness in as little as 15 seconds. This .5 ounce hand sanitizer from Purell® with our added leash is ideal for everyday use. Kills the germs that may make you and your family sick, no matter where you are!

Model	Description	Price
MT-3020	Hand Sanitizer (box of 125)	\$470.48 box





Pre-Filled Fire Extinguisher

The Model FE-2000 fire extinguisher is designed to be the ideal fire extinguisher for hazard areas that require a reliable and safe non-magnetic fire extinguisher. The agent is proven to be electrically non-conductive (will not damage sensitive electrical equipment like computers), leaves no residue, and requires no clean-up. The agent has no ozone depletion potential and is lower in toxicity compared to comparable extinguishing agents.

The UL/ULC approved fire extinguisher is completely non-magnetic stainless steel shell and non-magnetic valve, hose and nozzle. The unit was tested by Magnetic Resonance Safety Testing Services and shown to be non-magnetic to 3.0 Tesla. The extinguisher can safely be mounted within the MRI room.

Features

- Pre-Filled for A, B and C fires
- No Clean Up - No Residue
- Non-Magnetic (tested up to 3 Tesla)
- Safe for sensitive electronic equipment
- Safe for exposure to occupants
- Zero ozone depletion
- Six-year limited warranty from date of delivery to original end-user purchaser
- FAA approved
- Rechargeable
- UL/ULC Listed
- Meets or exceeds requirements of ANSI/UL 2129 and 711
- Temperature range is -40° F to 120° F (-40° C to 49° C)
- Mounting bracket included
- MSDS included in shipment

Model	Description	Price
FE-2000	English Labeling	\$870.00 ea.
FE-2001	Canadian Labeling (English & French)	\$870.00 ea.

More specifications available upon request or online at www.MRIequip.com



Custom Wall Mount Fire Extinguisher Cabinets

You determine the: Size, Finish, Color, Text, Window, Window Size, Window Placement, Locking or Non-Locking, Recessed, Surface Mounted, etc.

Please call for more information, we would be happy to do a quote for you today. Over 45 years of unsurpassed quality manufacturing for the building construction industry.

The most complete product line in the industry continues to grow through constant innovation — fire-rated cabinets, access panels, detention cabinets and a full line of fire hose and fire department valve cabinets.

Active participation in industry professional organizations such as Construction Specifications Institute, Door and Hardware Institute, American Institute of Architects, Fire Equipment Manufacturers' Association and National Association of Fire Equipment Distributors.

Model	Description	Price
FE-3000-CU ●	Fire Extinguisher Cabinets	Call for Pricing

For use outside MR Room

NOT MR SAFE

This product is magnetic



De-Ionized Water Fire Extinguisher

The ultimate extinguisher for class A fires, especially where a potential Class C (electrical) hazard exists. The fine spray from the unique misting nozzle provides safety from electrical shock, greatly enhances the cooling and soaking characteristics of the agent and reduces scattering of the burning materials. Unequivocally the best extinguisher for hospital environments, telecommunications facilities and "clean room" manufacturing facilities.

User Friendly :

Wand provides greatest operator safety
Maximum visibility during discharge
No electrical conductivity back to the operator
No thermal or static shock
Minimal training - simplest to use, least expensive agent available

Clean :

No toxicity
No respiratory concerns
De-ionized water extinguishing agent
Exclusive crevice free, butt welded cylinder
Antiseptic white paint finish (easy to clean interior & exterior)
No powder residue

Specifications:

	FE-1004	FE-1000
Capacity	1¾ Gallon	2½ Gallon
Application	Hose & Wand	Hose & Wand
U/L Rating	2A:C	2A:C
Extinguishing Height	18½"	24½"
Height with Hose & Wand, Installed in Clips	25"	27½"
Width with Hose & Wand, Installed in Clips	14"	14"
Depth (diameter)	7"	7"
Range (feet)	10'-12'	10'-12'
Discharge Time	72 Seconds	80 Seconds
Standard Bracket	Wall	Wall

Complete Kits

Includes:

2 Bottles De-Ionized Water
1 Fill Adapter

Model	Description	Price
FE-1006* ▲	1¾ Gallon Fire Extinguisher	\$244.00 ea.
FE-1003* ▲	2½ Gallon Fire Extinguisher	\$250.65 ea.

* Filling requires the use of a nitrogen pressurizing source set at no more than 125 PSI. The Fill Adapter is required for the refilling process and replaces the wand during such process. This will be needed by most professionals to handle the refilling process.

Extinguishers are not shipped filled.

Contact a local filling service near you.

Replacement Components

Model	Description	Price
FE-1004 ▲	1¾ Gallon Fire Extinguisher	\$200.35 ea.
FE-1000 ▲	2½ Gallon Fire Extinguisher	\$209.70 ea.
FE-1002	Fill Adapter	\$34.95 ea.
FE-1005 ■	Two 1¾ Gallon Bottles De-Ionized Water	\$25.20 ea.
FE-1001 ■	Two 2½ Gallon Bottles De-Ionized Water	\$28.20 ea.





NOT MR SAFE

This product is magnetic

Due to the strength of magnet, only regular ground shipping is available



NOT MR SAFE

This product is magnetic

Screening Magnet

Powerful Rare-Earth Handheld Magnet for Magnetic Resonance Screening/Sites:

Magnetic Resonance (MR) professionals must be able to screen patients, devices, and implants for possible ferromagnetic properties prior to allowing them to enter clinical and research magnetic resonance imaging (MRI) environments and magnet rooms (so-called, "Zone III and IV" regions). It has also been proposed that all MR sites have ready access to powerful, more than 1,000 gauss hand held magnets that they can use for such testing purposes. This magnet has been specially designed by experienced clinical and research MR practitioners and leaders in the field of MR safety to meet the needs of busy MR sites. The more than 1,200 gauss magnetic field of the magnet exerts a sufficiently powerful magnetic field to sensitively test even weakly ferromagnetic objects. The magnet itself is also composed of a neodymium - a so-called rare earth magnet - that permits it to lose its strength extremely slowly, at a rate of a 1% loss every 100 years. The recessed magnet in custom foam-fitted housing, the form fitted handle, and the surface ball bearing mounts allow for ready removal of the magnet - by anyone - should it become inadvertently attached to large ferromagnetic objects, such as oxygen tanks, wheelchairs, etc. The bright red anodized coloring of the magnet housing contrasts well with the silver labeling clearly stating that the magnet is not MR safe and should not be brought into the MR magnet room/Zone IV regions. This magnet was developed in part with the help of Dr. Emanuel Kanal, Department of Radiology, University of Pittsburgh Medical School, and is field tested to provide years of reliable MR safety and compatibility testing under typical clinical and research MRI settings. Also optionally available is a foam padded carrying and storage case to protect the magnet between uses.

Specifications:

Weight: 2 lbs

Color: Red, Silver lettering

Magnetic Field Strength: 1,200 gauss

Dimensions: 3½" height x 3" width x 5½" length

Model	Description	Price
MG-1000	● Screening Magnet	\$370.00 ea.
MG-1001	● Carrying Case	\$75.60 ea.

SuperWand

The SuperWand's exclusive 360° detection field provides uniform sensitivity and tip pinpointing to detect weapons and other metal objects with extreme accuracy. From U.S. airports to sporting events, the SuperWand is the hand-held of choice for today's security screener.

Accurate detection of all ferrous, non-ferrous and stainless steel and other metallic objects

Elongated profile offers 360° detection coverage plus tip pinpointing for easy scanning from head to foot

Easy and convenient to use with one-touch button and three-color LED lights

Clear audible or Silent / Vibrate LED alarm signal

Large 9.5" scan area allows thorough scanning with minimum passes

Ergonomically designed grip fits comfortably in virtually every hand

9V battery (included) operation

Model	Description	Price
MG-4000	● SuperWand	\$199.00 ea.



NOT MR SAFE

This product is magnetic



Reference Manual for Magnetic Safety, Implants and Devices

The Reference Manual for Magnetic Resonance Safety, Implants and Devices: 2011 Edition (ISBN-10, 0-9746410-7-3; ISBN-13, 978-0-9746410-7-2; Biomedical Research Publishing Group) is an indispensable textbook for radiologists, MRI technologists, and facility managers. This internationally acclaimed series includes updated guidelines and recommendations from the latest information in the peer-reviewed literature as well as documents developed by the International Society for Magnetic Resonance in Medicine (ISMRM), the American College of Radiology (ACR), the Food and Drug Administration (FDA), the National Electrical Manufacturers Association (NEMA), the International Electrotechnical Commission (IEC), the Medical Devices Agency (MDA), and the Institute for Magnetic Resonance, Safety, Education and Research (IMRSE). Features of the 2011 Edition include patient screening forms in English and Spanish and guidelines for scanning patients with electronically-activated devices.


-Section I presents fully updated safety guidelines and recommendations.

-Section II has the latest information for implants, devices, and materials tested for the MRI environment.

-Updated information pertaining to MRI Contrast Agents and Nephrogenic Systemic Fibrosis (NSF) is included along with websites to monitor emerging findings.

-An Appendix provides website information for more than 220 biomedical companies to facilitate finding the latest company information for implants and devices.

Each year a new edition comes out, upon ordering you will receive the most current addition available.



Model	Description	Price
MT-1108	2011 Edition Safety Manual	\$99.00 ea.

GaussAlert Magnetic Field Strength Alarm System


GaussAlert is designed to help keep MR Conditional equipment outside of the MRI exclusion zone. GaussAlert is programmed to alarm when the preset magnetic field strength is exceeded. It produces an attention-getting audio alarm when MR conditional equipment is too close to the MRI magnet.

Medical equipment that is used in the MRI room, such as infusion pumps, contrast injectors, patient monitors, anesthesia machines and others, may malfunction and potentially injure the patient if not kept outside of a specific magnetic field strength. In addition, MR Conditional equipment can degrade the image quality of a scan if brought too close to the MRI magnet. With GaussAlert, the possibility of equipment malfunction and image degradation is dramatically reduced.

GaussAlert can also be useful as an alert for MR unsafe equipment, such as a crash cart or an emergency resuscitator, which can be kept near, but should not be taken into the MRI magnet room.

The GaussAlert is available in three factory preset gauss thresholds: 10 (MG-2504), 30 (MG-2505) and 100 (MG-2506) gauss, based on the equipment manufacturer's specifications. GaussAlert easily and quickly attaches to any piece of equipment and requires no periodic maintenance for up to five years. It is internally-powered and continuously monitors the ambient magnetic field strength.

GaussAlert compliments the function of our FerrAlert Halo Entry detector by monitoring equipment that has passed through the portal but still must be monitored inside the magnet room.




Model	Description	Price
MG-2504	10 Gauss Threshold	\$950.00 ea.
MG-2505	30 Gauss Threshold	\$950.00 ea.
MG-2506	100 Gauss Threshold	\$950.00 ea.



Pre-screening and Entry Ferromagnetic Detection Systems

Entry System



FerrAlert™ HALO ENTRY – Ferromagnetic Detector complies with ACR and JCAHO standards to protect the magnet room against moderate to large size ferromagnetic projectiles. The device can be installed freestanding at the approach to the MR magnet room, or unobtrusively onto the door frame outside or inside the magnet room, depending on the swing of the door. This product significantly helps to protect the investment in the magnet and to avoid a catastrophic accident. FerrAlert™ HALO ENTRY has been certified by Siemens Medical Solutions for installation and use inside the magnet room for both 1.5T and 3.0T MRI magnets.

OPERATION: When the GREEN HALO array of READY indicator lights, located in the vertical columns of the portal, are illuminated (ON), the system is working. If a ferromagnetic object is detected, an AUDIO ALARM will sound at the portal (and operator's display module, if required). The GREEN READY HALO indicator lights will extinguish and AMBER LOCATION indicator lights will illuminate. The location of the ferromagnetic object in the Sensor Portal is indicated by an illumination of the group of amber LEDs closest to the ferromagnetic object. For example, if the bottom left amber LEDs are activated, the ferromagnetic object is near the floor on the left side. The NIGHT ALARM switch on the optional OPERATOR'S DISPLAY MODULE should be turned ON when there are no MR personnel present.

Pre-screening System



FerrAlert™ HALO PRESCREEN is a highly sensitive patient and personnel screener that detects VERY SMALL ferrous objects down to the size of a typical hairpin. It is up to 25 times more sensitive than our FerrAlert™ HALO ENTRY detector. The recommended placement for this device is in the patients' screening or locker-room area, away from the magnet room. This product helps increase patient throughput and avoid potential misdiagnosis. It is also very useful in screening staff from other departments who are visiting the MRI area.

OPERATION: When the GREEN HALO array of READY indicator lights, located in the vertical columns of the portal are illuminated (ON), the system is working. If a ferromagnetic object is detected, an AUDIO ALARM will sound, the GREEN READY HALO indicator lights will extinguish and AMBER LOCATION indicator lights will illuminate. The location of the ferromagnetic object in the Sensor Portal is indicated by an illumination of the group of amber LEDs closest to the ferromagnetic object. For example, if the top left amber LEDs are activated, the ferromagnetic object is near the left side of the portal at head level.

Solo Detector



Ferromagnetic screener for patients and personnel - The latest advancement in Pre-Screening Technology for MRI Safety. **OPERATION** APPROACH FERRALERT SOLO DETECTOR. When close to the detector, 6 green READY TO SCAN lights in the array will turn on. MAKE A COMPLETE ROTATION. If any ferromagnetic hazard is detected, the green READY TO SCAN lights will turn off. An audio alert will sound briefly and amber LOCATION DETECTION lights will flash on the array at the location of the detected ferromagnetic hazard. Additionally, an indicator at the top of the array aids the screening supervisor to see the location of the ferromagnetic hazard. IF THE GREEN READY TO SCAN LIGHTS STAY LIT, NO FERROMAGNETIC HAZARDS WERE DETECTED. After moving away, the detector will return to standby mode. The purple POWER indicator will be illuminated. The detector constantly monitors the ambient magnetic field and the blue bar graph at the top of the array displays the degree of magnetic disturbance.

Model	Description	Price
MG-2500	Entry System	Call for Pricing
MG-2501	Pre-Screening System	Call for Pricing
MG-2507	Solo Detector	Call for Pricing

Visit www.MRIequip.com
 for more specifications, or call for complete brochure.

See Page 1 and 2 for Non-Ferromagnetic Stretcher and Wheelchair



Caution, Restricted to Screened MRI Patients & MRI Personnel

Black / Yellow, 4" x 11"

Model	Description	Price
MT-1166	Sticker	\$16.50 ea.
MT-1165	Plastic Sign	\$19.95 ea.



Danger, Must Be Accompanied by Trained MRI Personnel

Black, Red and White, 4" x 11"

Model	Description	Price
MT-1168	Sticker	\$16.50 ea.
MT-1167	Plastic Sign	\$19.95 ea.



Notice, Approaching MRI Restricted Area Sign

Blue and White, 4" x 11"

Model	Description	Price
MT-1170	Sticker	\$16.50 ea.
MT-1169	Plastic Sign	\$19.95 ea.



Yellow Insert Signs

Thumb notch on pocket allows for easy removal of inserts (inserts sold separately)

Brass grommets for easy mounting

Can hold one header and multiple inserts of your choice

Header Insert - 3½" x 8"

Inserts - 1¾" x 8"

Model	Description	Price
MT-1180	3 Pocket, 12½" x 8"	\$38.40 ea.
MT-1181	5 Pocket, 16½" x 8"	\$49.30 ea.



Insert Sign Headers

3½" x 8"


Model	Description	Price
MT-1182	CAUTION	\$4.70 ea.
MT-1183	DANGER	\$4.70 ea.




Insert Sign Inserts

1¾" x 8"

Model	Description	Price
MT-1184	"High Radiation Area"	\$3.70 ea.
MT-1185	"No Entry by Unauthorized Personnel"	\$3.70 ea.
MT-1186	"Radiation Area"	\$3.70 ea.
MT-1187	"Restricted Area"	\$3.70 ea.
MT-1188	Customize your own	\$12.50 ea.





MRI Zone Sign, Zone 1

10" x 14"

Aluminum

General Public - White with Black lettering

Model	Description	Price
MT-1152	■ Zone 1	\$30.00 ea.



MRI Zone Sign, Zone 2

10" x 14"

Rigid Plastic

Unscreened MRI Patients - White/Blue and Black lettering

Model	Description	Price
MT-1153	■ Zone 2	\$28.25 ea.



MRI Zone Sign, Zone 3

10" x 14"

Rigid Plastic

Caution - Screened MRI Patients and Personnel Only

Yellow with Black lettering

Model	Description	Price
MT-1154	■ Zone 3	\$28.25 ea.



MRI Zone Sign, Zone 4

10" x 14"

Rigid Plastic

Danger - Screened MRI Patients Under Constant Direct

Supervision of Trained MRI Personnel Only

White/Black/Red with White and Black lettering

Model	Description	Price
MT-1155	■ Zone 4	\$33.75 ea.



MRI Zone Sign, Complete Set

Set of one each of all four Zone Signs

Model	Description	Price
MT-1160	■ Zone Sign Set of 4	\$112.50 set



MRI Zone Sign General Diagram

12" High x 18" Wide

White background with black lettering

Rigid plastic with foam tape for hanging

Model	Description	Price
MT-1161	■ Zone General Diagram	\$42.00 ea.



MRI Zone IV 17" Floor Sticker

Reads: DANGER! MRI Zone IV, Screened MRI Patients under constant direct supervision of trained MRI personnel only
 17" Diameter
 White / Red / Black

Model	Description	Price
MT-1923	Zone Floor Sticker	\$49.75 ea.



8" Round "Danger Restricted Access" Gauss Line Sticker

8" Diameter
 White / Red / Black

Model	Description	Price
MT-2500	5 Gauss Line	\$34.95 ea.
MT-2501	10 Gauss Line	\$34.95 ea.
MT-2502	30 Gauss Line	\$34.95 ea.
MT-2503	100 Gauss Line	\$34.95 ea.



Floor Sticker Warning Sign "Magnet Always On"

Size: 3' x 5'
 Pressure Sensitive Vinyl, Self adhesive
 Laminate is Flextrack, Non-skid

Model	Description	Price
MT-1128	Floor Sticker	\$240.00 ea.



Floor Sticker Warning Sign "Magnet Always On"

Model	Description	Price
MT-1018*	8" Diameter, 2 Color Custom Slip-Guard Floor Sticker	Call for Pricing
MT-1019*	17" Diameter, 2 Color Custom Slip-Guard Floor Sticker	Call for Pricing
MT-1061	17" Diameter Floor Sticker - Magnet Always On	\$56.50 ea. Call for Quantity Discounts.



MT-1061
 17" Diameter Floor Sticker Sign

MT-1018* & MT-1019*

Customizable 8" and 17" Floor Stickers

* Slip-Guard floor signs are available with your message and symbol/logo. When ordering, specify size, message, design layout, symbol/logo, background color, 2 color choices and quantity.





MRI Warning Signs



MT-1025



MT-1026



MT-1033



MT-1034

12" x 18" Signs in Spanish and English

Custom Signs Available - Call for More Details

New codes are requiring facilities to use **Reflective** Warning Stickers and Signs in case of fire or power outage*.

Plastic Signs

Hardware is NOT included

Model	Description	Price
MT-1025	12" x 18" English - Items Not Allowed	\$34.50 ea.
MT-1118	12" x 18" Reflective , English - Items Not Allowed	\$57.20 ea.
MT-1026	12" x 18" Spanish - Items Not Allowed	\$34.50 ea.
MT-1119	12" x 18" Reflective , Spanish - Items Not Allowed	\$57.20 ea.
MT-1033	12" x 18" English - No Unauthorized Entry	\$34.50 ea.
MT-1120	12" x 18" Reflective , English - No Unauthorized Entry	\$57.20 ea.
MT-1034	12" x 18" Spanish - No Unauthorized Entry	\$34.50 ea.
MT-1121	12" x 18" Reflective , Spanish - No Unauthorized Entry	\$57.20 ea.
MT-1041	4" x 8" English - Magnet Is Always On	\$19.50 ea.
MT-1114	4" x 8" Reflective , English - Magnet Is Always On	\$26.25 ea.
MT-1042	4" x 8" Spanish - Magnet Is Always On	\$19.50 ea.
MT-1115	4" x 8" Reflective , Spanish - Magnet Is Always On	\$26.25 ea.

Warning Stickers

Model	Description	Price
MT-1043	4" x 8" English - Magnet Is Always On	\$15.50 ea.
MT-1116	4" x 8" Reflective , English - Magnet Is Always On	\$24.75 ea.
MT-1044	4" x 8" Spanish - Magnet Is Always On	\$15.50 ea.
MT-1117	4" x 8" Reflective , Spanish - Magnet Is Always On	\$24.75 ea.



MT-1041 & MT-1043

4" x 8" Stickers and Signs in Spanish and English

*Codes are effective in varying states, contact your local fire dept. for more information.

■ MR Safe ▲ MR Conditional ● Not MR Safe



MRI Carpeted Floor Mat and Custom Carpet

MRI floor mat - Magnet is Always On, as shown in picture to the right. Made of indoor material. Carpeted floor mats are now available with your message and symbol/logo. When ordering, specify message, background color, if for indoor or indoor/outdoor use and symbol/logo to be shown on carpet.

Model	Description	Price
MT-1045	3' x 5' Carpeted Floor Mat - Magnet Always On	\$195.00 ea.
MT-1122	Custom Carpeted Floor Mat	Call for Pricing

Indoor Outdoor Carpeted Floor Mat

MRI floor mat - Magnet is Always On, as shown in picture to the right. Made of indoor/outdoor material. Made per order, please allow 4-6 weeks for delivery.

Model	Description	Price
MT-1123	Indoor/Outdoor Carpeted Floor Mat	Call for Pricing

MRI Warning Carpets

MT-1045
3' x 5' Carpeted Floor Mat



Floor Mat Carpet Warning Sign

"Caution! Entering MRI Restricted Area"

3' x 5' Yellow background with black lettering.

Model	Description	Price
MT-1124	Yellow Floor Mat	\$200.00 ea.



Specialized MRI Catalogs!

Call to request your copy today or
e-mail your request to
cs@mrlequip.com.

Breast Imaging Products



Pediatric Products



Disposable Supplies





- Die Cut
- Permanent Adhesive
- Laminated

MR Safe Stickers

Model	Description	Price
MT-8006-50	■ ½" x ¾", 50 Pack	\$19.75 pk.
MT-8006-100	■ ½" x ¾", 100 Pack	\$34.50 pk.
MT-8001-10	■ 1½" x 2", 10 Pack	\$14.50 pk.
MT-8001-50	■ 1½" x 2", 50 Pack	\$66.50 pk.
MT-8001-100	■ 1½" x 2", 100 Pack	\$124.00 pk.
MT-8000-10	■ 3½" x 4", 10 Pack	\$15.00 pk.
MT-8000-50	■ 3½" x 4", 50 Pack	\$67.50 pk.
MT-8000-100	■ 3½" x 4", 100 Pack	\$125.00 pk.

MR Conditional Stickers

Model	Description	Price
MT-8007-50	■ ½" x ¾", 50 Pack	\$19.75 pk.
MT-8007-100	■ ½" x ¾", 100 Pack	\$34.50 pk.
MT-8003-10	■ 1½" x 2", 10 Pack	\$14.50 pk.
MT-8003-50	■ 1½" x 2", 50 Pack	\$66.50 pk.
MT-8003-100	■ 1½" x 2", 100 Pack	\$124.00 pk.
MT-8002-10	■ 3½" x 4", 10 Pack	\$15.00 pk.
MT-8002-50	■ 3½" x 4", 50 Pack	\$67.50 pk.
MT-8002-100	■ 3½" x 4", 100 Pack	\$125.00 pk.

Not MR Safe Stickers

Model	Description	Price
MT-8008-50	■ ½" x ¾", 50 Pack	\$19.75 pk.
MT-8008-100	■ ½" x ¾", 100 Pack	\$34.50 pk.
MT-8005-10	■ 1½" x 2", 10 Pack	\$14.50 pk.
MT-8005-50	■ 1½" x 2", 50 Pack	\$66.50 pk.
MT-8005-100	■ 1½" x 2", 100 Pack	\$124.00 pk.
MT-8004-10	■ 3½" x 4", 10 Pack	\$15.00 pk.
MT-8004-50	■ 3½" x 4", 50 Pack	\$67.50 pk.
MT-8004-100	■ 3½" x 4", 100 Pack	\$125.00 pk.

Multi-Packs


Multi-packs include many different sizes and styles of stickers.

Model	Description	Price
MT-8010-SM	■ 3 of Each Style, 18 Total	\$25.00 pk.
MT-8010-MD	■ 6 of Each Style, 36 Total	\$50.00 pk.
MT-8010-LG	■ 10 of Each Style, 60 Total	\$70.00 pk.
MT-8010-XL	■ 20 of Each Style, 120 Total	\$135.00 pk.
MT-8009	■ ½" x ¾", 50 of Each Style	\$51.75 pk.

■ MR Safe ▲ MR Conditional ● Not MR Safe



MRI Safe - Do Not Remove From MRI Area - Warning Stickers



- Die Cut
- Permanent Adhesive
- Laminated

Model	Description	Price
MT-1065	■ Red - 5 Large, 5 Medium, 5 Small	\$19.95 pk.
MT-1092	■ Green - 5 Large, 5 Medium, 5 Small	\$19.95 pk.
MT-1066	■ Red - 6 Large, 9 Medium, 15 Small	\$28.50 pk.
MT-1093	■ Green - 6 Large, 9 Medium, 15 Small	\$28.50 pk.

Pack of 10

Pack of 50

Pack of 100

Description	Model	Price	Model	Price	Model	Price
1" x 3" Red	■ MT-1064	\$11.10 pk.	MT-1071	■ \$44.40 pk.	MT-1072	■ \$77.70 pk.
1" x 3" Green	■ MT-1083	\$11.10 pk.	MT-1086	■ \$44.40 pk.	MT-1089	■ \$77.70 pk.
1 1/2" x 6" Red	■ MT-1063	\$13.80 pk.	MT-1069	■ \$47.20 pk.	MT-1070	■ \$80.50 pk.
1 1/2" x 6" Green	■ MT-1084	\$13.80 pk.	MT-1087	■ \$47.20 pk.	MT-1090	■ \$80.50 pk.
4" x 6" Red	■ MT-1062	\$16.65 pk.	MT-1067	■ \$49.95 pk.	MT-1068	■ \$88.80 pk.
4" x 6" Green	■ MT-1085	\$16.65 pk.	MT-1088	■ \$49.95 pk.	MT-1091	■ \$88.80 pk.

DANGER! NOT MRI Safe Warning Stickers

- Die Cut
- Permanent Adhesive
- Laminated



Model	Description	Price
MT-1105	■ Red - 5 Large, 5 Medium, 5 Small	\$19.95 pk.
MT-1106	■ Red - 6 Large, 9 Medium, 15 Small	\$28.50 pk.

Pack of 10

Pack of 50

Pack of 100

Description	Model	Price	Model	Price	Model	Price
1" x 3" Red	■ MT-1099	\$11.10 pk.	MT-1100	■ \$44.40 pk.	MT-1101	■ \$77.70 pk.
1 1/2" x 6" Red	■ MT-1077	\$13.80 pk.	MT-1078	■ \$47.20 pk.	MT-1079	■ \$80.50 pk.
4" x 6" Red	■ MT-1102	\$16.65 pk.	MT-1103	■ \$49.95 pk.	MT-1104	■ \$88.80 pk.



Split, Swinging Stop Sign

Mounts separately from the door. Even with the door open, the sign continues to block the entrance and then can be pushed through separately. Hinges are different than pictured

Model	Description	Price
MT-5000	English	\$1,105.00 ea.
MT-5001	Spanish	\$1,105.00 ea.

For use outside MR Room



NOT MR SAFE

This product is magnetic



For use outside MR Room

Visit
www.mriequip.com
for more specifications.

Lighted Signs

Choose from "MAGNET IN USE", "MAGNET ALWAYS ON" or create a Custom Message.

Sign must be hard wired in and comes with a battery backup.

Size of Sign:

11" W x 9" T x 2 1/4" D

Mounting Bracket:

5 3/4" x 4 1/2" x 1 1/8" Tall

NOT MR SAFE

This product is magnetic

Two Sided and Custom Signs Available - call for pricing.

Model	Description	Price
MT-1006	Red, Magnet In Use	\$199.45 ea.
MT-1008	Red, Magnet Always On	\$199.45 ea.



"Do Not Enter" Floor Cones

Complete kit contains:

Set of 2 Cones, 10' of Chain, Lock in Caution Sign (10" x 10") with Exclamation Point

Model	Description	Price
MT-1612	Complete Kit	\$219.00 set
MT-1605	"Do Not Enter" Floor Cones	\$83.00 ea.
MT-1606	Barrier Chain	\$45.00 ea.
MT-1611	Lock in Sign w/ Exclamation Point	\$47.00 ea.

For more stanchion posts and options please visit www.MRIEquip.com



Plastic Name Badge Clip

Set of 25, of one color.
 Choose yellow, white, red, green or royal blue

Model	Description	Price
BD-2000-CL	■ Name Badge Clip (case of 25)	\$39.95 ea.



Badge ID Holder / Convention Pouch

12 per case
 4" x 2 1/2"
 48" Cord
 2 pockets: one open and one zippered
 Nylon

Model	Description	Price
MT-1803	■ Badge ID Holder (12)	\$42.85 cs.



Lanyard with Breakaway Cord/Bulldog Clip

10 per case
 Have your name custom imprinted on each lanyard with a minimum order of 100 each and \$50.00 set-up fee.

Model	Description	Price
MT-1800	■ Lanyard (case of 100)	\$52.60 cs.



Brass Pin for Name Badges

Model	Description	Price
BD-1002	▲ Brass Pin Only (package of 10)	\$20.00 pk.



One Piece Belt with Safety Buckle

80" belt in 2" black webbing with plastic safety buckle.

Model	Description	Price
GB-1003	■ Safety Belt	\$19.95 ea.





Breast Shield

Today, CT Scans play a critical role in the diagnosis and management of most patients with complex thoracic and cardiopulmonary diseases. As the use of CT imaging increases, however, so does the risk for latent carcinogenic effects of ionizing radiation on radiosensitive tissue including the female breast. The recently redesigned AttenuRad Breast Shield System consists of a thin (1 mm) piece of bismuth impregnated synthetic rubber covering that mounts firmly on an offset quarter inch foam base. The breast shield system offers up to a 57% reduction of radiation exposure during a CT scan without significant changes in image quality, including artifacts and beam hardening effects. The single-piece AttenuRad CT Breast Shield is designed for multiple uses.

Model	Description	Price
MO-3002-MD	■ Medium Breast Shield	\$87.95 ea.
MO-3002-LG	■ Large Breast Shield	\$87.95 ea.



Adult Breast Shield Covers - Case of 10

Fits both MO-3002-MD & MO-3002-LG.

Plastic, Latex Free Covers comes in a Package of 10

Introducing MRlequip.com's new disposable covering for the Adult Breast Shield System. The new plastic covering is cost-effective alternative to the single-use Adult Breast Shield, reducing the risk of contamination while still providing an artifact and streak-free image.

2 ml thick


Model	Description	Price
MO-3014	■ Breast Shield Covers	\$20.00 cs.



CT Breast Shield System

Studies have indicated that infants and children are as much as 10 times more susceptible to carcinogenesis from radiation than adults. Breast tissue exposed to radiation from CT procedures involving the chest or abdomen is an area of particular concern for young females. The CT Breast Shield System provides significant in-plane breast shielding (57%) for pediatric MDCT, without qualitative or quantitative changes in the CT image. The use of the system's offset foam base gives the shields additional substance and durability, making them easy to use. The use of the elevating foam also decreases the patient's exposure to scattered radiation.

Model	Description	Price
MO-3010	■ Neonatal, Box of 10	\$119.05 box
MO-3011	■ 1 Year Old, Box of 8	\$119.05 box
MO-3012	■ 5 Year Old, Box of 6	\$119.05 box
MO-3013	■ 10 Year Old, Box of 4	\$119.05 box





Thyroid Shields - Case of 10

Positioned over the Adam's Apple, the AttenuRad CT Thyroid shield is designed to provide maximum protection to the radio sensitive thyroid gland. To ensure tight positioning, a curved notch in the center of the shield is placed around the Adam's Apple and an adhesive backing holds the shield firm as the patient breathes and swallows. The AttenuRad CT Thyroid shield, when used for the soft tissues of the neck, is cut to a thickness of 1 mm allowing for a 60% dose savings to the Thyroid.

Model	Description	Price
MO-3001	■ Thyroid Shield	\$109.90 cs.



CT MODALITY

Eye Shields - Case of 20

The lens of the eye is sensitive to radiation during routine Cranial CT. Patient undergoing routine CT scans of the head are especially vulnerable to complications. The AttenuRad CT Eye Shield is a single-piece bismuth impregnated synthetic rubber covering designed to provide significant protection. The shield molds over the patients nasal bridge with the shielding portions laying flat over the patients closed eyes. Adhesive material ensures a tight seal with minimal air gaps. This shield is cut to a thickness of 1 mm, allowing for 50% dose savings to the eye lens. The AttenuRad CT eye shield, which can be used both for adults and pediatrics, is intended to be a single use product because of potential biological contamination.

Model	Description	Price
MO-3000	■ Eye Shield	\$109.90 cs.



Multi-Modality Markers

Product Benefits:

- Provides clear reference points
- Shows on all sequences
- Does not indent into soft tissue
- No artifacts
- 50 Per Container

For further details visit www.mriequip.com or call customer service



Model	Description	Price
BP-1015	■ Multi-Modality Markers	\$100.50 cs.
BP-1037	■ Multi-Modality Markers, Solid	\$105.60 cs.
BP-1030	■ Multi-Modality Markers for Nuclear Medicine/Pet	\$122.10 cs.

"Great Markers! Perfect size and show on both T1 and T2 weighted images. Also, the markers do not indent the skin." - Open MRI Specialists of Stockbridge





MRlequip.com proudly introduces a revolutionary patent pending draping system for Stereotactic Tables.

These drapes will minimize the daily cleaning routine of coils as well as ease the mind of clientele. The drapes provide comfort to patients, while providing protection from contagious diseases such as MRSA. Finally, your staff will be delighted with your decision to simplify clean up between patients, thus achieving a greater throughput.



Accufit Fischer/Siemens Table Drape

25 per Case

CE Certified

Model	Description	Price
BS-1500	Table Drape (25 per case)	\$148.75 cs.



Accufit Fischer/Siemens Autoguide Drape

25 per Case

CE Certified

Model	Description	Price
BS-1501	Autoguide Drape (25 per case)	\$376.25 cs.



Accufit Lorad Table Drape System

20 per Case

Includes table cover and footrest drape

CE Certified

Model	Description	Price
BS-1502	Lorad Table Drape System (20 per case)	\$200.00 cs.




Accufit Lorad Stage Drape System

20 per Case

Includes needle guidance drape, stage drape and control panel cover. The Lorad® Two Piece Stage Drape offers complete coverage including the camera face, stage and needle guide. This option completely protects the patient and the equipment. Worried about MRSA or other bodily fluids? This system is the complete coverage option.

CE Certified

Model	Description	Price
BS-1503	Lorad Stage Drape System (20 per case)	\$220.00 cs.



Accufit Lorad One Piece Stage Drape System

20 per Case

Designed to accommodate the technologist for an easy, quick and cost effective alternative to our BS-1503 Lorad Two Piece Drape System. May be used as a two piece for older systems by tearing the drape at the perforations and attaching on the back of the paddle and the front of the paddle. System also includes the control panel cover.

CE Certified

Model	Description	Price
BS-1506	Lorad One Piece Stage Drape System (20 per case)	\$220.00 cs.





MagnaWand™ Cleaning Wand

A non-magnetic tool designed specifically to clean and disinfect MRI and CAT scan machines. Finally Hospitals will be able to easily clean between patients without losing valuable scanning time.

The cleaning wand reaches inside the MRI magnet bore or CAT scan gantry to clean and disinfect those difficult to reach places in seconds.

Introducing our new wand hook.

Our **Wand hook**, is designed to hold our MagnaWand, it is self adhesive and can be mounted on the MRI magnet, or on the side of a counter, next to our Acrylic dispenser.

Introducing our new MagnaWand Acrylic Dispenser.

The MagnaWand **Acrylic disposable dispenser** is self adhesive, can be mounted or stand alone atop a counter. Easily dispenses one disposable at a time.

Comes with 30 disposables.

Disinfectant cleaner protects against:
MRSA, S.aures, and more



Model	Description	Price
MT-4500	Kit includes one wand, 25 pads and one 24 oz. bottle of spray disinfectant	\$155.00 ea.
MT-4506	Cleaning Pad Dispenser w/30 Pads	\$64.99 ea.
MT-4501	100 Cleaning pads	\$155.00 cs.
MT-4504	400 Cleaning pads	\$580.00 cs.
MT-4505	800 Cleaning pads	\$1,070.00 cs.
MT-4507	MagnaWand Wand Hook	\$4.50 ea.
MT-4503	Cleaning wand only	\$130.00 ea.

SpectraSan 24-Broad Spectrum Disinfectant

SpectraSan 24 disinfects for 24 hours, not 10 minutes like most disinfectants!

Effective against H1N1 Flu Virus (Swine)

Kills MRSA in 2 minutes, Non-Flammable, Odorless, colorless

Disinfects without bleaching, SpectraSan 24 has been one of the most effective and safe disinfectants to use to kill influenza A (H1N1).

Disinfectant, Fungicide, & Virucide

SpectraSan 24 is an EPA registered, broad spectrum, antibacterial, antifungal and antiviral disinfectant that can provide up to 24-hour residual effectiveness. It kills the most dangerous pathogens as well as everyday household germs.

Model	Description	Price
MT-4508	SpectraSan 24, 32oz Spray Bottle	\$15.00 ea.



■ **MR Safe** ▲ **MR Conditional** ● **Not MR Safe**



Maxima 128 Disinfectant Cleaner - Sanitizer

6 products in 1 – Disinfectant Cleaner, Lime & Soap Scum Control, Glass Cleaner, Bowl Cleaner, Brightwork Cleaner and Laundry Sanitizer. Unique “soft-acid”-based, cutting-edge disinfectant cleans hard, non-porous surfaces better than conventional disinfectant cleaners. Unsurpassed efficacy - HBV effective. Staphylocidal, Pseudomonacidal, Virucidal, Fungicidal and Mildewcidal.

Features:

Disinfectant cleaner
Simplifies Cleaning – 6 Products in 1 Bottle
Lemon-Like pH Cleans Surfaces Better than Conventional Disinfectants
Hospital-Grade Disinfectant
Unsurpassed Efficacy

Uses

Hospitals
Nursing Homes
Medical Offices
Veterinary Clinics

Typical Specifications

Typical Dilution: 1:128
Working pH: 4.1–5.5
Color: Green
Odor: Sassafras
Hlth Rtnq (In-Use): 1
USDA: C1
Case of 4 one gallon containers

Model	Description	Price
MT-4509	Maxima 128 (case of 4)	\$174.00 cs.

Janitorial Cart with Mop Bucket and Wringer

This Janitor's cart was specifically designed to be utilized in areas where MRI equipment is being used. The cart is completely non-ferrous and is non-magnetic, so it will not interfere with MRI equipment usage. All components are made from plastics and/or non-ferrous metals(aluminum and stainless steel) to ensure top quality and to be completely nonmagnetic.

Janitorial Cart Specifications:

55" x 30" x 38"
Grey with yellow 25 gallon vinyl bag
MR Conditional at 3-Tesla or less

Mop Bucket and Wringer Specifications:

35 Quart
Yellow
4 - 3" Non-Locking casters



Model	Description	Price
MT-1920	▲ Janitorial Cart Complete	\$504.20 ea.



Multi-Use Floor Sign

This unique and patented design utilizes four separate multilingual message centers built into a single sign, simply flip to the appropriate message center. Messages are imprinted in red and black for increased visibility. Caution, Wet Floor, Area Closed and Out of Order. Will hang on Janitorial Cart MT-1922

Specifications:

2 Per case

26" High

Model	Description	Price
MT-1903	▲ Multi-Use Floor Sign	\$69.00 cs.



Caution Floor Sign

Twenty-six inch high, easel-type, yellow with a two sided imprint. Accepts barrier chain and common barrier tape for connecting multiple signs and cordoning off an area. Will hang on Janitorial Cart MT-1922

Specifications:

6 Per case

26" High

Model	Description	Price
MT-1904	▲ Caution Floor Sign	\$139.08 cs.



Mop Bucket with Wringer

Can be used on Janitorial Cart MT-1922

Specifications:

35 Quart

Yellow

4 - 3" Non-Locking casters

Model	Description	Price
MT-1902	▲ Mop bucket with Wringer	\$144.20 ea.



Janitorial Cart

This Janitor's cart was specifically designed to be utilized in areas where MRI equipment is being used. The cart is completely non-ferrous and is non-magnetic, so it will not interfere with MRI equipment usage. All components are made from plastics and/or non-ferrous metals (aluminum and stainless steel) to ensure top quality and to be completely nonmagnetic.

Specifications:

55" x 30" x 38"

Grey with yellow 25 gallon vinyl bag

MR Conditional at 3-Tesla or less

Model	Description	Price
MT-1922	▲ Janitorial Cart	\$360.00 ea.





MT-1144



MT-1917



MT-1916

MT-1048



MT-1049

MT-7011



MT-7012
MT-7013



MT-1146
MT-1147



Dust Mop Press Wringer

Will not wring 24" mop head

Specifications

Height w/casters: 15½"
Width: 13½"
Length: 23½"

Wringer Specifications

Length: 21½"
Width: 6"

Model	Description	Price
MT-1144	Mop Bucket Press Wringer	\$216.00 ea.

Broom/Mop Holders



Model	Description	Price
MT-1916	Stick-On Holder, Single	\$5.95 ea.
MT-1917	Green Plastic Holder, Multi	\$35.00 ea.

Mop Handle and Heads

Model	Description	Price
MT-1048	Mop Handle Only	\$34.30 ea.
MT-1049	Mop Head Only (6 per case)	\$144.15 cs.

Mop/Duster Frame & Handle

70" Extension Pole

Model	Description	Price
MT-7010	18" Frame	\$87.90 ea.
MT-7011	24" Frame	\$91.00 ea.

Rubbermaid Mop/Duster Frame & Handle

58" Pole, 24" Head has a squeegee

Model	Description	Price
MT-1146	18" Frame	\$80.00 ea.
MT-1147	24" Frame	\$99.75 ea.

Mop/Duster Pads


6 Per Case

Model	Description	Price
MT-7012	6" x 20" (6 per case)	\$54.95 cs.
MT-7013	6" x 26" (6 per case)	\$79.50 cs.



Pliers & Cutters

Fig	Model	Description	Price
A	TL-1759	7" Diagonal Cutting Pliers	\$212.20 ea.
B	TL-1029	7" Long Nose Cutter Pliers	\$175.75 ea.
C	TL-1031	6½" - 5/16" Max Pliers	\$105.65 ea.
C	TL-1032	8" - 1/2" Max Pliers	\$118.55 ea.
D	TL-1027	10" Groove Joint Pliers, 1½" Max	\$199.75 ea.
D	TL-1028	12" Groove Joint Pliers, 2½" Max	\$239.75 ea.
E	TL-1521	8½" Side Cutting Pliers, Lineman's	\$195.65 ea.



Hammers

Model	Description	Price
TL-1040	14" Ball Peen, 1.5 lbs	\$107.25 ea.
TL-1041	14¼" Ball Peen, 2 lbs	\$125.20 ea.
TL-1042	14½" Ball Peen, 2.25 lbs	\$142.80 ea.
TL-1043	14½" Ball Peen, 2.75 lbs	\$145.20 ea.
TL-1037	14" Claw, .75 lbs	\$112.25 ea.
TL-1038	14" Claw, 1 lbs	\$113.75 ea.
TL-1039	14" Claw, 1.25 lbs	\$171.00 ea.



10 Piece Allen Wrench Set

Set includes one each of the following:

1/16", 3/32", 5/32", 7/32", 5/16", 5/64", 1/8", 3/16", 1/4", 3/8".

Model	Description	Price
TL-1013	10 Piece Allen Wrench Set	\$209.50 set



7 Piece Combination Wrench Set

Set includes one each of the following:

3/8", 7/16", 1/2", 9/16", 5/8", 3/4", 7/8".

Model	Description	Price
TL-1012	7 Piece Combination Set	\$495.00 set



Allen Wrench Sets

Model	Description	Price
TL-1722	5 Piece Standard Set, Includes: 9/16", 5/8", 11/16", 3/4", 7/8"	\$588.00 set
TL-1721	13 Piece Standard Set, Includes: 1/16", 5/64", 3/32", 1/8", 5/32", 3/16", 7/32", 1/4", 5/16", 3/8", 7/16", 1/2", 9/16"	\$337.45 set
TL-1746	4 Piece Metric Set, Includes: 12 mm, 14 mm, 17 mm, 19 mm	\$472.75 set
TL-1745	9 Piece Metric Set, Includes: 2 mm, 3 mm, 4 mm, 5 mm, 6 mm, 8 mm, 10 mm, 12 mm, 14 mm	\$355.00 set



3T Safe

Titanium Combination Wrench Sets

Metric

Set includes one each of the following:
 10 mm, 13 mm, 14 mm, 15 mm, 17 mm, 19 mm.

S.A.E.

Set includes one each of the following:
 $\frac{3}{8}$ ", $\frac{7}{16}$ ", $\frac{1}{2}$ ", $\frac{9}{16}$ ", $\frac{5}{8}$ ", $\frac{11}{16}$ ", $\frac{3}{4}$ ".

Individual pieces and other sizes available upon request.

Minimum purchase applies. Please call for pricing.

Model	Description	Price
TL-2231	▲ 6 Piece, Metric Set	\$1,010.00 set
TL-2230	▲ 7 Piece, S.A.E Set	\$1,020.00 set



3T Safe

Titanium Wrench Sets

Metric

Set includes one each of the following:
 6 x 7 mm, 8 x 9 mm, 10 x 12 mm, 11 x 13 mm, 13 x 15 mm, 14 x 16 mm, 17 x 19 mm.

S.A.E.

Set includes one each of the following:
 $\frac{7}{32}$ " x $\frac{1}{4}$ ", $\frac{3}{8}$ " x $\frac{7}{16}$ ", $\frac{1}{2}$ " x $\frac{9}{16}$ ", $\frac{1}{2}$ " x $\frac{11}{16}$ ".

Individual pieces and other sizes available upon request.

Minimum purchase applies. Please call for pricing.

Model	Description	Price
TL-2232	▲ 7 Piece, Metric Set	\$875.00 set
TL-2233	▲ 4 Piece, S.A.E Set	\$477.00 set



Titanium Hex Key Sets

Metric

Set includes one each of the following:
 1.5 mm, 2 mm, 2.5 mm, 3 mm, 4 mm, 5 mm, 6 mm, 8 mm, 10 mm.


S.A.E.

Set includes one each of the following:
 $\frac{1}{16}$ ", $\frac{5}{64}$ ", $\frac{3}{32}$ ", $\frac{7}{64}$ ", $\frac{1}{8}$ ", $\frac{9}{64}$ ", $\frac{5}{32}$ ", $\frac{3}{16}$ ", $\frac{7}{32}$ ", $\frac{1}{4}$ ".

Individual pieces and other sizes available upon request.

Minimum purchase applies. Please call for pricing.

Model	Description	Price
TL-2234	▲ 9 Piece, Metric Set	\$418.65 set
TL-2235	▲ 10 Piece, S.A.E Set	\$325.25 set



Individual Titanium

Pliers, Cutter, Wrench & Wire Stripper

Fig	Model	Description	Price
A	TL-2122	▲ Needle nose Pliers	\$367.20 ea.
B	TL-2167	▲ Side Cutters	\$352.65 ea.
C	TL-2052	▲ Adjustable Wrench	\$448.20 ea.
D	TL-2168	▲ Wire Strippers	\$382.30 ea.



Titanium Phillips 18 Piece T-5 Gyroscan Tool Kit

Includes:

Hex Key Wrenchs: 1.5mm, 3mm, 4mm & 5mm
Open End Wrenchs: 5x5.5mm, 6x7mm, 8x9mm,
Combination Wrenchs: 10mm, 13mm, 17mm & 19mm
Flatside Screwdrivers: 3/32 x 3", 3/16" x 6" & 1/4" x 5"
Also: Needle Nose Pliers, Folding Ruler, Tweezers and a Case

Model	Description	Price
TL-2118	▲ 18 Piece Phillips Tool Kit	\$1,979.65 set



Titanium Siemens 32 Piece Standard Tool Kit

Includes:

Hex Key Wrenchs: 1/16", 3/32", 2mm, 2.5mm, 3mm, 4mm, 5mm, 6mm, 8mm & 10mm
Open End Wrenchs: 6x7mm, 8x9mm, 9x11mm, 10x12mm, 13x15mm, 14x16mm,
17x19mm, 18x20mm, 21x22mm, 24x27mm & 30x32mm
Flatside Screwdrivers: 3/32 x 3", 1/8" x 6", 3/16" x 6" & 5/16" x 6"
Phillips Screwdrivers: #1 x 3", #2 x 6" & #3 x 6"
Also: Needle Nose Pliers, Side Cutter, Tweezers and a Case

Model	Description	Price
TL-2243	▲ 32 Piece Siemens Tool Kit	\$3,566.75 set



Titanium Standard Screwdriver Set

Set includes one each of the following:

3/32" w/ 6" blade, 1/8" w/ 6" blade, 3/16" w/ 6" blade, 1/4" w/ 8" blade.

Titanium Phillips Screwdriver Set

Set includes one each of the following:

3" blade w/ #0 tip, 3" blade w/ #1 tip, 6" blade w/ #2 tip, 6" blade w/ #3 tip, 1 1/2" blade w/ #2 tip


Model	Description	Price
TL-2236	▲ 4 Piece Standard Set	\$186.90 set
TL-2237	▲ 5 Piece Phillips Set	\$315.40 set



Adjustable Pipe Wrench

Aluminum, 3T Safe


Model	Description	Price
TL-1044	▲ Fits Pipe up to 1"	\$310.20 ea.
TL-1045	▲ Fits Pipe up to 1 1/2"	\$360.70 ea.
TL-1046	▲ Fits Pipe up to 2"	\$413.60 ea.
TL-1047	▲ Fits Pipe up to 2 1/2"	\$377.85 ea.
TL-1048	▲ Fits Pipe up to 3 1/2"	\$600.80 ea.
TL-1049	▲ Fits Pipe up to 5"	\$1,161.90 ea.





Titanium Tweezers, 3T Safe

Model	Description	Price
TL-2222	▲ Titanium Tweezers	\$114.30 ea.



Titanium Ruler Set

2 Piece set includes: 6" Ruler and 78" (2m) Folding Fiberglass Ruler

Model	Description	Price
TL-2242	▲ Ruler Set	\$64.75 set



Titanium Speed Handle

Model	Description	Price
TL-2221	▲ 3/8" Drive	\$241.35 ea.



Bronze Sockets for 1/4" Drive, 1.5T Safe Metric

Set includes one each of the following:
 7 mm, 8 mm, 9 mm, 10 mm, 11 mm, 12 mm, 13 mm.

S.A.E.

Set includes one each of the following:

3/16", 7/32", 1/4", 9/32", 5/16", 11/32", 3/8", 7/16".

Individual pieces and other sizes available upon request.
 Minimum purchase applies. Please call for pricing.

Model	Description	Price
TL-2239	▲ 7 Piece, Metric Set	\$287.20 set
TL-2238	▲ 8 Piece, S.A.E Set	\$288.95 set



Bronze Sockets for 3/8" Drive, 1.5T Safe Metric

Set includes one each of the following:
 8 mm, 9 mm, 10 mm, 11 mm, 12 mm, 13 mm, 14 mm, 15 mm.

S.A.E.

Set includes one each of the following:

5/16", 3/8", 7/16", 1/2", 9/16", 5/8", 11/16", 3/4".

Individual pieces and other sizes available upon request. Minimum purchase applies. Please call for pricing.

Model	Description	Price
TL-2241	▲ 8 Piece, Metric Set	\$287.20 set
TL-2240	▲ 8 Piece, S.A.E Set	\$339.70 set



Bronze Universal Joints

Model	Description	Price
TL-2215	▲ 1/4" Drive	\$237.75 ea.
TL-2216	▲ 3/8" Drive	\$268.40 ea.



100% non-magnetic • Hypoallergenic • Non-sparking • Rustproof • 45% Lighter than steel
• Stronger than Stainless Steel • More resistant than chrome Molybdenum • Extremely acid resistant

Titanium Combination Pliers

Model	Description	Price
TL-2503	▲ Combination Pliers, 160 mm	\$364.70 ea.
TL-2504	▲ Combination Pliers, 200 mm	\$418.75 ea.



Titanium Long Nose Pliers

Model	Description	Price
TL-2505	▲ Long Nose Pliers, 180 mm	\$355.25 ea.
TL-2506	▲ Long Nose Pliers, 200 mm	\$388.75 ea.



Titanium Diagonal Cutting Pliers

Model	Description	Price
TL-2507	▲ Cutting Pliers, 160 mm	\$321.80 ea.



Titanium Crimping Pliers

Model	Description	Price
TL-2509	▲ Crimping Pliers, 200 mm	\$418.75 ea.



Titanium Multi-Position Automatic Pliers

Model	Description	Price
TL-2510	▲ Multi-Position Pliers, 250 mm	\$388.75 ea.
TL-2511	▲ Multi-Position Pliers, 375 mm	\$285.75 ea.
TL-2512	▲ Multi-Position Pliers, 180 mm	\$355.25 ea.



Titanium Slip Joint Pliers

Model	Description	Price
TL-2513	▲ Slip Joint Pliers, 200 mm	\$388.75 ea.



Titanium Circlips Pliers

Model	Description	Price
TL-2515	▲ Circlips Pliers, 180 mm	\$321.80 ea.





Titanium Reversible Ratchet Handle

100% non-magnetic • Hypoallergenic • Non-sparking • Rustproof • 45% Lighter than steel • Stronger than Stainless Steel • More resistant than chrome Molybdenum • Extremely acid resistant

Model	Description	Price
TL-2500 ▲	1/4" Ratchet Handle	\$643.60 ea.
TL-2501 ▲	3/8" Ratchet Handle	\$703.65 ea.
TL-2502 ▲	1/2" Ratchet Handle	\$1,174.60 ea.



Titanium Ball Peen Hammer, 450 gr.

100% non-magnetic • Hypoallergenic • Non-sparking • Rustproof • 45% Lighter than steel • Stronger than Stainless Steel • More resistant than chrome Molybdenum • Extremely acid resistant

Model	Description	Price
TL-2514 ▲	Ball Pein Hammer	\$392.10 ea.



4 Piece Set

Screw Driver - 5/16" Tip, 6 1/2" Combination Pliers, 8" Adjustable End Wrench, 9" Crate Opener

Model	Description	Price
TL-1007 ▲	4 Piece Set	\$335.35 set



6 Piece Set

Screw Driver - 5/16" Tip, 6 1/2" Combination Pliers, 8" Adjustable End Wrench, 9" Crate Opener, #2 Phillips Screwdriver, #10 Groove Joint Pliers


Model	Description	Price
TL-1008 ▲	6 Piece Set	\$405.00 set



6 Piece Screwdriver Set

Standard Screwdriver, 5/16" Tip, 3" Blade, Standard Screwdriver, 7/16" Tip, 6" Blade, Standard Screwdriver, 25/64" Tip, 8" Blade, #1 Phillips Screwdriver, #2 Phillips Screwdriver, #3 Phillips Screwdriver

Model	Description	Price
TL-1010 ▲	6 Piece Screwdriver Set	\$229.75 set



Trash Can w/Foot Pedal

Model	Description	Price
MT-1003 ■	8 Gallon, Red - 16 1/4" x 15 3/4" x 17 1/8"	\$146.50 ea.
MT-1004 ■	12 Gallon, Red - 16 1/4" x 15 3/4" x 23 5/8"	\$177.00 ea.
MT-1005 ■	18 Gallon, Red - 19 3/4" x 16 1/8" x 26 1/2"	\$204.00 ea.
MT-1111 ■	8 Gallon, Beige - 16 1/4" x 15 3/4" x 17 1/8"	\$146.50 ea.
MT-1112 ■	12 Gallon, Beige - 16 1/4" x 15 3/4" x 23 5/8"	\$175.00 ea.
MT-1113 ■	18 Gallon, Beige - 19 3/4" x 16 1/8" x 26 1/2"	\$207.30 ea.

■ MR Safe ▲ MR Conditional ● Not MR Safe



17 Piece Toolkit with Bucket

8" Combination Pliers, Shoe Handle Brush, Putty Knife (stiff), Common Knife 5-3/4" Blade, 12" Groove Joint Pliers, 7" Long Nose Pliers w/Cutters, #3 Phillips Screwdriver, Std. Screwdriver, 7/16" Tip, 16 oz Claw Hammer, 8" Adjustable Wrench, 12" Adjustable Wrench, 14" Adjustable Pipe Wrench, Deck Scraper, 1-1/2" x 15", Putty Knife, 9" Crate Opener, 12" Tin Snips, 12" Bung Wrench

Model	Description	Price
TL-1002	▲ 17 Piece Set	\$1,351.15 set



16 Piece Toolkit with Toolbox

8" Combination Pliers, Shoe Handle Brush, Putty Knife (stiff), Common Knife 3-1/8" Blade, #10 Groove Joint Pliers, 6-1/2" Combination Pliers, #2 Phillips Screwdriver, Std. Screwdriver, 7/16" Tip, 16 oz Claw Hammer, 10" Adjustable Wrench, 14" Adjustable Pipe Wrench, Deck Scraper, 2" x 12", 5/16" x 3/8" Box Wrench, 7/16" x 1/2" Box Wrench, 9/16" x 5/8" Box Wrench, 3/4" x 7/8" Box Wrench

Model	Description	Price
TL-1001	▲ 16 Piece Set	\$906.75 set



11 Piece Toolkit with Toolbox

8" Combination Pliers, Shoe Handle Brush, Putty Knife (stiff), Common Knife 3-1/8" Blade, #10 Groove Joint Pliers, #2 Phillips Screwdriver, Std. Screwdriver, 7/16" Tip, 16 oz Claw Hammer, 6" Adjustable Wrench, 10" Adjustable Wrench, 8" Adjustable Pipe Wrench

Model	Description	Price
TL-1000	▲ 11 Piece Set	\$799.50 set



Socket Sets, Ratchets and Accessories

Ratchets: 1/4", 1/2" and 3/4"
Hinged Wrenches: 1/2" and 3/4"
Extensions: 4" and 17"
Universal Swivel Joints for 1/2" Drive

Please call MRlequip for pricing and more Information at 866-529-3526.

SAE Sizes: 3/16", 7/32", 1/4", 9/32", 5/16", 3/8", 7/16", 1/2", 9/16", 19/32", 5/8", 11/16", 3/4", 25/32", 13/16", 7/8", 15/16", 1", 1 1/16", 1 1/8", 1 1/4", 1 5/16", 1 3/8", 1 7/16", 1 1/2", 1 9/16", 1 5/8", 1 11/16", 1 3/4", 1 13/16", 1 7/8" & 2"

Metric Sizes: 4mm, 4.5mm, 5mm, 5.5mm, 6mm, 7mm, 8mm, 9mm, 10mm, 11mm, 12mm, 13mm, 14mm, 15mm, 16mm, 17mm, 18mm, 19mm, 20mm, 21mm, 22mm, 23mm, 24mm, 25mm, 26mm, 27mm, 28mm, 29mm, 30mm, 32mm, 36mm, 38mm, 41mm, 46mm, 50mm, 55mm & 60mm



Aluminum



Fiberglass



MRI Double Sided Step Stool Ladder

Duty Rating: IAA 375, ANSI Sticker: A14.2

Aluminum

375 lb Weight capacity

Model	Description	Price
LA-1012	▲ 4 Foot	\$280.00 ea.
LA-1000	▲ 5 Foot	\$329.80 ea.
LA-1001	▲ 7 Foot	\$440.35 ea.
LA-1002	▲ 6 Foot	\$384.20 ea.
LA-1013	▲ 8 Foot	\$515.80 ea.

Fiberglass

375 lb Weight capacity

Model	Description	Price
LA-1008	▲ 4 Foot	\$324.60 ea.
LA-1009	▲ 5 Foot	\$385.00 ea.
LA-1010	▲ 6 Foot	\$433.00 ea.
LA-1011	▲ 7 Foot	\$469.00 ea.
LA-1014	▲ 8 Foot	\$561.00 ea.

Non-Folding Stand, Double Climbing Side

Aluminum

Model	Description	Price
LA-1005	▲ 16" Non-Folding, 8½ lbs	\$138.00 ea.
LA-1006	▲ 24" Non-Folding, 10 lbs	\$162.00 ea.

18" Folding Stand, Single Climbing Side

18" x 8" platform size

Aluminum

Model	Description	Price
LA-1003	▲ 18" Folding, 9 lbs	\$162.00 ea.
LA-1004	▲ 24" Folding, 11½ lbs	\$209.00 ea.

SafeFlex Stand

16" High

8 1/2 lbs

Aluminum

No Intermediate Steps

Flexible

Model	Description	Price
LA-1007	▲ SafeFlex Stand	\$147.40 ea.



Office Hours: Monday thru Friday 7:00 am to 4:30 pm (CST)

TERMS: Our terms are NET 45 DAYS on all new and current accounts. C.O.D. if not current. Payments become due 45 days after the invoice date, if the order is complete or not. Shipments will be held if an account contains invoices over 45 days past due. The account will go on permanent C.O.D. when there is a history of past dues. Third party collection fees and legal fees must be paid in full before an account can be reopened.

FREIGHT: Freight will be prepaid and added to your invoice. Freight charges are part of your invoice and are due accordingly. Freight charges are not discountable.

PROOF OF DELIVERY: We will provide proof of delivery within 30 days of shipment upon request.

CLAIMS FOR DAMAGE OR SHORTAGE: Must be made within 48 hours of receiving. Shipping containers must be kept if shipment is damaged or short.

SHIPMENTS ARRIVING VIA LTL TRUCK (not UPS or Fedex ground): **The packages must be opened and inspected at time of delivery for any concealed damage.**

PRICES: All prices are subject to change without notice. Prices reflect single unit cost, except where noted.

MINIMUM ORDER: No Order is too small. Excludes tools.

RETURNS: No item will be accepted for return after 30 days from date of purchase. All returns must be labeled with a return authorization number which can be obtained from our customer service department. **Claims for shortage, breakage or any other reason must be made within 48 hours of receipt.** Shipping containers must be kept if shipment is damaged or short. All returns are subject up to a 25% restocking fee. The customer is responsible for all freight charges incurred. **Custom Products are NOT Returnable.**

CANCELLATION CHARGES: Orders may be canceled without charge providing your request is received prior to shipment except custom orders. If the merchandise has been shipped there will be up to a 25% restocking fee, plus all freight charges incurred. Please get cancellation confirmation to avoid these charges.

OTHER ORDER INFORMATION: We reserve the right to split shipments. We sell on first come first serve basis. Where available, substitutions will be shipped if the product ordered is not available or discontinued.



DUNS # 140786828
Cage: 3M5A0

6248 Bittersweet Lane
Nisswa, MN 56468

Toll Free 866-529-3526
Fax Free 888-433-5174

Purchase Orders Accepted



To:

Presorted Standard
U.S. Postage
Paid
Ripon, WI
Permit No. 100

Attn: Director of Radiology or MRI Technologist

**Low Price
Guarantee**



Non-Ferromagnetic Folding Gurney

NEW

Passes through
ferromagnetic
detection
systems.
See page 1
for more details



NEW

MRI Step Stools - Bariatric and Narrow

See pages 17-18 for specifications
on the NEW Step Stools



Fixed Gurney

Now available with green upholstery
See page 9 for more details

NEW

Now Available with
GREEN Upholstery



Visit www.mriequip.com for many more products and options.